


Prosthetics

Upper Limb



PROSTHETICS

Upper Limb

1	1 Information about the catalog	4
	1.1 Otto Bock: About Us	4
	1.2 Information about the catalog	5
	1.3 Ordering	7
2	1.4 Order sheets prosthetic upper limbs	8
	1.5 Arm Prostheses System Overview	20
	1.6 Contact	21
3	2 Axon-Bus® Prosthetic System	22
3	3 MyoBock® Prosthetic system	32
	3.1 Myo Terminal Devices	32
4	3.2 Myo Prosthetic Gloves	74
	3.3 Myo Flexion	82
	3.4 Myo Rotation	88
	3.5 Myo Elbow Components	94
5	3.6 Myo Battery Manager	106
	3.7 Myo Control Elements	118
	3.8 Myo Master	128
6	3.9 Myo Software	132
4	4 Movo Exoskeletal	138
	4.1 Movo Terminal Devices	138
7	4.2 Movo Prosthetic Gloves	170
	4.3 Movo Wrist Units	176
	4.4 Movo Elbow Joints	186
	4.5 Movo Shoulder Joints	202
8		

5	Movo Endoskeletal	206
5.1	Movo Terminal Devices	206
5.2	Movo Prosthetic Gloves	212
5.3	Movo Wrist Units	218
5.4	Movo Elbow Components	222
5.5	Movo Arm Bars	238
5.6	Movo Shoulder Joints	248
6	Physo	254
6.1	Physo Terminal Devices	254
6.2	Physo Wrist Units	304
7	Into	308
7.1	Liner	308
7.2	Body harnesses and accessories	314
8	Miscellaneous	332
8.1	Socket Comfort, Materials & Workplace Arm Prosthetics	332
8.2	Test, Inspection and Adjustment Equipment	350
8.3	Myo Service Parts	358

1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



"Our objective is to offer maximum mobility, independence and normality to people with physical disabilities. User functionality is therefore the most important criterion for measuring our products."

Professor Hans Georg Näder
Chairman and CEO

The name Otto Bock has stood for technology, innovation, quality and customer orientation for more than 90 years. Driven by a pioneering spirit, courage and decisiveness, the prosthetist Otto Bock founded Otto Bock Orthopädische Industrie GmbH in Berlin in 1919. He had the courage to break new ground and set higher standards which would revolutionise an entire industry. Under the leadership of Dr. Max Näder, Otto Bock became a company of international standing. Thanks to his creativity and inventive talent, Max Näder continued to set standards in orthopaedic technology with the development of products, such as the modular leg prosthesis system or myoelectric arm prosthesis. The company began to establish an international sales structure in 1958, when the first foreign subsidiary was founded. After years of consistent and dynamic expansion, Otto Bock is now a true global player and a strong corporate brand. Today our name stands for high-quality, functional and technologically outstanding products and services in orthopaedic and rehabilitation technology around the world. Whatever we do, people are always our number one priority: We are committed to helping them achieve maximum mobility, independence and normality.

Otto Bock HealthCare is a modern, customer and success-oriented company with a long-standing tradition – a global player with local roots. With 45 sales and service companies and export activities in over 140 countries around the world, we are constantly in close contact with our customers. Thanks to this intimate relationship with the market, we understand user needs and customer requirements and integrate them into the products we develop. Although we have created a future-proof organisation with our global network of development, manufacturing and production sites, we remain committed to the German market and our local roots. Duderstadt, located in Germany's Eichsfeld region, is not only where our largest development and production sites are situated, it is also home to the Otto Bock HealthCare headquarters.

We will continue to use our experience and expertise responsibly in order to improve the quality of life of disabled people by providing functional and technologically outstanding solutions in the future: "Quality for life - made by Otto Bock".

Otto Bock offers a comprehensive range of practice-oriented products for upper extremity applications. The proven and tested MyoBock system is available as a 4.8 volt system for children and as a 7.2 volt system for children, youths and adults. New controls and the lithiumion charging technology enlarge the offer for the different applications in fitting practice.

The MyoBoy allows realistic on-site simulation of all control variants. This facilitates the entire fitting process.

MyoBock systems offer dual advantages of high function and natural appearances.

We are in a continuing process of increasing our offer for you. When ordering new components you need to install into a prosthesis in combination with older products, please check whether they are compatible.

Our ordering overview is designed to assist you with the order process. It is organised according to the desired function of the arm prosthesis and designed to correspond with your daily practice. An alphabetically organised "Table of Contents" as well as an index arranged in order number make it easier to use this catalogue.

The various control systems for our System Electric Hands and System Electric Greifers and their corresponding service parts have been printed in different colours to further assist you.

Our Myo specialists are available on our OT hotline' at +49 (0) 5527 848-3333 and will gladly provide you with information and technical tips on other technology or service questions.

Technical changes reserved.



Zeichenerklärung

	Informationsblatt, Poster
	Informationsmaterial
	Verarbeitungs-/ Bedienungsanleitung
	Diese gekennzeichneten Produkte sind in ausgewählten Ländern als Warenzeichen eingetragen.
	Verkleinerung
	Vergrößerung
	Mischungsverhältnis
	Zwei-Komponenten: A-Komponente, B-Komponente
	Selbstklebend
	Beidseitig klebend
	Waschbar bei 40°C im Schonwaschgang
	Waschbar bei 60°C
	646A309 Wissen und Anwendung

Gefahrstoff-Symbolerläuterungen (R-/S-Sätze)*

	Xi Reizend
	Xn Gesundheitsschädlich
	O Brandfördernd
	F Leichtentzündlich
	F+ Hochentzündlich
	N Umweltgefährlich
	C Ätzend
	T Giftig

Gefahrstoff-Symbolerläuterungen (P-/H-Sätze)*

Gefahrenklassen	Gefahrenkategorien
	Entzündbare Gase 1
	Entzündbare Aerosole 1, 2
	Entzündbare Flüssigkeiten 1, 2, 3
	Entzündbare Feststoffe 1, 2
	Selbsterzetzende Stoffe und Gemische Typen B, C, D, E, F
	Pyrophore Flüssigkeiten 1
	Pyrophore Feststoffe 1
	Selbsterhitzungsfähige Stoffe und Gemische 1, 2
	Stoffe und Gemische, die bei Berührung mit Wasser entzündbare Gase abgeben 1, 2, 3
	Organische Peroxide Typen B, C, D, E, F
	Oxidierende Gase 1
	Oxidierende Flüssigkeiten 1, 2, 3
	Oxidierende Feststoffe 1, 2, 3
	Auf Metalle korrosiv wirkend 1
	Hautätzend 1A, 1B, 1C
	Schwere Augenschädigung 1
	Akute Toxizität (oral, dermal, inhalativ) 1, 2, 3
	Akute Toxizität (oral, dermal, inhalativ) 4
	Reizung der Haut 2
	Augenreizung 2
	Sensibilisierung der Haut 1
	Spezifische Zielorgan-Toxizität (einmalige Exposition) 3
	Atemwegsreizung
	Narkotisierende Wirkungen
	Sensibilisierung der Atemwege 1
	Keimzellmutagenität 1A, 1B, 2
	Karzinogenität 1A, 1B, 2
	Reproduktionstoxizität 1A, 1B, 2
	Spezifische Zielorgan-Toxizität (einmalige Exposition) 1, 2
	Spezifische Zielorgan-Toxizität (wiederholte Exposition) 1, 2
	Aspirationsgefahr 1
	Gewässergefährdend
	- Akut gewässergefährdend 1
	- Chronisch gewässergefährdend 1,2

* Die abgedruckten Gefahrstoffsymbole (R-/S-Sätze und P-/H-Sätze) im Katalog entsprechen der Gefahrstoffkennzeichnung zur Zeit der Drucklegung. Diese beziehen sich auf die Materialien im Rohzustand. Änderungen vorbehalten.


Bitte beachten Sie, dass die im Katalog aufgeführten Grundfarben nicht zwingend mit der tatsächlichen Farbwirkung übereinstimmen.

Information about the catalog

6	Physo	286
6.1	Physo Greifgeräte	286
6.2	Physo Handgelenke	336
7	Into	340
7.1	Line	340
7.2	Verhaltenshilfen und Zubehör	346
8	Sonstiges	366
8.1	Schaffensort, Materialien & Arbeitsplatz Anpassbarkeit	366
8.2	Tea-, Prof- und Einzelgeräte	384
8.3	Service-Info Myo	392


Der Katalog ist in **8 Kapitel** gegliedert.

Über Register an der rechten Seite und den Seitentiteln oben auf jeder Seite können Sie sich schnell und einfach zurechtfinden.



Stichwortverzeichnis	Index
A	1
B	2
C	3
D	4
E	5
F	6
G	7
H	8
I	9
J	10
K	11
L	12
M	13
N	14
O	15
P	16
Q	17
R	18
S	19
T	20
U	21
V	22
W	23
X	24
Y	25
Z	26

Das Stichwortverzeichnis enthält sämtliche Produkte in alphabetischer Reihenfolge. Ebenso können Sie über das Indexverzeichnis mittels Kennzeichen die Seite zum gesuchten Produkt auffinden.



Bestellblätter für verschiedene Produkte mit den jeweiligen Faxnummern finden Sie im Kapitel Bestellinformation / Stichwortverzeichnis / Index.

User information about the products

Application examples, product characteristics and advantages

Name of product

8E38=8 SensorHand Speed

With quick-disconnect wrist

Suitable for all amputation levels, except wrist disarticulation. Passive wrist rotation with ratchet lock (can be replaced with the 11S30 Friction Ring). The SensorHand Speed features the automatic Auto-Grasp SUVA sensor technology, the FlexiGrip function as well as different control modes for the fitting with one or two electrodes.

13E184=* Coding Plugs of various colours or the 757T13 MyoSelect are used to select the desired control mode. Please also note the table in the user manual for this.

13E184=1 (white: DMC plus control)

13E184=2 (red: AutoControl LowInput control)

13E184=3 (green: AutoControl)

13E184=4 (blue: VarioControl)

13E184=5 (yellow: VarioDual Control)

13E184=6 (violet: DMC plus Control; SUVA sensors and FlexiGrip can be switched off)

The SensorHand Speed can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

With central coaxial plug connection, automatic shut-off electronics and integrated on-off switch, low-friction bevel gear, positive back lock and system inner hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number

Article number	Side	Size	Inner hand
8E38=8-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E38=8-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E38=8-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E38=8-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E38=8-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E38=8-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4

Technical Data

Reference Number	8E38=8	8E38=8	8E38=8
Size	7 1/4	7 3/4	8 1/4
Operating voltage	6/7.2 V	6/7.2 V	6/7.2 V
Opening width	100 mm	100 mm	100 mm
Proportional gripping force	15-300 mm/sec	15-300 mm/sec	15-300 mm/sec
Proportional speed	0-100 N	0-100 N	0-100 N
Weight with system inner hand for	462 g Women, Adolescents	462 g Men	462 g Men

Technical data such as size, weight, etc.

Image of product



646D165 647H495

- The electrode(s) must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see Page 63, 64, 65, 66

Note

1
2

3

4

5

6

7

8


Below-Elbow to Above-Elbow Residual Limb

Function: Opening and closing using 2 electrodes

SAMPLE

Sample order for a system with EnergyPack 757B20/757B21

For the assembly of the prosthesis, please select only **one** item per box!




System: DMC
VariPlus control unit

8E33=9

System: Digital Twin
control unit

8E33=7



System: VariPlus Speed

L 8E38=9 R

7 1/4

7 3/4

8 1/4

System: SensorHand Speed

L 8E38=8 R

7 1/4

7 3/4

8 1/4

System: DMC plus
control unit

L 8E38=6 R

7

7 1/4

7 3/4

8 1/4

System: Digital Twin
control unit

L 8E38=7 R

7

7 1/4

7 3/4

8 1/4

Coding plug

13E184=2* (red)


13E184=5** (yellow)

13E184=6*** (purple)

*for SensorHand mode:
Auto Control Low Input

** for VarioDual

*** for DMC plus
(sensors can be switched off)




Prosthetic glove
(For order no.,
see catalogue) Quantity 1

8S11N=210x78R4

Electrode cable

13E129=G100

13E129=G300

13E129=G600 2


13E129=G1000

Electrode

13E202 13E200


2 50

60




Tube valve for suction
socket – if 13E202

12V10



Coaxial plug

9E169



Coupling piece


10S4

Lamination ring

10S1=45


10S1=50

10S1=54



Connection cable


757P41



Battery mounting set
screwable

757Z184=2


757Z190=2



EnergyPack


757B20 Quantity 2

757B21 Quantity



Battery charger


757L20



Battery connecting
cable

13E188=200

13E188=600



Battery mounting set

757Z184=1

757Z190=1


MyoEnergy Integral

757B35=1

757B35=3

MyoCharge Integral


757L35



Battery connecting
cable

13E51=2


13E51=4



Battery mounting set


757Z103=1

757Z103=2



X-Change Pack

757B15 Quantity



Pulse battery charger

757L14

Frame Legend

- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology

Sender/customer number


.....

.....

.....

Customer no.:

Signature



System: DMC
VariPlus control unit

8E33=9


System: Digital Twin
control unit

8E33=7

Below-Elbow to Above-Elbow Residual Limb

Function: Opening and closing using 2 electrodes

For the assembly of the prosthesis,
please select only **one** item per box!



System: VariPlus Speed

L	8E38=9	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: SensorHand Speed


L	8E38=8	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: DMC plus control unit

L	8E38=6	R
<input type="checkbox"/>	7	<input type="checkbox"/>
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: Digital Twin control unit

L	8E38=7	R
<input type="checkbox"/>	7	<input type="checkbox"/>
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>




Electrode cable

13E129=G100

13E129=G300


13E129=G600

13E129=G1000




Electrode

13E202	13E200
<input type="checkbox"/>	50 <input type="checkbox"/>
<input type="checkbox"/>	60 <input type="checkbox"/>




Tube valve for suction socket - if 13E202

12V10




Coaxial plug

9E169



Coupling piece

10S4




Lamination ring

10S1=45


10S1=50

10S1=54



Connection cable


757P41



Battery mounting set screwable

757Z184=2


757Z190=2



EnergyPack


757B20 Quantity

757B21 Quantity



Battery charger


757L20



Battery connecting cable

13E188=200


13E188=600



Battery mounting set

757Z184=1


757Z190=1



MyoEnergy Integral

757B35=1

757B35=3



MyoCharge Integral

757L35

Coding plug

13E184=2* (red)



13E184=5** (yellow)

13E184=6*** (purple)

*for SensorHand mode:
Auto Control Low Input

** for VarioDual


*** for DMC plus
(sensors can be switched off)

Battery connecting cable

13E51=2


13E51=4



Battery mounting set


757Z103=1

757Z103=2



X-Change Pack

757B15 Quantity



Pulse battery charger

757L14

Frame Legend

- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology


Sender/customer number

.....

.....

Customer no.:

Signature



Prosthetic glove
(For order no., see catalogue) Quantity

.....

1

Transcarpal residual limb

Function: Opening and closing using 2 electrodes

For the assembly of the prosthesis, please select only **one** item per box!

2

3


4

5

6

7

8




System: DMC plus control unit

L	8E44=6	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: Digital Twin control unit

L	8E44=7	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>




Electrode cable

13E129=G100

13E129=G300


13E129=G600

13E129=G1000



Electrode

13E202	13E200
<input type="checkbox"/>	50 <input type="checkbox"/>
<input type="checkbox"/>	60 <input type="checkbox"/>




Extension cable

9E185=30

9E185=40


9E185=50



Battery mounting set

757Z184=1


757Z190=1



EnergyPack


757B20 Quantity

757B21 Quantity



Battery charger


757L20



Battery connecting cable

13E51=2


13E51=4



Battery mounting set


757Z103=1

757Z103=2




X-Change Pack

757B15 Quantity




Pulse battery charger

757L14




Prosthetic glove
(For order no., see catalogue) Quantity



Distributor


13E190=150

oder



Distributor


13E190



MyoEnergy Integral

757B35=1

757B35=3



MyoCharge Integral

757L35

Frame Legend


- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology

Sender/customer number

.....
.....
.....

Customer no.:

Signature



System: DMC
VariPlus control unit

8E33=9


System: Digital Twin
control unit

8E33=7

Below-Elbow to Above-Elbow Residual Limb

Function: Opening and closing using 1 electrode

For the assembly of the prosthesis,
please select only **one** item per box!



System: VariPlus Speed

L	8E38=9	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: SensorHand Speed

L	8E38=8	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: Digital Twin control unit


L	8E38=7	R
<input type="checkbox"/>	7	<input type="checkbox"/>
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

Coding plug

13E184=2* (red)

13E184=5** (yellow)

*for AutoControl-LowInput
** for VarioDual



Prosthetic glove
(For order no.,
see catalogue) Quantity

Electrode cable

13E129=G100

13E129=G300

13E129=G600


13E129=G1000

Electrode

13E202 13E200


50

60




Tube valve for suction
socket – if 13E202

12V10



Coaxial plug

9E169



Coupling piece


10S4

Lamination ring

10S1=45


10S1=50

10S1=54



Connection cable


757P41



Battery mounting set
screwable

757Z184=2


757Z190=2



EnergyPack


757B20 Quantity

757B21 Quantity



Battery charger


757L20



Battery connecting
cable

13E188=200


13E188=600



Battery mounting set


757Z184=1

757Z190=1




X-Change Pack

757B15 Quantity



Pulse battery charger


757L14



Battery connecting
cable

13E51=2


13E51=4



Battery mounting set

757Z103=1


757Z103=2



MyoEnergy Integral

757B35=1

757B35=3



MyoCharge Integral

757L35

Frame Legend

- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology

Sender/customer number

.....
.....

Customer no.:

Signature

1

2

3

4


5

6

7

8

1



System: DMC
VariPlus control unit

8E34=9

System: Digital Twin
control unit


8E34=7

Wrist disarticulation

Function: Opening and closing using 2 electrodes

For the assembly of the prosthesis,
please select only **one** item per box!

2



System: VariPlus Speed

L 8E39=9 R

7 1/4

7 3/4

8 1/4

System: SensorHand Speed

L 8E39=8 R

7 1/4

7 3/4

8 1/4

System: DMC plus control unit

L 8E39=6 R

7

7 1/4

7 3/4

8 1/4

System: Digital Twin control unit

L 8E39=7 R

7

7 1/4

7 3/4

8 1/4


Electrode cable

13E129=G100

13E129=G300

13E129=G600

13E129=G1000




Electrode

13E202 13E200

50


60



Tube valve for suction
socket – if 13E202

12V10

3




Extension cable

9E185=30

9E185=40


9E185=50



Battery mounting set
screwable

757Z184=2


757Z190=2



EnergyPack

757B20 Quantity


757B21 Quantity



Battery charger

757L20


4



Battery connecting
cable

13E51=2


13E51=4



Battery mounting set


757Z103=1

757Z103=2



X-Change Pack


757B15 Quantity



Pulse battery charger

757L14


5



Distributor


13E190=150

oder



Distributor


13E190



MyoEnergy Integral

757B35=1

757B35=3



MyoCharge Integral

757L35

6

Coding plug

13E184=2* (red)


13E184=5** (yellow)

13E184=6*** (purple)

*for SensorHand mode:
Auto Control Low Input

** for VarioDual

*** for DMC plus
(sensors can be switched off)



Frame Legend

- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology

Sender/customer number


.....

.....

Customer no.:

Signature

8




Prosthetic glove
(For order no.,
see catalogue) Quantity

Transcarpal residual limb


Function: Opening and closing using 1 electrode

For the assembly of the prosthesis,
please select only **one** item per box!



System: Digital Twin control unit

L	8E44=7	R
<input type="checkbox"/>	7	<input type="checkbox"/>
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>




Electrode cable

13E129=G100

13E129=G300


13E129=G600

13E129=G1000



Electrode

13E202	13E200
<input type="checkbox"/>	50 <input type="checkbox"/>
<input type="checkbox"/>	60 <input type="checkbox"/>




Extension cable

9E185=30

9E185=40


9E185=50



Battery mounting set
screwable

757Z184=2


757Z190=2



EnergyPack


757B20 Quantity

757B21 Quantity




Battery charger

757L20




Prosthetic glove
(For order no.,
see catalogue) Quantity



Battery connecting
cable

13E51=2


13E51=4



Battery mounting set


757Z103=1

757Z103=2




X-Change Pack

757B15 Quantity



Pulse battery charger


757L14



Distributor


13E190=150

oder



Distributor


13E190



MyoEnergy Integral

757B35=1

757B35=3



MyoCharge Integral

757L35

Frame Legend

- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology

Sender/customer number

.....

.....

.....

Customer no.:

Signature

1

2

3

4


5

6

7

8

1




System: DMC
VariPlus control unit

8E34=9

System: Digital Twin
control unit

8E34=7

2



System: VariPlus Speed

L	8E39=9	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: SensorHand Speed

L	8E39=8	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: Digital Twin control unit

L	8E39=7	R
<input type="checkbox"/>	7	<input type="checkbox"/>
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

3

4

5

6

7


8

Coding plug

13E184=3* (grün)

13E184=4** (blau)

*for SensorHand mode:
Auto Control Low Input
** for VarioDual



Prosthetic glove
(For order no., see catalogue) Quantity

Electrode cable

13E129=G100

13E129=G300

13E129=G600

13E129=G1000

Extension cable

9E185=30

9E185=40

9E185=50

Battery connecting cable


13E51=2

13E51=4

Distributor

13E190=150

oder



Distributor

13E190

Electrode

13E202	13E200
<input type="checkbox"/>	50 <input type="checkbox"/>
<input type="checkbox"/>	60 <input type="checkbox"/>

Battery mounting set screwable


757Z184=2

757Z190=2

Battery mounting set


757Z103=1

757Z103=2



Tube valve for suction socket – if 13E202


12V10



EnergyPack

757B20 Quantity

757B21 Quantity




X-Change Pack

757B15 Quantity

MyoEnergy Integral


757B35=1

757B35=3



Battery charger

757L20



Pulse battery charger

757L14

MyoCharge Integral

757L35

Wrist disarticulation

Function: Opening and closing using 1 electrode

For the assembly of the prosthesis, please select only **one** item per box!

Frame Legend

- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology


Sender/customer number

.....

.....

Customer no.:

Signature



System: DMC
VariPlus control unit

8E33=9


System: Digital Twin
control unit

8E33=7

Below-Elbow to Above-Elbow Residual Limb

Function: Opening and closing using 2 electrodes
and pronation and supination or MyoRotronic

For the assembly of the prosthesis,
please select only **one** item per box!



System: VariPlus Speed

L 8E38=9 R

7 1/4

7 3/4

8 1/4

System: SensorHand Speed

L 8E38=8 R

7 1/4

7 3/4

8 1/4

System: DMC plus control unit

L 8E38=6 R

7

7 1/4

7 3/4

8 1/4

System: Digital Twin control unit

L 8E38=7 R

7

7 1/4

7 3/4

8 1/4

Electrode cable

13E129=G100

13E129=G300

13E129=G600


13E129=G1000

Electrode

13E202 13E200

50

60



Electric wrist rotator


10S17

Lamination ring

10S1=45


10S1=50

10S1=54




MyoRedronic

13E205




Tube valve for suction socket – if 13E202

12V10



Connection cable


757P41



Battery mounting set screwable

757Z184=2


757Z190=2



EnergyPack


757B20 Quantity

757B21 Quantity



Battery charger


757L20



Battery connecting cable

13E188=200


13E188=600



Battery mounting set


757Z184=1

757Z190=1



X-Change Pack

757B15 Quantity



Pulse battery charger



757L14

Coding plug

13E184=5* (yellow)

13E184=6** (purple)


*for SensorHand mode: VarioDual
** DMC plus (sensors can be switched off)

Battery connecting cable

13E51=2

13E51=4



Battery mounting set

757Z103=1

757Z103=2

MyoEnergy Integral

757B35=1

757B35=3

MyoCharge Integral

757L35

Frame Legend

- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology


Sender/customer number

.....

.....

Customer no.:

Signature



Prosthetic glove
(For order no., see catalogue) Quantity

1

2

3


4

5

6

7

8




System: DMC
VariPlus control unit

8E33=9

System: Digital Twin
control unit

8E33=7



System: VariPlus Speed

L	8E38=9	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: SensorHand Speed

L	8E38=8	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: DMC plus control unit

L	8E38=6	R
<input type="checkbox"/>	7	<input type="checkbox"/>
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: Digital Twin control unit

L	8E38=7	R
<input type="checkbox"/>	7	<input type="checkbox"/>
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

Coding plug

13E184=2* (red)



13E184=5** (yellow)

13E184=6*** (purple)

* for SensorHand mode:
AutoControl-LowInput

** for VarioDual


*** for DMC plus
(sensors can be switched off)

Prosthetic glove
(For order no.,
see catalogue) Quantity

Below-Elbow to Above-Elbow Residual Limb

Function: Opening and closing using 2 electrodes and
pronation and supination using switch




Electrode cable

13E129=G100

13E129=G300


13E129=G600

13E129=G1000




Electrode

13E202	13E200
<input type="checkbox"/>	50 <input type="checkbox"/>
<input type="checkbox"/>	60 <input type="checkbox"/>




Tube valve for suction
socket – if 13E202


12V10




Harness pull switch
9X14




Cable pull switch
9X18



Rocker switch
9X25




Electric wrist rotator
10S17




Lamination ring

10S1=45	<input type="checkbox"/>
10S1=50	<input type="checkbox"/>
10S1=54	<input type="checkbox"/>




Connecting cable

13E50=250	<input type="checkbox"/>
13E50=1200	<input type="checkbox"/>




Connecting cable

13E97=250	<input type="checkbox"/>
13E97=1200	<input type="checkbox"/>




Connection cable
757P41




Battery mounting set
screwable

757Z184=2	<input type="checkbox"/>
757Z190=2	<input type="checkbox"/>




EnergyPack

757B20 Quantity	<input type="checkbox"/>
757B21 Quantity	<input type="checkbox"/>




Battery charger
757L20




Battery connecting
cable

13E188=200	<input type="checkbox"/>
13E188=600	<input type="checkbox"/>




Battery mounting set

757Z184=1	<input type="checkbox"/>
757Z190=1	<input type="checkbox"/>




X-Change Pack

757B15 Quantity	<input type="checkbox"/>
-----------------	--------------------------




Pulse battery charger
757L14




Battery connecting
cable

13E51=2	<input type="checkbox"/>
13E51=4	<input type="checkbox"/>




Battery mounting set

757Z103=1	<input type="checkbox"/>
757Z103=2	<input type="checkbox"/>



MyoEnergy Integral

757B35=1	<input type="checkbox"/>
757B35=3	<input type="checkbox"/>



MyoCharge Integral
757L35

Frame Legend

- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology


Sender/customer number

.....

.....

Customer no.:

Signature




System: DMC
VariPlus control unit
8E33=9

System: Digital Twin
control unit
8E33=7

Below-Elbow to Above-Elbow Residual Limb

Function: Opening and closing using 1 electrode and
pronation and supination using switch



System: VariPlus Speed

L	8E38=9	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: SensorHand Speed

L	8E38=8	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: DMC plus control unit

System: Digital Twin control unit



L	8E38=7	R
<input type="checkbox"/>	7	<input type="checkbox"/>
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

Coding plug


13E184=3* (red)

13E184=4** (yellow)


* for AutoControl
** for VarioControl


Prosthetic glove
(For order no.,
see catalogue) Quantity




Electrode cable
13E129=G100
13E129=G300
13E129=G600
13E129=G1000




Electrode
13E202 13E200
 50
 60




Tube valve for suction
socket - if 13E202
12V10




Harness pull switch
9X14




Electric wrist rotator
10S17




Lamination ring
10S1=45
10S1=50
10S1=54




Connecting cable
13E50=250
13E50=1200




Connecting cable
13E97=250
13E97=1200




Cable pull switch
9X18




Rocker switch
9X25




Connection cable
757P41




Battery mounting set
screwable
757Z184=2
757Z190=2




EnergyPack
757B20 Quantity
757B21 Quantity




Battery charger
757L20




Battery connecting
cable
13E188=200
13E188=600




Battery mounting set
757Z184=1
757Z190=1




X-Change Pack
757B15 Quantity




Pulse battery charger
757L14




Battery connecting
cable
13E51=2
13E51=4



Battery mounting set
757Z103=1
757Z103=2



MyoEnergy Integral
757B35=1
757B35=3



MyoCharge Integral
757L35

Frame Legend

- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology

Sender/customer number

.....
.....

Customer no.:

Signature

1

2

3

4

5

6

7

8

1

2

3

4

5


6

7

8


Below-Elbow to Above-Elbow Residual Limb

Function: Opening and closing using switch



System: DMC
VariPlus control unit
8E33=9

System: Digital Twin
control unit
8E33=7



System: VariPlus Speed

L	8E38=9	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: SensorHand Speed

L	8E38=8	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: DMC plus control unit

System: Digital Twin control unit

L	8E38=7	R
<input type="checkbox"/>	7	<input type="checkbox"/>
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>


Coding plug

13E184=2* (red)


13E184=3** (grün)

13E184=4*** (blau)

* for AutoControl-LowInput
** for AutoControl
*** for VarioControl



Prosthetic glove
(For order no.,
see catalogue) Quantity




Electrode cable

13E129=G100

13E129=G300


13E129=G600

13E129=G1000




Connecting cable


13E99=1200




Linear transducer
9X50




Linear transducer
9X52




Harness pull switch
9X14




Cable pull switch
9X18



Rocker switch
9X25




Pressure switch
9X37




Coaxial plug

9E169



Coupling piece

10S4




Lamination ring

10S1=45


10S1=50

10S1=54



Connection cable


757P41



Battery mounting set
screwable

757Z184=2


757Z190=2




EnergyPack

757B20 Quantity

757B21 Quantity




Battery charger
757L20



Battery connecting
cable

13E188=200


13E188=600



Battery mounting set


757Z184=1

757Z190=1




X-Change Pack

757B15 Quantity




Pulse battery charger
757L14



Battery connecting
cable

13E51=2


13E51=4



Battery mounting set

757Z103=1


757Z103=2



MyoEnergy Integral

757B35=1

757B35=3



MyoCharge Integral
757L35

Frame Legend


- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology

Sender/customer number

.....
.....

Customer no.:

Signature




System: DMC
VariPlus control unit

8E33=9

System: Digital Twin
control unit

8E33=7




Electrode cable

13E129=G100

13E129=G300


13E129=G600

13E129=G1000




Connecting cable

13E99=1200




Harness pull switch

9X14




Cable pull switch

9X18




Rocker switch

9X25




Linear transducer

9X50




Linear transducer

9X52



Tube valve for suction
socket - if 13E202

12V10



System: VariPlus Speed

L	8E38=9	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>


System: SensorHand Speed

L	8E38=8	R
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>

System: DMC plus control unit


System: Digital Twin control unit

L	8E38=7	R
<input type="checkbox"/>	7	<input type="checkbox"/>
<input type="checkbox"/>	7 1/4	<input type="checkbox"/>
<input type="checkbox"/>	7 3/4	<input type="checkbox"/>
<input type="checkbox"/>	8 1/4	<input type="checkbox"/>



Electric wrist rotator

10S17




Lamination ring

10S1=45

10S1=50


10S1=54



Connecting cable

13E50=250


13E50=1200



Connecting cable


13E97=250

13E97=1200




Harness pull switch

9X14




Cable pull switch

9X18




Rocker switch

9X25



Connection cable


757P41



Battery mounting set
screwable

757Z184=2


757Z190=2



EnergyPack


757B20 Quantity

757B21 Quantity



Battery charger


757L20



Battery connecting
cable

13E188=200


13E188=600



Battery mounting set


757Z184=1

757Z190=1



X-Change Pack

757B15 Quantity



Pulse battery charger

757L14

Coding plug

13E184=2* (red)



13E184=3** (grün)

13E184=4*** (blau)

* for SensorHand mode:
AutoControl-LowInput

** for AutoControl


*** for VarioControl

Battery connecting
cable

13E51=2


13E51=4



Battery mounting set

757Z103=1


757Z103=2



MyoEnergy Integral


757B35=1

757B35=3



MyoCharge Integral

757L35



Prosthetic glove
(For order no.,
see catalogue) Quantity

Frame Legend

- black frame suitable for Li-Ion technology
- light brown frame suitable for Ni-Mh technology
- dark brown frame suitable for Li-Pol technology

Sender/customer number

.....

.....

.....

Customer no.:

Signature

1

2

3

4

5


6

7

8

Arm Prostheses System Overview

1
2
3
4
5
6
7
8





CUSTOMER SERVICE is capitalized at Otto Bock. Our competent staff are available to you in person with extensive technical know-how at your disposal, you will learn about the latest developments and provide advice on all matters relating to our products. For complex issues, our product experts and specialists will be in the processing technology. Highly qualified sales representatives help you with specific technical solutions and their implementation on the ground. In addition, we offer comprehensive service and marketing concepts. On www.ottobock.de you can always find the latest product information.

1

2

3

4

5

6

7

8

1

2

3

4

5

6

7

8



1

2

Axon-Bus Prosthetic System

3

4

5

6

7

8

1
2
3
4
5
6
7
8



8K500 Michelangelo Hand with Axon-Bus Prosthetic System

The Axon-Bus System is a new system for transradial fittings. Axon stands for Adaptive exchange of neuroplacement data. The Axon-Bus itself is a new development by Ottobock for the field of exoprosthetics. It was derived from safety-related bus systems in the aviation and motor vehicle industries and is a true innovation of our Research and Development department.

The advantage is that it constitutes a self-contained data transmission system in which all components are optimally adapted to each other. The individual components communicate with each other perfectly, eliminating losses in terms of data transmission, speed and functionality. This results in a clear plus in safety and reliability for the user. Compared to conventional systems, this is achieved by a significant reduction in sensitivity to outside interference.

In combination with the Michelangelo Hand, the Axon-Bus System offers more degrees of freedom than ever before. The user benefits from enhanced hand functionality. The modular prosthetic system can be expanded with additional components in the future. The adaptation of the DynamicArm and ErgoArm, additional hand sizes, electric rotation and flexion for the wrist joint, new electrodes and a new gripper or hook along with the realisation of additional features are all in progress. These components will be harmonised with the Axon-Bus system.

The current Axon-Bus System is suitable for transradial fittings only.

It consists of the following components:

- Michelangelo Hand
- AxonSkin Natural
- AxonWrist
- AxonEnergy Integral
- AxonCharge Integral
- AxonMaster
- AxonSoft

Article number	Side	Size	for
8K500=R-M	Right	M (7 3/4)	Women, men
8K500=L-M	Left	M (7 3/4)	Women, men

⚠ **The individual components of the Axon-Bus prosthetic system cannot be ordered separately at this time, but are only available as the all-inclusive package 8K500=***.

8E500 Michelangelo Hand

The 8E500 Michelangelo Hand is to be used exclusively for exoprosthetic fittings of the upper limbs and serves as an anatomical replacement for the human hand.


The 8E500 Michelangelo Hand can be used for unilateral or bilateral amputees from a transradial amputation level or in the case of dysmelia with forearm fittings.

The Michelangelo Hand cannot be combined with components of the Ottobock MyoBock system or components from other manufacturers.

Article number	Side	Size
8E500=L-M	Left (L)	M (7 3/4)
8E500=R-M	Right (R)	M (7 3/4)
Reference Number 8E500		
Size	M (7 3/4)	
Weight	420 g	
Operating voltage	11,1 V	
Operating temperature	-10 bis +60 °C	
Storage temperature	(-)20 - (+) 40 °C	
Relative humidity	max. 80% non-condensing	
Opening width	120 mm	
Gripping force in Opposition Mode	70 N	
Gripping force in Lateral Mode	60 N	
Gripping force in Neutral Mode	15 N	
Average speed	325 mm/sec	
for	Women, men	



 646D501=EN

 647G587

AxonSkin Natural

The AxonSkin glove serves as structural protection for the Michelangelo prosthetic hand and creates a natural appearance for the user.

It is intended exclusively for use with the Ottobock Michelangelo Hand.


Use the colour pattern set N 646M47 to determine the colour.

The Michelangelo Hand is locked with the help of the hand holding tool (711M64) and hand holding tool (711M1) in order to apply and remove the prosthesis glove.

Article number	Side	Colour	For
8S501=L-M2	Left (L)	2	Men
8S501=R-M2	Right (R)	2	Men
8S501=L-M4	Left (L)	4	Men
8S501=R-M4	Right (R)	4	Men
8S501=L-M6	Left (L)	6	Men
8S501=R-M6	Right (R)	6	Men
8S501=L-M8	Left (L)	8	Men
8S501=R-M8	Right (R)	8	Men
8S501=L-M11	Left (L)	11	Men
8S501=R-M11	Right (R)	11	Men
8S501=L-M16	Left (L)	16	Men
8S501=R-M16	Right (R)	16	Men



 646D501=EN

 647G596

1

2

3

4

5

6

7

8

Article number	Side	Colour	For
8S502=L-M2	Left (L)	2	Women
8S502=R-M2	Right (R)	2	Women
8S502=L-M4	Left (L)	4	Women
8S502=R-M4	Right (R)	4	Women
8S502=L-M6	Left (L)	6	Women
8S502=R-M6	Right (R)	6	Women
8S502=L-M8	Left (L)	8	Women
8S502=R-M8	Right (R)	8	Women
8S502=L-M11	Left (L)	11	Women
8S502=R-M11	Right (R)	11	Women
8S502=L-M16	Left (L)	16	Women
8S502=R-M16	Right (R)	16	Women



8S500 AxonSkin Visual / AxonSkin Black

The AxonSkin Visual and AxonSkin Black are two additional PVC gloves for the Michelangelo Hand.

The Michelangelo Hand is locked with the help of the hand holding tool (711M64) and hand holding tool (711M1) in order to apply and remove the prosthesis glove.

Article number	Side	Colour
8S500=R-M0	Right (R)	Visual
8S500=L-M0	Left (L)	Visual
8S500=R-M20	Right (R)	Black
8S500=L-M20	Left (L)	Black

 646D501=EN

 647G596

10V500 AxonWrist

The 10V500=2 AxonWrist in combination with the Michelangelo Hand facilitates a fitting within the framework of the modular Ottobock Axon-Bus prosthetic system.

The 10V500=2 AxonWrist allows passive flexion/extension and passive pronation/supination. These functions support the user during his or her daily activities and provide the highest rehabilitation value. The 10V500=2 AxonWrist promotes a physiologically correct body posture and reduces the need for unnatural, compensating body movements.

The 10V500=2 AxonWrist is to be used exclusively for exoprosthetic fittings of the upper limbs and enables passive rotation and flexion of the Michelangelo Hand.

In combination with myo-controlled gripping prostheses of the Ottobock Axon-Bus modular prosthetic system, the 10V500=2 AxonWrist can be applied to a transradial amputation level and to cases of dysmelia forearm fittings.

Article number	Side
10V500=2L	Left (L)
10V500=2R	Right (R)

Reference Number	10V500
Weight	ca. 180 g
Operating temperature	-10 bis +60 °C
Storage temperature	(-) 20 - (+) 40 °C
Transport temperature	(-) 20 - (+) 40 °C
Relative humidity	max. 80 %, non-condensing
Rotation, flexion	75° / 4 ratchet positions
Rotation, extension	45° / 3 ratchet positions
Pronation / supination	360 °C



 646D501=EN

 647G589


757B500 AxonEnergy Integral


The 757B500 AxonEnergy Integral serves exclusively to provide power to the Axon-Bus prosthetic system. The battery consists of 3 Li-Ion cells. The integrated electronics protect the battery against short circuits, overvoltage, undervoltage and charging outside the allowable temperature range.


The Axon-Bus cable with the three-pin receptacle is used to exchange data and connects the respective prosthesis components to the battery.

The 757B500 AxonEnergy Integral may only be used in conjunction with components of the Ottobock Axon-Bus prosthetic system. Only the 757L500 AxonCharge Integral may be used for charging.

Article number	757B500
Weight	142 g
Operating temperature	0 bis +60 °C
Storage temperature	(-) 20 - (+) 40 °C
Transport temperature	(-) 20 - (+) 40 °C
Relative humidity	max. 80% non-condensing
Capacity	1500 mAh
Dimensions LxWxH	75 x 60 x 21 mm
Approx. output voltage	11.1 V
Approx. charging time	3.5 h



 646D501=EN

 647G592

1



 646D501=EN

 647G593

757L500 AxonCharge Integral

The 757L500 AxonCharge Integral is designed to charge the 757B500 AxonEnergy Integral of the Ottobock Axon-Bus prosthetic system.

Charging is performed automatically after the charging plug has been connected to the charging receptacle of the 757B500 AxonEnergy Integral. The charging plug is secured to the charging receptacle by the integrated magnet. The special contour of the charging receptacle and charging plug ensures quick, reliable positioning of the two components. LEDs indicate when the charger is ready and the current battery capacity.

The 757L500 AxonCharge Integral is intended exclusively for the purpose of charging the 757B500 AxonEnergy Integral.


Article number	757L500
Operating temperature	0 bis +60 °C
Storage temperature	(-) 20 - (+) 40 °C
Transport temperature	(-) 20 - (+) 40 °C
Relative humidity	max. 80% non-condensing


2

3

4



 646D501=EN

 647G590

13E500 AxonMaster

The 13E500 AxonMaster is the central control unit in the Axon-Bus prosthetic system. The AxonMaster samples the user's control signals and routes them to the corresponding prosthesis components via the Axon-Bus (joint control and power supply connection). This allows the user to control the prosthesis and switch between the prosthesis components. The AxonMaster also controls the data communication process of the Axon-Bus. The 13E500 AxonMaster is intended exclusively for exoprosthetic fittings of the upper limbs and serves to control and coordinate all active Ottobock Axon-Bus prosthesis components within the system.

5 control programmes are available:

- MultiGrip
- DMC LowInput
- Digital
- VarioControl
- DoubleChannel

Article number	13E500
Weight	15 g
Operating voltage	11,1 V
Operating temperature	0 bis +60 °C
Storage temperature	(-) 20 - (+) 60 °C
Transport temperature	(-) 20 - (+) 60 °C
Relative humidity	max. 80% non-condensing
Power supply	757B500
Dimensions	53 x 28 x 9 mm

5

6

7

8

560X500=V1.0 AxonSoft

The 560X500=* AxonSoft software is intended exclusively as a tool for adjusting the 13E500 AxonMaster within the framework of unilateral or bilateral prosthetic fittings with the Michelangelo Hand and the Axon-Bus system components for the upper limb.

The 60X5 BionicLink PC is the only permissible method of data transfer between the system components.

Recommended hardware requirements for a PC with 32-bit or 64-bit platforms

- At least 1 GB free hard drive space
- Graphics card with Open GL support
- Minimum resolution 1024 x 768 at 96 DPI (higher DPI settings can result in errors in the display, depending on the resolution)
- 32-bit colour depth (16.7 million colours)
- CD-ROM drive or DVD-ROM drive
- 1 available USB port (if applicable)
- Mouse and keyboard (if applicable)

Minimum PC hardware requirements

- PC with a Pentium III/1 GHz processor, 32-bit (x86)
- 512 MB RAM (memory)
- 1 GB free hard drive space
- Graphics card with Open GL support
- Resolution 1024 x 768 at 96 DPI
- 32-bit colour depth (16.7 million colours)
- CD-ROM drive
- 1 available USB port (if applicable)
- Mouse and keyboard



646D501=EN



647G591

1

2

3

4

5

6

7

8

1

Accessories



711M64 Donning Tool for Michelangelo

The Donning Tool is essential for the correct application of the AxonSkin glove in order to protect the mechanism of the Michelangelo Hand and for easier handling by the technician.

2

3



711M1 Assembly Tool

One side with M12x1.5 outer thread and other side with M12x1.5 outer thread for Otto Bock System Electric Hands and Michelangelo Hand


4

5

6

7

8



1
2
3
4
5
6
7
8

1
2
3
4
5
6
7
8



Myo Terminal Devices

1

2

3

4

5

6

7

8



Electric Hand 2000

The Otto Bock hand system can be operated for 1 1/2 to 13-year-old children. In the Electric Hand 2000, the thumb and fingers rotate around the same axis, thereby executing a more compact rotation. Therefore, objects can be gripped without unnatural compensating movements of the forearm and upper arm.

As a result, the design of the fingers enables the user to hold writing or eating utensils in a more natural manner. In order to achieve a low weight, the chassis components were covered with flexible plastic covers. Therefore, an inner hand was not required.

For the **sizes 5 1/2, 6 and 6 1/2**, the drive consists of **two motors** with a multi-stage overriding drive. The 1st motor opens and closes the hand at high speed and low strength. Therefore, the object being gripped is handled with a gentle touch. If the object is to be gripped tightly, the 2nd motor is activated and the required gripping strength is established. A slip clutch makes it possible to passively open the hand in an emergency, and protects the mechanical components against overloads.

In contrast to the other hand sizes, **size 5** is operated with **one motor**. Furthermore, the hand can also be passively opened via a spring. This makes gripping training for small children easier. The **Electronic Voluntary Opening Control** was designed for ease in control of the hand for those children who are only able to generate one muscle signal. The **EVO Control** is operated by a single electrode. The controller for the children's system allows for selection between digital, proportional and **EVO-** control options. Its special design even allows for the fitting of children with very long below-elbow residual limbs, while maintaining a fitting length that matches the length of the sound arm.

The 757M11 MyoBoy assists in placing the electrodes and ensuring ideal system arrangement.

8E51 Electric Hand 2000

For 4.8/7.4 Volt MyoBock System

Suitable for virtually all levels of amputation

With 9E347=* Finger Filler, passive hand rotation and central and centralised friction contact, **without** cosmetic glove, **without** controller and **without** lamination ring.

Available in four sizes:

Size 5 for children between 1 1/2 and 3 years old

Size 5 1/2 for children between 3 and 6 years old

Size 6 for children between 5 and 10 years old

Size 6 1/2 for children between 8 and 13 years old


The Electric Hand 2000 can be used for children fittings with 4.8 V and 7.4 V depending on the controller used.

Article number	Side	Size	Lamination Ring	Prosthetic Glove
8E51=L5	Left (L)	5	10S16=34	8S20(N)=136x41L
8E51=L5 1/2	Left (L)	5 1/2	10S16=34	8S20(N)=147x45L
8E51=L6	Left (L)	6	10S16=38	8S20(N)=162x56L
8E51=L6 1/2	Left (L)	6 1/2	10S16=38	8S20(N)=177x64L
8E51=R5	Right (R)	5	10S16=34	8S20(N)=136x41R
8E51=R5 1/2	Right (R)	5 1/2	10S16=34	8S20(N)=147x45R
8E51=R6	Right (R)	6	10S16=38	8S20(N)=162x56R
8E51=R6 1/2	Right (R)	6 1/2	10S16=38	8S20(N)=177x64R

Technical Data

Reference Number	8E51	8E51	8E51	8E51
Size	5	5 1/2	6	6 1/2
Operating voltage	4.8/7.4 V	4.8/7.4 V	4.8/7.4 V	4.8/7.4 V
Approx. average current consumption	200 mA	200 mA	200 mA	200 mA
Electrical "Close" switch at approx.	400 mA	400 mA	400 mA	400 mA
Electrical "Open" switch at approx.	400 mA	250 mA	250 mA	250 mA
Opening width	28 mm	37 mm	52 mm	58 mm
Approx. minimum grip force	8 N	25 N	35 N	35 N
Weight (without prosthetic glove)	86 g	115 g	125 g	130 g

• For suitable system prosthetic glove, see pages 76-77, 79



646D326
646D442

647H58

Spare parts for Children's Hand 2000

1



9E342 Finger Cover

Article number	For hand size
9E342=5	5
9E342=5 1/2	5 1/2
9E342=6	6
9E342=6 1/2	6 1/2

2

Technical Data

Reference Number	9E342
Material	Made of thermoplastic material
Consists of:	consists of thumb, index finger and middle finger cover

3



9E347 Finger Filler

Used to fill the little and ring fingers in the 8S20 Prosthetic Glove.

Article number	Side	For size
9E347=L5	Left (L)	5
9E347=R5	Right (R)	5
9E347=L5 1/2	Left (L)	5 1/2
9E347=R5 1/2	Right (R)	5 1/2
9E347=L6	Left (L)	6
9E347=R6	Right (R)	6
9E347=L6 1/2	Left (L)	6 1/2
9E347=R6 1/2	Right (R)	6 1/2

4


5

6



7

8





SensorHand Speed

When it comes to gripping speed and the ability to control it, the new SensorHand Speed beats all Electric Hands available up to now by lengths: it offers a maximum gripping speed of 300 mm/sec (up to now: 130 mm/sec)! Nevertheless, the SensorHand Speed's control is particularly easy and precise thanks to new intelligent software and modified signal processing.

The SensorHand Speed with SUVA* Sensor Technology improves secure grasp and the ability to hold on to objects. Via the sensor integrated into the thumb and a measurement bracket, the SensorHand Speed recognises when the gripped object changes position and then automatically readjusts the gripping force. Therefore, the user does not need to constantly monitor the gripped object in order to readjust manually. This improvement in gripping security provides effective everyday relief.

Various control options using one or two electrodes are available. Simply exchange the 13E184=* Coloured Coding Plugs or use the 757T13 MyoSelect to determine the control mode. There's peace of mind knowing the Otto Bock SensorHand Speed makes it easier to grasp fragile objects or liquid-filled containers.

Developed in cooperation with the **Swiss Insurance Agency**

8E38=8 SensorHand Speed

With quick-disconnect wrist

Suitable for all amputation levels, except wrist disarticulation. Passive wrist rotation with ratchet lock (can be replaced with the 11S30 Friction Ring). The SensorHand Speed features the automatic Auto-Grasp SUVA sensor technology, the FlexiGrip function as well as different control modes for the fitting with one or two electrodes.

13E184=* Coding Plugs of various colours or the 757T13 MyoSelect are used to select the desired control mode. Please also note the table in the user manual for this.

13E184=1 (white: DMC plus control)

13E184=2 (red: AutoControl LowInput control)

13E184=3 (green: AutoControl)

13E184=4 (blue: VarioControl)


13E184=5 (yellow: VarioDual Control)


13E184=6 (violet: DMC plus Control; SUVA sensors and FlexiGrip can be switched off)

The SensorHand Speed can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

With central coaxial plug connection, automatic shut-off electronics and integrated on-off switch, low-friction bevel gear, positive back lock and system inner hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.



 646D165

 647H495



Article number	Side	Size	Inner hand
8E38=8-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E38=8-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E38=8-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E38=8-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E38=8-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E38=8-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4

Technical Data

Reference Number	8E38=8	8E38=8	8E38=8
Size	7 1/4	7 3/4	8 1/4
Operating voltage	6/7.2 V	6/7.2 V	6/7.2 V
Opening width	100 mm	100 mm	100 mm
Proportional gripping force	15-300 mm/sec	15-300 mm/sec	15-300 mm/sec
Proportional speed	0-100 N	0-100 N	0-100 N
Weight with system inner hand	462 g	462 g	462 g
for	Women, Adolescents	Men	Men

- The electrode(s) must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80



 646D165  647H495

8E39=8 SensorHand Speed

With lamination ring

For wrist disarticulation. Passive wrist rotation with friction. The SensorHand Speed features the automatic Auto-Grasp SUVA sensor technology, the FlexiGrip function as well as different control modes for the fitting with one or two electrodes.

13E184=* Coding Plugs of various colours or the 757T13 MyoSelect are used to select the desired control mode. Please also note the table in the user manual for this.

13E184=1 (white: DMC plus control)

13E184=2 (red: AutoControl LowInput control)

13E184=3 (green: AutoControl)

13E184=4 (blue: VarioControl)

13E184=5 (yellow: VarioDual Control)

13E184=6 (violet: DMC plus Control; SUVA sensors and FlexiGrip can be switched off)

The SensorHand Speed can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a central flat cable, automatic shut-off electronics and integrated on-off switch, low friction bevel gear, positive back lock, and System Inner Hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand
8E39=8-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E39=8-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E39=8-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E39=8-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E39=8-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E39=8-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4

Technical Data

Reference Number	8E39=8	8E39=8	8E39=8
Size	7 1/4	7 3/4	8 1/4
Operating voltage	6/7.2 V	6/7.2 V	6/7.2 V
Opening width	100 mm	100 mm	100 mm
Proportional gripping force	0-100 N	0-100 N	0-100 N
Proportional speed	15-300 mm/sec	15-300 mm/sec	15-300 mm/sec
Weight with system inner hand	462 g	462 g	462 g
for	Women, Adolescents	Men	Men

- The electrode(s) must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80

8E41=8 SensorHand Speed

With threaded stud M12x1.5

Suitable for all amputation levels, except wrist disarticulation.

The SensorHand Speed features the automatic Auto-Grasp SUVA sensor technology, the FlexiGrip function as well as different control modes for the fitting with one or two electrodes.

13E184=* Coding Plugs of various colours or the 757T13 MyoSelect are used to select the desired control mode. Please also note the table in the user manual for this.

13E184=1 (white: DMC plus control)

13E184=2 (red: AutoControl LowInput control)

13E184=3 (green: AutoControl)

13E184=4 (blue: VarioControl)

13E184=5 (yellow: VarioDual Control)

13E184=6 (violet: DMC plus Control; SUVA sensors and FlexiGrip can be switched off)

The SensorHand Speed can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.


It features a central cable outlet with automatic shut-off electronics and integrated on-off switch, low friction bevel gear, positive back lock, and System Inner Hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand
8E41=8-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E41=8-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E41=8-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E41=8-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E41=8-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E41=8-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4

Technical Data

Reference Number	8E41=8	8E41=8	8E41=8
Size	7 1/4	7 3/4	8 1/4
Operating voltage	6/7.2 V	6/7.2 V	6/7.2 V
Opening width	100 mm	100 mm	100 mm
Proportional gripping force	0-100 N	0-100 N	0-100 N
Proportional speed	15-300 mm/sec	15-300 mm/sec	15-300 mm/sec
Weight with system inner hand	462 g	462 g	462 g
for	Women, Adolescents	Men	Men

- The electrode(s) must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80



646D165



647H495

1

2

3

4

5

6

7

8

Control Modes for the SensorHand Speed

Program 1	OPEN	CLOSE	Indication
DMC plus Sensor White coding plug Two electrodes	Myo-signal to the electrode Proportional speed	Myo-signal to the electrode Proportional speed	Preferred program for amputees with two strong muscle signals
Program 2	OPEN	CLOSE	Indication
AutoControl LowInput Red coding plug Two electrodes	Myo-signal to the electrode Proportional speed	Myo-signal to the electrode digitally (short signal at any level) Constant speed	For amputees with two weak muscle signals
AutoControl LowInput Red coding plug One electrode & one optional Otto Bock Switch	Myo-signal to the electrode Proportional speed	Signal from the switch Constant speed	For amputees with only 1 muscle and weak muscle signal
AutoControl LowInput Red coding plug Any Otto Bock switch	Hand opens as long as opening side of the switch is activated Constant speed	Signal from closing side of switch: Hand closes Constant speed	For patients with too weak or no muscle signal
Program 3	OPEN	CLOSE	Indication
AutoControl Green coding plug One electrode	Quick, sustained myo-signal to the electrode Constant speed	Very slow muscle relaxation: Hand stays open: Quick muscle relaxation: Hand closes Constant speed	For amputees with only 1 muscle and very weak muscle signal
AutoControl Green coding plug Any Otto Bock switch	Hand opens as long as switch is activated Constant speed	Hand closes automatically as soon as switch is released Constant speed	For amputees with weak or no muscle signals
Program 4	OPEN	CLOSE	Indication
VarioControl Blue coding plug One electrode	Speed and strength of muscle tension to the electrode Proportional speed	Speed and strength of muscle relaxation to electrode Proportional speed	For amputees with 1 muscle and strong uscle signal or tendency to co-contraction
Program 4	OPEN	CLOSE	Indication
VarioControl Blue coding plug One linear transducer	Speed and strength of pull to linear transducer Proportional speed	Speed of release of the linear transducer Proportional speed	For amputees with weak or no muscle signal
Program 5	OPEN	CLOSE	Indication
VarioDual Yellow coding plug Two electrodes	Speed and strength of muscle tension to first electrode Proportional speed	Speed and strength of muscle relaxation to first electrode Proportional speed. Grip force proportional to strength of muscle signal to second electrode.	Control for amputees with 2 strong muscle signals
Program 6	OPEN	CLOSE	Indication
DMC plus Sensor Purple coding plug Two electrodes	Myo-signal to electrode Proportional speed	Myo-signal to electrode Proportional speed	Program for amputees with 2 strong muscle signals SUVA Sensor and FlexiGrip can be switched off.

8E38=9 MyoHand VariPlus Speed

With quick-disconnect wrist

The MyoHand VariPlus Speed is a new development from Ottobock. It combines the mechanical characteristics of the SensorHand Speed and the control options of the DMC VariPlus System Electric Greifer. Thanks to the high gripping strength (approx. 100 N) and speed (up to 300 mm/s), objects can be gripped quickly and precisely. With the 757T13 MyoSelect, a total of 6 different programmes can be selected and adapted according to patient indications. This permits optimum adaptation to the needs and abilities of the prosthesis wearer.

The following control modes can be selected with the 757T13 MyoSelect:

DMC plus
AutoControl LowInput
VarioControl
VarioDual
DigitalControl
DoubleChannel Control

The MyoHand VariPlus Speed can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

With central coaxial plug connection, automatic shut-off electronics and integrated on-off switch, low-friction bevel gear, positive back lock and system inner hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand	For
8E38=9-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4	Women, Adolescent
8E38=9-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4	Men
8E38=9-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4	Men
8E38=9-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4	Women, Adolescent
8E38=9-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4	Men
8E38=9-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4	Men


Technical Data

Reference Number	8E38=9	8E38=9	8E38=9
Size	7 1/4	7 3/4	8 1/4
Static current	1 mA	1 mA	1 mA
Operating temperature	0 bis +70 °C	0 bis +70 °C	0 bis +70 °C
Opening width	100 mm	100 mm	100 mm
Proportional gripping force	0 - 100 N	0 - 100 N	0 - 100 N
Proportional speed	15 - 300 mm/sec	15 - 300 mm/sec	15 - 300 mm/sec
Weight with system inner hand	460 g	125 g, 460 g	460 g

- The electrode(s) must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80



 646D165

 647G504



646D165 647G504

8E39=9 MyoHand VariPlus Speed

With lamination ring

The MyoHand VariPlus Speed is a new development from Ottobock. It combines the mechanical characteristics of the SensorHand Speed and the control options of the DMC VariPlus System Electric Greifer. Thanks to the high gripping strength (approx. 100 N) and speed (up to 300 mm/s), objects can be gripped quickly and precisely. With the 757T13 MyoSelect, a total of 6 different programmes can be selected and adapted according to patient indications. This permits optimum adaptation to the needs and abilities of the prosthesis wearer.

The following control modes can be selected with the 757T13 MyoSelect:

DMC plus
AutoControl LowInput
VarioControl
VarioDual
DigitalControl
DoubleChannel Control

The MyoHand VariPlus Speed can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a central flat cable, automatic shut-off electronics and integrated on-off switch, low friction bevel gear, positive back lock, and System Inner Hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand	For
8E39=9-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4	Women, Adolescent
8E39=9-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4	Men
8E39=9-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4	Men
8E39=9-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4	Women, Adolescent
8E39=9-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4	Men
8E39=9-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4	Men

Technical Data

Reference Number	8E39=9	8E39=9	8E39=9
Size	7 1/4	7 3/4	8 1/4
Static current	1 mA	1 mA	1 mA
Operating temperature	0 bis +70 °C	0 bis +70 °C	0 bis +70 °C
Opening width	100 mm	100 mm	100 mm
Proportional gripping force	0 - 100 N	0 - 100 N	0 - 100 N
Proportional speed	15 - 300 mm/sec	15 - 300 mm/sec	15 - 300 mm/sec
Weight with system inner hand	460 g	460 g	460 g

- The electrode(s) must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80

8E41=9 MyoHand VariPlus Speed

With threaded stud M12x1.5

The MyoHand VariPlus Speed is a new development from Ottobock. It combines the mechanical characteristics of the SensorHand Speed and the control options of the DMC VariPlus System Electric Greifer. Thanks to the high gripping strength (approx. 100 N) and speed (up to 300 mm/s), objects can be gripped quickly and precisely. With the 757T13 MyoSelect, a total of 6 different programmes can be selected and adapted according to patient indications. This permits optimum adaptation to the needs and abilities of the prosthesis wearer.

The following control modes can be selected with the 757T13 MyoSelect:

DMC plus
AutoControl LowInput
VarioControl
VarioDual
DigitalControl
DoubleChannel Control

The MyoHand VariPlus Speed can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a central cable outlet with automatic shut-off electronics and integrated on-off switch, low friction bevel gear, positive back lock, and System Inner Hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.



646D165



647G504

Article number	Side	Size	Inner hand	For
8E41=9-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4	Women, Adolescent
8E41=9-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4	Men
8E41=9-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4	Men
8E41=9-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4	Women, Adolescent
8E41=9-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4	Men
8E41=9-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4	Men

Technical Data

Reference Number	8E41=9	8E41=9	8E41=9
Size	7 1/4	7 3/4	8 1/4
Static current	1 mA	1 mA	1 mA
Operating temperature	0 bis +70 °C	0 bis +70 °C	0 bis +70 °C
Opening width	100 mm	100 mm	100 mm
Proportional gripping force	0 - 100 N	0 - 100 N	0 - 100 N
Proportional speed	15 - 300 mm/sec	15 - 300 mm/sec	15 - 300 mm/sec
Weight with system inner hand	460 g	460 g	460 g

- The electrode(s) must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80

Control Modes for the MyoHand VariPlus Speed®

Control mode 1: DMC plus®

Control with 2 electrodes

Control mode 1	OPEN	CLOSE	Indication
DMC plus®	Sustained electrode signal Speed: proportional Adjustment regulator A	Sustained electrode signal Grip force: proportional After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed® with unwanted electrode signals is prevented. Speed: proportional Adjustment regulator A	For patients with two strong electrode signals

Control mode 2: AutoControl – Low Input

Control with 2 electrodes

Control mode 2	OPEN	CLOSE	Indication
AutoControl LowInput	Sustained electrode signal Reduced range of proportionality: Maximum speed once the LOW-threshold has been reached Speed: proportional Adjustment regulator A	Sustained electrode signal Grip force: Time proportional After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed® with unwanted electrode signals is prevented. Speed: constant Adjustment regulator B	For amputees with two weak electrode signals

Control with 1 electrode and 1 switch

Control mode 2	OPEN	CLOSE	Indication
AutoControl LowInput	Sustained electrode signal Reduced range of proportionality: Maximum speed once the LOW-threshold has been reached Speed: proportional Adjustment regulator A	Signal from the switch Grip force: Time proportional After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed® with unwanted electrode signals is prevented. Speed: constant Adjustment regulator B	For amputees with only one muscle and weak electrode signal

Control with 1 switch

Control mode 2	OPEN	CLOSE	Indication
AutoControl LowInput	MyoHand VariPlus Speed® opens as long as the OPEN-side of the switch is operated. Speed: constant Adjustment regulator A	MyoHand VariPlus Speed® closes as long as the CLOSE-side of the switch is operated. Grip force: Time proportional Speed: constant Adjustment regulator B	For amputees with weak or no electrode signal

Control mode 3: VarioControl

Control with 1 electrode

Control mode 3	OPEN	CLOSE	Indication
VarioControl	<p>Increasing electrode signal through muscle contraction</p> <p>Speed and strength of muscle contraction to the electrode</p> <p>Speed: proportional Adjustment regulator A</p>	<p>Declining electrode signal through muscle relaxation</p> <p>Grip force: Proportional to the decline of the electrode signal After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed® with unwanted electrode signals is prevented.</p> <p>Speed: proportional Adjustment regulator A</p>	For amputees with one strong electrode signal or tendency to co-contraction

Control with 1 linear transducer

Control mode 3	OPEN	CLOSE	Indication
VarioControl	<p>Speed and strength of pull on the linear transducer</p> <p>Speed: proportional Adjustment regulator A</p>	<p>Speed of release of pull on the linear transducer</p> <p>Grip force: Proportional to the release of the pull on the linear transducer After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed® with unwanted electrode signals is prevented.</p> <p>Speed: proportional Adjustment regulator A</p>	For amputees with weak or no electrode signal

Control mode 4: VarioDual

Control with 2 electrodes

Control mode 4	OPEN	CLOSE	Indication
VarioDual	<p>Increasing electrode signal through muscle contraction to the first electrode</p> <p>Speed and strength of muscle contraction to the electrode</p> <p>Speed: Proportional to the increase of the electrode signal Adjustment regulator A</p>	<p>Declining electrode signal through muscle relaxation to the first electrode or: Sustained electrode signal to the second electrode</p> <p>Grip force: Proportional to the signal strength to the second electrode. If only the first electrode is used, the MyoHand VariPlus Speed® closes up to the minimal grip force.</p> <p>Grip force increase: By another, stronger electrode signal to the second electrode. After gripping once with maximum force, the EMG signal required to OPEN the Hand will be set to a higher value. Opening the MyoHand VariPlus Speed® with unwanted electrode signals is prevented.</p> <p>Speed: Proportional to the decline of the electrode signal to the first electrode/proportional to the signal strength to the second electrode Adjustment regulator A</p>	For amputees with two strong electrode signals

1

Control mode 5: Digital Control

Control with 2 electrodes

Control mode 5	OPEN	CLOSE	Indication
Digital Control	Sustained electrode signal Adjustment regulator B	Sustained electrode signal Grip force: Duration of the signal Adjustment regulator B	For amputees with two weak electrode signals

2

Control with 1 electrode and 1 switch

Control mode 5	OPEN	CLOSE	Indication
Digital Control	Sustained electrode signal Adjustment regulator B	Signal from the switch Grip force: Duration of the signal Adjustment regulator B	For amputees with only one muscle and weak electrode signal

3

Control with 1 switch

Control mode 5	OPEN	CLOSE	Indication
Digital Control	MyoHand VariPlus Speed® opens as long as the OPEN-side of the switch is operated. Adjustment regulator B	MyoHand VariPlus Speed® closes as long as the CLOSE-side of the switch is operated. Grip force: Duration of the signal Adjustment regulator B	For amputees with weak or no electrode signals

4

Control mode 6: Double Channel Control

Control with 1 electrode


Control mode 6	OPEN	CLOSE	Indication
Double Channel Control	Quick, strong electrode signal that reaches the upper threshold within 80 ms and is maintained above the upper threshold for at least 30 ms Adjustment regulator B	Slow, gentle electrode signal that does not reach the upper threshold within 80 ms Grip force: Duration of the signal Adjustment regulator B	For amputees with one strong electrode signal


5

6

7

8





System Electric Hand DMC plus

Myoelectrically-controlled arm prostheses are more powerful thanks to **microcontroller technology**.

In the **Dynamic Mode Control plus**, **two independent measurement and regulating systems** proportionally control gripping speed and force at all opening widths.

Gripping speed and gripping force are determined by the strength of the muscle signal.

As soon as the strength of the muscle signal changes, gripping speed and gripping force immediately adapt to the new muscle signal.

The gripping speed is variable from 15 to 130 mm/sec. The high maximum gripping speed, the quickly learned system and the optimal balance of gripping force and speed make it possible for the patient to experience **directly controlled, physiological gripping**. The DMC plus control features both the well-known DMC control and the DMC plus mode. The new DMC plus mode also includes an integrated virtual "hand switch": after gripping once with maximum, a slightly higher signal is required to open the hand again. This improves gripping security and, for example, makes eating easier.

The DMC plus mode is easily activated by removing the function plug integrated in the electronics.

8E38=6 System Electric Hand DMC plus

With quick-disconnect wrist

Suitable for all amputation levels, except wrist disarticulation.

Passive wrist rotation with ratchet lock (can be replaced by the 11S30 Friction Ring).

The DMC plus control features a DMC and DMC plus control mode. The desired control mode is selected with the integrated 13E185 Function Plug. In DMC plus control mode, after gripping once with maximum grip force, a higher signal is required to open the hand. This reduces the risk of opening the hand with undesired muscle signals. In this system, two independent measurement and control systems proportionally control gripping speed as well as gripping force. Gripping speed and gripping force are determined by the strength of the muscle signal. The System Electric Hand DMC plus can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a central coaxial plug connection, automatic shut-off electronics and integrated on-off switch, low-friction bevel gear, positive back lock and system inner hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand
8E38=6-L7	Left (L)	7	8X18=L7
8E38=6-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E38=6-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E38=6-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E38=6-R7	Right (R)	7	8X18=R7
8E38=6-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E38=6-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E38=6-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4


Technical Data

Reference Number	8E38=6	8E38=6
Size	7	7 1/4, 7 3/4, 8 1/4
Operating voltage	6/7.2 V	6/7.2 V
Opening width	79 mm	100 mm
Proportional gripping force	0-90 N	0-90 N
Proportional speed	15-130 mm/sec	15-130 mm/sec
Weight with system inner hand	355 g	457 g
for	Women, Adolescents	Women, Adolescents, Men

- The electrodes must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80



 646D44

 647H326

Particularly short and light: the new Electric Hand size 7.

It closes the gap between the Electric Hand 2000 hand system for children and the well-known System Electric Hands for adults. It is particularly suitable for the fitting of adolescents or women with small, dainty hands.

The new Electric Hand size 7 is offered with the well-known Digital Twin and DMC plus control systems.

1
2
3
4
5
6
7
8



646D44 647H326

8E39=6 System Electric Hand DMC plus

With lamination ring

Suitable for wrist disarticulation.

Passive wrist rotation with friction.

The DMC plus control features a DMC and DMC plus control mode. The desired control mode is selected with the integrated 13E185 Function Plug. In DMC plus control mode, after gripping once with maximum grip force, a higher signal is required to open the hand. This reduces the risk of opening the hand with undesired muscle signals. In this system, two independent measurement and control systems proportionally control gripping speed as well as gripping force. Gripping speed and gripping force are determined by the strength of the muscle signal. The System Electric Hand DMC plus can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a central flat cable, automatic shut-off electronics and integrated on-off switch, low friction bevel gear, positive back lock, and System Inner Hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand
8E39=6-L7	Left (L)	7	8X18=L7
8E39=6-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E39=6-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E39=6-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E39=6-R7	Right (R)	7	8X18=R7
8E39=6-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E39=6-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E39=6-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4

Technical Data

Reference Number	8E39=6	8E39=6
Size	7	7 1/4, 7 3/4, 8 1/4
Operating voltage	6/7.2 V	6/7.2 V
Opening width	79 mm	100 mm
Proportional gripping force	0-90 N	0-90 N
Proportional speed	15-130 mm/sec	15-130 mm/sec
Weight with system inner hand	355 g	457 g
for	Women, Adolescents	Women, Adolescents, Men

- The electrodes must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80

Particularly short and light: the new Electric Hand size 7.

It closes the gap between the Electric Hand 2000 hand system for children and the well-known System Electric Hands for adults. It is particularly suitable for the fitting of adolescents or women with small, dainty hands.

The new Electric Hand size 7 is offered with the well-known Digital Twin and DMC plus control systems.

8E41=6 System Electric Hand DMC plus

With threaded stud M12x1.5

Suitable for all amputation levels, except wrist disarticulation.

The DMC plus control features a DMC and DMC plus control mode. The desired control mode is selected with the integrated 13E185 Function Plug. In DMC plus control mode, after gripping once with maximum grip force, a higher signal is required to open the hand. This reduces the risk of opening the hand with undesired muscle signals. In this system, two independent measurement and control systems proportionally control gripping speed as well as gripping force. Gripping speed and gripping force are determined by the strength of the muscle signal. The System Electric Hand DMC plus can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a central cable outlet with automatic shut-off electronics and integrated on-off switch, low friction bevel gear, positive back lock, and System Inner Hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand
8E41=6-L7	Left (L)	7	8X18=L7
8E41=6-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E41=6-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E41=6-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E41=6-R7	Right (R)	7	8X18=R7
8E41=6-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E41=6-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E41=6-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4


Technical Data

Reference Number	8E41=6	8E41=6
Size	7	7 1/4, 7 3/4, 8 1/4
Operating voltage	6/7.2 V	6/7.2 V
Opening width	100 mm	100 mm
Proportional gripping force	0-90 N	0-90 N
Proportional speed	15-130 mm/sec	15-130 mm/sec
Weight with system inner hand	355 g	457 g
for	Women, Adolescents	Women, Adolescents, Men

- The electrodes must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80



 646D44

 647H326

Particularly short and light: the new Electric Hand size 7.

It closes the gap between the Electric Hand 2000 hand system for children and the well-known System Electric Hands for adults. It is particularly suitable for the fitting of adolescents or women with small, dainty hands.

The new Electric Hand size 7 is offered with the well-known Digital Twin and DMC plus control systems.

1

2

3

4

5



System Electric Hand Digital Twin

The System Electric Hand Digital Twin combines both the classic digital control and double-channel control in a single hand.

The orthopaedic technician decides whether one or two electrodes will be used to control the prosthesis. The desired control mode is selected with the integrated function plug. This way the orthopaedic technician may decide during trial fitting which control mode is best for the patient.

The System Electric Hand Digital Twin can be operated with various interchangeable battery technologies.

6

7

8

8E38=7 System Electric Hand Digital Twin

With quick-disconnect wrist

Suitable for all amputation levels, except wrist disarticulation.

Passive wrist rotation with ratchet lock (can be replaced by the 11S30 Friction Ring). The Digital Twin control features a digital and dual channel control mode. The desired control mode is selected with the integrated 13E185 function plug. The System Electric Hand Digital Twin can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a central coaxial plug connection, automatic shut-off electronics and integrated on-off switch, low-friction bevel gear, positive back lock and system inner hand. With potentiometer for adjusting the switching threshold (double-channel control). An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand
8E38=7-L7	Left (L)	7	8X18=L7
8E38=7-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E38=7-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E38=7-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E38=7-R7	Right (R)	7	8X18=R7
8E38=7-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E38=7-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E38=7-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4


Technical Data

Reference Number	8E38=7	8E38=7
Size	7	7 1/4, 7 3/4, 8 1/4
Operating voltage	6/7.2 V	6/7.2 V
Opening width	79 mm	100 mm
Approx. maximum gripping strength.	90 N	90 N
Average speed	110 mm/sec	110 mm/sec
Weight with system inner hand	355 g	457 g
for	Women, Adolescents	Women, Adolescents, Men

- The electrodes must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80




 646D44

 647H327

Particularly short and light: the new Electric Hand size 7

It closes the gap between the Electric Hand 2000 hand system for children and the well-known System Electric Hands for adults. It is particularly suitable for the fitting of adolescents or women with small, dainty hands.

The Electric Hand size 7 is offered with the well-known Digital Twin and DMC plus control systems.



 646D44

 647H327

8E39=7 System Electric Hand Digital Twin

With lamination ring

Suitable for wrist disarticulation. Passive wrist rotation with friction.

The Digital Twin control features a digital and dual channel control mode. The desired control mode is selected with the integrated 13E185 function plug. The System Electric Hand Digital Twin can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a centrally guided flat cable with automatic shut-off electronics and integrated on-off switch, low-friction bevel gear, positive back lock and system inner hand. With potentiometer for adjusting the switching threshold (double-channel control). An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand
8E39=7-L7	Left (L)	7	8X18=L7
8E39=7-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E39=7-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E39=7-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E39=7-R7	Right (R)	7	8X18=R7
8E39=7-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E39=7-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E39=7-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4

Technical Data

Reference Number	8E39=7	8E39=7
Size	7	7 1/4, 7 3/4, 8 1/4
Operating voltage	6/7.2 V	6/7.2 V
Opening width	79 mm	100 mm
Approx. maximum gripping strength.	90 N	90 N
Average speed	110 mm/sec	110 mm/sec
Weight with system inner hand	355 g	457 g
for	Women, Adolescents	Women, Adolescents, Men

- The electrodes must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80

Particularly short and light: the new Electric Hand size 7

It closes the gap between the Electric Hand 2000 hand system for children and the well-known System Electric Hands for adults. It is particularly suitable for the fitting of adolescents or women with small, dainty hands.

The Electric Hand size 7 is offered with the well-known Digital Twin and DMC plus control systems.

8E41=7 System Electric Hand Digital Twin

With threaded stud M12x1.5

Suitable for all amputation levels, except wrist disarticulation.

The Digital Twin control features a digital and dual channel control mode. The desired control mode is selected with the integrated 13E185 function plug. The System Electric Hand Digital Twin can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a central cable outlet with automatic battery-saving circuit and integrated on-off switch, low-friction bevel gear, positive back lock and system inner hand. With potentiometer for adjusting the switching threshold (double-channel control). An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand
8E41=7-L7	Left (L)	7	8X18=L7
8E41=7-L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E41=7-L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E41=7-L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E41=7-R7	Right (R)	7	8X18=R7
8E41=7-R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E41=7-R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E41=7-R8 1/4	Right (R)	8 1/4	8X18=R8 1/4



 646D44

 647H327

Technical Data

Reference Number	8E41=7	8E41=7
Size	7	7 1/4, 7 3/4, 8 1/4
Operating voltage	6/7.2 V	6/7.2 V
Opening width	79 mm	100 mm
Approx. maximum gripping strength.	90 N	90 N
Average speed	110 mm/sec	110 mm/sec
Weight with system inner hand	355 g	457 g
for	Women, Adolescents	Women, Adolescents, Men

- The electrodes must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80

Particularly short and light: the new Electric Hand size 7

It closes the gap between the Electric Hand 2000 hand system for children and the well-known System Electric Hands for adults. It is particularly suitable for the fitting of adolescents or women with small, dainty hands.

The Electric Hand size 7 is offered with the well-known Digital Twin and DMC plus control systems.



Transcarpal Hand

For the first time, many patients who in the past could only be fitted with cosmetic prostheses due to their residual limb length can now benefit from the functional advantages of a System Electric Hand thanks to the new Transcarpal Hand.

Short: Many of the components in the Transcarpal Hand have been redesigned to make it especially compact: - 37 mm in comparison with the 8E39 Hands, which up to now had been the shortest ones, and - 49 mm in comparison with the 8E38 Hands. The new anchoring technique, which replaces the normal lamination ring and a newly developed drive unit, are particularly helpful in saving space.

Lightweight: Furthermore, the Transcarpal Hand is approximately 150 g lighter (= - 30%) than comparable hand types for wrist disarticulations.

The Transcarpal Hand is available with **DMC plus**[®] control or **Digital Twin**[®] control variants.

Due to the special anchoring technique used, it is essential that the patient is able to actively pronate and supinate the hand. Passive or electromotive pronation and supination are not possible with the Transcarpal Hand.

The Transcarpal Hand makes active gripping possible for many patients with very low amputation levels.

8E44=6 Transcarpal-Hand DMC plus

With lamination plate

For residual limb lengths from wrist disarticulation to transcarpal levels. Without wrist rotation, which means that active pronation and supination are absolutely required. (Fine adjustment of the basic position is possible, however, only after the socket has been completed.)

The DMC plus control features a DMC and DMC plus control mode. In DMC plus control mode, after gripping once with maximum grip force, a higher signal is required to open the hand. This reduces the risk of opening the hand with undesired muscle signals. The desired control mode is selected with the integrated 13E185 Function Plug. In the DMC plus system, two independent measurement and control systems proportionally control gripping speed as well as gripping force. The Transcarpal-Hand DMC plus can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a centrally guided flat cable, automatic shut-off electronics and integrated on-off switch, low-friction bevel gear, positive back lock and system inner hand. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand	For
8E44=6-L7 1/4	Left (L)	7 1/4	8X24=L7 1/4 mm	Women, Adolescent
8E44=6-L7 3/4	Left (L)	7 3/4	8X24=L7 3/4	Men
8E44=6-L8 1/4	Left (L)	8 1/4	8X24=L8 1/4	Men
8E44=6-R7 1/4	Right (R)	7 1/4	8X24=R7 1/4	Women, Adolescent
8E44=6-R7 3/4	Right (R)	7 3/4	8X24=R7 3/4	Men
8E44=6-R8 1/4	Right (R)	8 1/4	8X24=R8 1/4	Men

Technical Data

Reference Number	8E44=6	8E44=6	8E44=6
Size	7 1/4	7 3/4	8 1/4
Operating voltage	6/7,2 V	6/7,2 V	6/7,2 V
Operating temperature	0 bis +70 °C	0 bis +70 °C	0 bis +70 °C
Opening width	100 mm	100 mm	100 mm
Proportional gripping force	0 - 90 N	0 - 90 N	0 - 90 N
Proportional speed	15 - 130 mm/sec	15 - 130 mm/sec	15 - 130 mm/sec
Weight with system inner hand	308 g	308 g	308 g

- The electrodes must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80



646D44

647H398

1

2

3

4

5

6

7

8



 646D44

 647H398

8E44=7 Transcarpal-Hand Digital Twin

With lamination plate

For residual limb lengths from wrist disarticulation to transcarpal levels. Without wrist rotation, which means that active pronation and supination are absolutely required. (Fine adjustment of the basic position is possible, however, only after the socket has been completed.)

The Digital Twin control features a digital and dual channel control mode. The desired control mode is selected with the integrated 13E185 function plug. The Transcarpal-Hand Digital Twin can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a centrally guided flat cable, automatic shut-off electronics and integrated on-off switch, low-friction bevel gear, positive back lock and system inner hand. With potentiometer for adjusting the switching threshold. An integrated slip clutch allows the hand to be opened in case of power supply or myoelectric control failure.

Article number	Side	Size	Inner hand	For
8E44=7-L7 1/4	Left (L)	7 1/4	8X24=L7 1/4	Women, Adolescent
8E44=7-L7 3/4	Left (L)	7 3/4	8X24=L7 3/4	Men
8E44=7-L8 1/4	Left (L)	8 1/4	8X24=L8 1/4	Men
8E44=7-R7 1/4	Right (R)	7 1/4	8X24=R7 1/4	Women, Adolescent
8E44=7-R7 3/4	Right (R)	7 3/4	8X24=R7 3/4	Men
8E44=7-R8 1/4	Right (R)	8 1/4	8X24=R8 1/4	Men

Technical Data

Reference Number	8E44=7	8E44=7	8E44=7
Size	7 1/4	7 3/4	8 1/4
Operating voltage	6/7,2 V	6/7,2 V	6/7,2 V
Operating temperature	0 bis +70 °C	0 bis +70 °C	0 bis +70 °C
Opening width	100 mm	100 mm	100 mm
Approx. maximum gripping strength.	90 N	90 N	90 N
Average speed	110 mm/sec	110 mm/sec	110 mm/sec
Weight with system inner hand	308 g	308 g	308 g

- The electrodes must be adjusted with the 757M11 MyoBoy!
- For the suitable system prosthetic glove, see pages 77-80

8E12 System Electric Hand for other systems

For switch control or other manufacturer's systems, 6 Volt, with M12x1.5 threaded stud


Suitable for all amputation levels, except wrist disarticulation.

It features a central cable outlet and on-off switch, positive back lock and System Inner Hand. Without electronic control.

Article number	Side	Size	Inner hand
8E12=L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E12=L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E12=L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E12=R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E12=R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E12=R8 1/4	Right (R)	8 1/4	8X18=R8 1/4

Technical Data

Reference Number	8E12	8E12	8E12
Size	7 1/4	7 3/4	8 1/4
Operating voltage	6 V	6 V	6 V
Approx. average current consumption	200 mA	200 mA	200 mA
Approx. max current	800 mA	800 mA	800 mA
Opening width	100 mm	100 mm	100 mm
Approx. maximum gripping strength.	90 N	90 N	90 N
Average speed	110 mm/sec	110 mm/sec	110 mm/sec
Weight with system inner hand	460 g	460 g	460 g
for	Women, Adolescents	Men	Men



647G133

- ⚠ **Not suitable for operation with EnergyPack!**
- ⚠ For the suitable system prosthetic glove, see pages 77-80

1

2

3

4


5

6

7

8



 647G133

8E37 System Electric Hand for other systems

For other manufacturer's systems, 6 Volt, with Quick Disconnect Wrist

Suitable for all amputation levels, except wrist disarticulation.

With central coaxial plug connection, on-off switch, positive back lock and System Inner Hand. Without electronic control.

Article number	Side	Size	Inner hand
8E37=L7 1/4	Left (L)	7 1/4	8X18=L7 1/4
8E37=L7 3/4	Left (L)	7 3/4	8X18=L7 3/4
8E37=L8 1/4	Left (L)	8 1/4	8X18=L8 1/4
8E37=R7 1/4	Right (R)	7 1/4	8X18=R7 1/4
8E37=R7 3/4	Right (R)	7 3/4	8X18=R7 3/4
8E37=R8 1/4	Right (R)	8 1/4	8X18=R8 1/4

Technical Data

Reference Number	8E37	8E37	8E37
Size	7 1/4	7 3/4	8 1/4
Operating voltage	6 V	6 V	6 V
Approx. average current consumption	200 mA	200 mA	200 mA
Approx. max current	800 mA	800 mA	800 mA
Opening width	100 mm	100 mm	100 mm
Approx. maximum gripping strength.	90 N	90 N	90 N
Average speed	110 mm/sec	110 mm/sec	110 mm/sec
Weight with system inner hand	420 g	420 g	420 g
for	Women, Adolescents	Men	Men

- ⦿ **Not suitable for operation with EnergyPack!**
- ⦿ Not suitable for Ottobock switch control.
- ⦿ Specially adapted for other manufacturer's systems, such as INAIL, Utah.
- ⦿ For the suitable system prosthetic glove, see pages 77-80

Spare parts for System Electric Hands

8X18 System Inner Hand

For Ottobock System Electric Hands in sizes 6 3/4, 7, 7 1/4, 7 3/4 and 8 1/4. Energy saving, lightweight plastic design, partially reinforced, with wire inserts in the fingers and sealed retainer ring.

Article number	Side	Size	Retaining ring
8X18=L6 3/4	Left (L)	6 3/4	9S15=42
8X18=L7	Left (L)	7	9S15=48
8X18=L7 1/4	Left (L)	7 1/4	9S187=7 1/4
8X18=L7 3/4	Left (L)	7 3/4	9S187=7 3/4
8X18=L8 1/4	Left (L)	8 1/4	9S187=8 1/4
8X18=R6 3/4	Right (R)	6 3/4	9S15=42
8X18=R7	Right (R)	7	9S15=48
8X18=R7 1/4	Right (R)	7 1/4	9S187=7 1/4
8X18=R7 3/4	Right (R)	7 3/4	9S187=7 3/4
8X18=R8 1/4	Right (R)	8 1/4	9S187=8 1/4



• For the suitable system prosthetic glove, see pages 77-80

9S187 Retaining ring, wide

Article number	For hand size
9S187=7 1/4	7 1/4
9S187=7 3/4	7 3/4
9S187=8 1/4	8 1/4

Technical Data

Article number	9S187=7 1/4	9S187=7 3/4	9S187=8 1/4
for	8X18 System Inner Hands	8X18 System Inner Hands	8X18 System Inner Hands



1



8X24 System Inner Hand

Article number	Side	Size
8X24=L7 1/4	Left (L)	7 1/4
8X24=L7 3/4	Left (L)	7 3/4
8X24=L8 1/4	Left (L)	8 1/4
8X24=R7 1/4	Right (R)	7 1/4
8X24=R7 3/4	Right (R)	7 3/4
8X24=R8 1/4	Right (R)	8 1/4

2

Technical Data

Article number	8X24=L7 1/4	8X24=L7 3/4	8X24=L8 1/4	8X24=R7 1/4	8X24=R7 3/4	8X24=R8 1/4
Size	7 1/4	7 3/4	8 1/4	7 1/4	7 3/4	8 1/4
for	Transcarpal Hand					

3

• For the suitable system prosthetic glove, see pages 77-80

4



624Z12 Special Thread

For sealing the 8X24 Inner Hand from the Transcarpal Hand

5

Accessories for System Electric Hands

6



8Y1 Tweezers

7



743F1 Pinch Gauge

For testing the grip force of System Electric Hands

8

453A1 Derma Protection ArmComfort

Improves the transition from the below-elbow socket to the upper arm. The special polymer gel layer, the high elasticity and the anatomic shape increase wearing comfort.

Article number	453A1	453A1=S
for	Adults	Children



1

2

3

4

5

6

7

8



Electric Greifer System

The Otto Bock Electric Greifer System increases the number of potential applications of the MYOBOCK System for work and special tasks. The Electric Greifer has a flexion joint as well as a hand joint with passive pronation and supination.

Gripping surfaces may be individually adjusted by adjusting the finger tips.

The Electric Greifer is especially useful for heavy and manual labour, and where precise gripping is required.

Using a quick disconnect wrist, the Greifer can be quickly replaced with a System Electric Hand

- > **Large opening width**
- > **Digital Twin Control**
- > **DMC VariPlus Control**
- > **Adjustable gripping tips**
- > **Safety lever for releasing the hand**
- > **Hand wheel allows manual operation**
- > **Energy-saving**

8E33=9 System Electric Greifer DMC VariPlus

With quick-disconnect wrist

Suitable for all amputation levels, except wrist disarticulation.

The System Electric Greifer DMC VariPlus is based on the Ottobock DMC system (DMC=Dynamic Mode Control): This system, which has been developed by Ottobock, uses two independent measurement and control systems in order to optimally control the grip speed and grip force in accordance with the patient's muscle signal.

6 different programs can be selected and adjusted using the 757T13 MyoSelect. They allow to optimally adapt the prosthesis to the amputee's individual abilities and requirements.

The System Electric Greifer DMC VariPlus can be operated with the 757B35=* MyoEnergy Integral, the 757B20 or 757B21 EnergyPack or the 757B15 X-ChangePack.

With central coaxial plug connection, energy-saving, automatic shut-off electronics and on-off switch. Metal finger tips and grip surfaces covered with rubber for normal gripping activities.

Technical Data

Article number	8E33=9
Weight	540 g
Operating voltage	6/7.2 V
Opening width	95 mm
Proportional gripping force	0-160 N
Proportional speed	8-200 mm/sec

 647H382

- The electrodes must be adjusted with the 757M11 MyoBoy!

8E34=9 System Electric Greifer DMC VariPlus

With lamination ring

Suitable for wrist disarticulation.

The System Electric Greifer DMC VariPlus is based on the Ottobock DMC system (DMC=Dynamic Mode Control): This system, which has been developed by Ottobock, uses two independent measurement and control systems in order to optimally control the grip speed and grip force in accordance with the patient's muscle signal.

6 different programs can be selected and adjusted using the 757T13 MyoSelect. They allow to optimally adapt the prosthesis to the amputee's individual abilities and requirements.

The System Electric Greifer DMC VariPlus can be operated with the 757B35=* MyoEnergy Integral, the 757B20 or 757B21 EnergyPack or the 757B15 X-ChangePack.

With central flat cable connection, energy-saving, automatic shut-off electronics and on-off switch. Metal finger tips and grip surfaces covered with rubber for normal gripping activities.

Technical Data

Article number	8E34=9
Weight	520 g
Operating voltage	6/7.2 V
Opening width	95 mm
Proportional gripping force	0-160 N
Proportional speed	8-200 mm/sec

 647H382

- The electrodes must be adjusted with the 757M11 MyoBoy!



1

2

3

4

5

6

7

8

Control Modes for the System Electric Greifer DMC VariPlus

Control mode 1: DMC plus®

Control with 2 electrodes

Control mode 1	OPEN	CLOSE	Indication
DMC plus®	Sustained electrode signal Speed: proportional 8 to 200 mm/s	Sustained electrode signal Grip force: proportional After gripping once with maximum force, the electrode signal required to open the Greifer will be set to a higher value. Opening the Greifer with unwanted electrode signals is prevented. Speed: proportional 8 to 200 mm/s	For patients with two strong electrode signals

Control mode 2: AutoControl – Low Input

Control with 2 electrodes

Control mode 2	OPEN	CLOSE	Indication
AutoControl-Low Input	Sustained electrode signal Reduced range of proportionality: Maximum speed once the LOW-threshold has been reached Speed: proportional 8 to 200 mm/s	Sustained electrode signal Grip force: Time proportional After gripping once with maximum force, the electrode signal required to open the Greifer will be set to a higher value. Opening the Greifer with unwanted electrode signals is prevented. Speed: constant 200 mm/s	For amputees with two weak electrode signals

Control with 1 electrode and 1 switch

Control mode 2	OPEN	CLOSE	Indication
AutoControl-Low Input	Sustained electrode signal Reduced range of proportionality: Maximum speed once the LOW-threshold has been reached Speed: proportional 8 to 200 mm/s	Signal from the switch Grip force: Time proportional After gripping once with maximum force, the electrode signal required to open the Greifer will be set to a higher value. Opening the Greifer with unwanted electrode signals is prevented. Speed: constant 200 mm/s	For amputees with only one muscle and weak electrode signal

Control with 1 switch

Control mode 2	OPEN	CLOSE	Indication
AutoControl-Low Input	Greifer opens as long as the OPEN-side of the switch is operated. Speed: constant 200 mm/s	Greifer closes as long as the CLOSE-side of the switch is operated. Grip force: Time proportional Speed: constant 200 mm/s	For amputees with weak or no electrode signal

Control mode 3: VarioControl

Control with 1 electrode

Control mode 3	OPEN	CLOSE	Indication
VarioControl	Increasing electrode signal through muscle contraction Speed and strength of muscle contraction to the electrode Speed: proportional 8 to 200 mm/s	Declining electrode signal through muscle relaxation Grip force: Proportional to the decline of the electrode signal After gripping once with maximum force, the electrode signal required to open the Greifer will be set to a higher value. Opening the Greifer with unwanted electrode signals is prevented. Speed: proportional 8 to 200 mm/s	For amputees with one strong electrode signal or tendency to co-contraction

Control with 1 linear transducer

Control mode 3	OPEN	CLOSE	Indication
VarioControl	Speed and strength of pull on the linear transducer Speed: proportional 8 to 200 mm/s	Speed of release of pull on the linear transducer Grip force: Proportional to the release of the pull on the linear transducer After gripping once with maximum force, the electrode signal required to open the Greifer will be set to a higher value. Opening the Greifer with unwanted electrode signals is prevented. Speed: proportional 8 to 200 mm/s	For amputees with weak or no electrode signal

Control mode 4: VarioDual

Control with 2 electrodes

Control mode 4	OPEN	CLOSE	Indication
VarioDual	Increasing electrode signal through muscle contraction to the first electrode Speed and strength of muscle tension to the electrode Speed: Proportional to the increase of the electrode signal 8 to 200 mm/s	Declining electrode signal through muscle relaxation to the first electrode or: Sustained electrode signal to the second electrode Grip force: Proportional to the signal strength to the second electrode. If only the first electrode is used, the Greifer closes up to the minimal grip force. Grip force increase: By another, stronger electrode signal to the second electrode. After gripping once with maximum force, the electrode signal required to open the Greifer will be set to a higher value. Opening the Greifer with unwanted electrode signals is prevented. Speed: Proportional to the decline of the electrode signal to the first electrode / proportional to the signal strength to the second electrode. 8 to 200 mm/s	For amputees with two strong electrode signals

Control mode 5: Digital Control

Control with 2 electrodes

Control mode 5	OPEN	CLOSE	Indication
Digital Control	Sustained electrode signal Speed: constant 200 mm/s	Sustained electrode signal Grip force: Duration of the signal Speed: constant 200 mm/s	For amputees with two weak electrode signals

Control with 1 electrode and 1 switch

Control mode 5	OPEN	CLOSE	Indication
Digital Control	Sustained electrode signal Speed: constant 200 mm/s	Signal from the switch Grip force: Duration of the signal Speed: constant 200 mm/s	For amputees with only one muscle and weak electrode signal

Control with 1 switch

Control mode 5	OPEN	CLOSE	Indication
Digital Control	Greifer opens as long as the OPEN-side of the switch is operated. Speed: constant 200 mm/s	Greifer closes as long as the CLOSE-side of the switch is operated. Grip force: Duration of the signal Speed: constant 200 mm/s	For amputees with weak or no electrode signals

Control mode 6: Double Channel Control

Control with 1 electrode

Control mode 6	OPEN	CLOSE	Indication
Double Channel Control	Quick, strong electrode signal that reaches the upper threshold within 80 ms and is maintained above the upper threshold for at least 30 ms Speed: constant 200 mm/s	Slow, gentle electrode signal that does not reach the upper threshold within 80 ms Grip force: Duration of the signal Speed: constant 200 mm/s	For amputees with one strong electrode signal

1

2

3

4

5

6

7

8



 647H382

8E33=7 System Electric Greifer Digital Twin

[With quick-disconnect wrist](#)

Suitable for all amputation levels, except wrist disarticulation.

The Digital Twin control features a digital and dual channel control mode. The desired control mode is selected with the integrated 13E189 function plug. The Electric Greifer Digital Twin can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack. Passive wrist rotation with ratchet lock (can be replaced by the 11S30 Friction Ring).

With central coaxial plug connection, automatic shut-off electronics and on-off switch. Metal finger tips and grip surfaces covered with rubber for normal gripping activities. With two potentiometers to adjust the switching threshold.

Technical Data

Article number	8E33=7
Weight	540 g
Operating voltage	6/7.2 V
Opening width	95 mm
Approx gripping strength	160 N
Approx. gripping speed	180 mm/sec

- ⦿ The electrodes must be adjusted with the 757M11 MyoBoy!



 647H382

8E34=7 System Electric Greifer Digital Twin

[With lamination ring](#)

Suitable for wrist disarticulation. Passive wrist rotation with friction.

The Digital Twin control features a digital and dual channel control mode. The desired control mode is selected with the integrated 13E189 function plug. The Electric Greifer Digital Twin can be operated with the 757B35=* MyoEnergy Integral, the 757B20/757B21 EnergyPack or the 757B15 X-ChangePack.

It features a central flat cable outlet with automatic shut-off electronics and on-off switch. Metal finger tips and grip surfaces covered with rubber for normal gripping activities. With two potentiometers to adjust the switching threshold.

Technical Data

Article number	8E34=7
Weight	520 g
Operating voltage	6/7.2 V
Opening width	95 mm
Approx gripping strength	160 N
Approx. gripping speed	180 mm/sec

- ⦿ The electrodes must be adjusted with the 757M11 MyoBoy!

8E32=6 Electric Greifers for other manufacturer's systems

For other manufacturer's systems, 6 Volt, with Quick Disconnect Wrist

Suitable for all amputation levels, except wrist disarticulation.

Passive wrist rotation with ratchet lock (can be replaced with the 11S30 Friction Ring).


With central coaxial plug connection and on-off switch. Without electronic control. Metal finger tips and grip surfaces covered with rubber for normal gripping activities.

Technical Data

Article number	8E32=6
Weight	540 g
Operating voltage	6 V
Approx. average current consumption	200 mA
Approx. max current	700 mA
Opening width	95 mm
Approx gripping strength	140 N
Average speed	120 mm/sec

- Not suitable for operation with EnergyPack.
- Not suitable for Ottobock switch control.
- Specially adapted for other manufacturer's systems, such as INAIL, Utah.



 647G50

1

2

3

4

5

6

7

8

Accessories for System Electric Greifer

1



9S138 Finger Tip Set

For use without tip padding,

Article number	9S138
Consists of:	Finger Tips (pair) Finger Tip Blanks (pair) 501S54=M3,3x8 Oval Head Countersunk Screw (2 pieces)

2

3



9S145 Rubber Gripping Pad Set

Article number	9S145
Consists of:	9S146 Tip Padding (pair) 9S147=1PAA Rubber pad (pair) 9S147=2PAA Rubber pad (pair) 9S148 Padded Stirrup

4

5



9S149 Finger Tip Blank Set

For subsequent adaptation to special working needs,

Article number	9S149
Consists of:	Spacer (2 pieces) Finger Tip Blanks (pair) 501S54=M3,3x8 Oval Head Countersunk Screw (2 pieces)

6

7



9S234=PAA Finger Tip (pair)

For use with tip padding,

Article number	9S234=PAA
with	501S54=M3,3x8 Oval Head Countersunk Screw (2 pieces)

8

9S278=PAA Finger Tip (pair)

For use without tip padding,

Article number	9S278=PAA
with	501S54=M3,3x8 Oval Head Countersunk Screw (2 pieces)



1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



Myo Prosthetic Gloves

1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



 646D423  647G571

8S20N MyolinoSkin Natural

In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.

To order them, please add an **N** before the =:

Example 8S20=136x41xL4 Standard glove

8S20N=136x41xL4 MyolinoSkin Natural

For available colours, please ask our customer service.

The PVC-based material ensures that Skin Natural is very robust compared to other materials and it has a relatively long life. This provides security when using the gloves. Due to a special modern surface treatment, the Skin Natural gloves are also easy to clean.

Article number	Side	Size	Sleeve length	Sleeve end circumference
8S20N=136x41L	Left (L)	5	210 mm	200 mm
8S20N=136x41R	Right (R)	5	210 mm	200 mm
8S20N=147x45L	Left (L)	5 1/2	215 mm	200 mm
8S20N=147x45R	Right (R)	5 1/2	215 mm	200 mm
8S20N=162x56L	Left (L)	6	220 mm	210 mm
8S20N=162x56R	Right (R)	6	220 mm	210 mm
8S20N=177x64L	Left (L)	6 1/2	240 mm	220 mm
8S20N=177x64R	Right (R)	6 1/2	240 mm	220 mm

The entire colour palette of the Skin Natural is produced with six colours. The new 646M47 colour scale helps when selecting the right pattern.

- For special cleaner for prosthetic gloves, see pages 81, 341

8S11N MyoSkin Natural

In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.

To order them, please add an **N** before the =:

Example 8S11=210x78xL4 Standard glove

8S11N=210x78xL4 MyoSkin Natural

For available colours, please ask our customer service.

The PVC-based material ensures that Skin Natural is very robust compared to other materials and it has a relatively long life. This provides security when using the gloves. Due to a special modern surface treatment, the Skin Natural gloves are also easy to clean.

Article number	Side	Size	Inner hand	Sleeve length	Sleeve end circumference	Package
8S11N=190x76L	Left (L)	7 1/4	8X18=L7 1/4 8X24=L7 1/4	300 mm	230 mm	Adolescents / Men
8S11N=190x76R	Right (R)	7 1/4	8X18=R7 1/4 8X24=R7 1/4	300 mm	230 mm	Adolescents / Men
8S11N=210x78L	Left (L)	7 3/4	8X18=L7 3/4 8X24=L7 3/4	320 mm	250 mm	Men
8S11N=210x78R	Right (R)	7 3/4	8X18=R7 3/4 8X24=R7 3/4	320 mm	250 mm	Men
8S11N=225x80L	Left (L)	8 1/4	8X18=L8 1/4 8X24=L8 1/4	350 mm	260 mm	Men
8S11N=225x80R	Right (R)	8 1/4	8X18=R8 1/4 8X24=R8 1/4	350 mm	260 mm	Men

The entire colour palette of the Skin Natural is produced with six colours. The new 646M47 colour scale helps when selecting the right pattern.



646D423



647G571

- For special cleaner for prosthetic gloves, see pages 81, 341

1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



646D423 647G571

8S12N MyoSkin Natural

In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.

To order them, please add an **N** before the =:

Example 8S12=190x78xL4 Standard glove

8S12N=190x78xL4 MyoSkin Natural

For available colours, please ask our customer service.

The PVC-based material ensures that Skin Natural is very robust compared to other materials and it has a relatively long life. This provides security when using the gloves. Due to a special modern surface treatment, the Skin Natural gloves are also easy to clean.

Article number	Side	Size	Inner hand	Sleeve length	Sleeve end circumference	Package
8S12N=190x78L	Left (L)	7 1/4	8X18=L7 1/4 8X24=L7 1/4	460 mm	250 mm	Women
8S12N=190x78R	Right (R)	7 1/4	8X18=R7 1/4 8X24=R7 1/4	460 mm	250 mm	Women

The entire colour palette of the Skin Natural is produced with six colours. The new 646M47 colour scale helps when selecting the right pattern.

For special cleaner for prosthetic gloves, see pages 81, 341



646D423 647G571

8S13N MyoSkin Natural

In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.

To order them, please add an **N** before the =:

Example 8S13=7L4 Standard glove

8S13N=7L4 MyoSkin Natural

For available colours, please ask our customer service.

The PVC-based material ensures that Skin Natural is very robust compared to other materials and it has a relatively long life. This provides security when using the gloves. Due to a special modern surface treatment, the Skin Natural gloves are also easy to clean.

Article number	Side	Size	Inner hand	Sleeve length	Sleeve end circumference	Package
8S13N=7L	Left (L)	7	8X18=L8 1/4	280 mm	242 mm	Adolescents / Women
8S13=7R	Right (R)	7	8X18=R7	280 mm	242 mm	Adolescents / Women

The entire colour palette of the Skin Natural is produced with six colours. The new 646M47 colour scale helps when selecting the right pattern.

For special cleaner for prosthetic gloves, see pages 81, 341

8S20 Prosthetic Glove for Children

The prosthetic glove for children impresses with its natural appearance, durability and flexibility.

Article number	Side	Size	Sleeve length	Sleeve end circumference
8S20=136X41L*	Left (L)	5	210 mm	200 mm
8S20=147X45L*	Left (L)	5 1/2	215 mm	200 mm
8S20=162X56L*	Left (L)	6	220 mm	210 mm
8S20=177x64L*	Left (L)	6 1/2	240 mm	220 mm
8S20=136X41R*	Right (R)	5	210 mm	200 mm
8S20=147X45R*	Right (R)	5 1/2	215 mm	200 mm
8S20=162X56R*	Right (R)	6	220 mm	210 mm
8S20=177X64R*	Right (R)	6 1/2	240 mm	220 mm

*Available in 18 different shades. When ordering, please add the colour code according to the 646M3 colour swatches to the end of the order number, e.g., **8S20=136x41L4**

- For special cleaner for prosthetic gloves, see pages 81, 341



646D49

647G468

8S11 Prosthetic Glove for Adolescents and Men

The Ottobock prosthetic glove for adolescents and men impresses with its natural appearance, durability and flexibility.

Article number	Side	Size	Inner hand	Sleeve length	Sleeve end circumference	Package
8S11=190X76L	Left (L)	7 1/4	8X18=L7 1/4 8X24=L7 1/4	300 mm	230 mm	Adolescents / Men
8S11=190X76R	Right (R)	7 1/4	8X18=R7 1/4 8X24=R7 1/4	300 mm	230 mm	Adolescents / Men
8S11=192X78L	Left (L)	7 1/4	8X18=L7 1/4 8X24=L7 1/4	440 mm	245 mm	Adolescents / Men
8S11=192X78R	Right (R)	7 1/4	8X24=R7 1/4	440 mm	245 mm	Adolescents / Men
8S11=210X78L	Left (L)	7 3/4	8X18=L7 3/4 8X24=R7 3/4	320 mm	250 mm	Men
8S11=210X78R	Right (R)	7 3/4	8X18=R7 3/4 8X24=R7 3/4	320 mm	250 mm	Men
8S4=220X80L	Left (L)	8	8X14=L8	225 mm	210 mm	Men
8S4=220X80R	Right (R)	8	8X14=R8	225 mm	210 mm	Men
8S11=225X80L	Left (L)	8 1/4	8X18=L8 1/4 8X24=L8 1/4	350 mm	260 mm	Men
8S11=225X80R	Right (R)	8 1/4	8X24=R8 1/4	350 mm	260 mm	Men

*Available in 18 different shades. When ordering, please add the colour code according to the 646M3 colour swatches to the end of the order number, e.g., **8S11=190x76L4**

- For special cleaner for prosthetic gloves, see pages 81, 341



646D49

647G468

1



8S12 Prosthetic Glove for Women

The prosthetic glove for women impresses with its natural appearance, durability and flexibility.

Article number	Side	Size	Inner hand	Sleeve length	Sleeve end circumference	Package
8S12=190X78L	Left (L)	7 1/4	8X18=L7 1/4 8X24=L7 1/4	460 mm	250 mm	Women
8S12=190X78R	Right (R)	7 1/4	8X18=R7 1/4 8X24=R7 1/4	460 mm	250 mm	Women

*Available in 18 different shades. When ordering, please add the colour code according to the 646M3 colour swatches to the end of the order number, e.g., **8S12=190x78L4**

- For special cleaner for prosthetic gloves, see pages 81, 341

2

3

4

646D49 647G468

5



8S13 Prosthetic Glove for Children and Adolescents

The prosthetic glove for children and adolescents impresses with its natural appearance, durability and flexibility.

Article number	Side	Size	Inner hand	Sleeve length	Sleeve end circumference	Package
8S13=7L	Left (L)	7	8X18=L7	280 mm	242 mm	Adolescents / Women
8S13=7R	Right (R)	7	8X18=R7	280 mm	242 mm	Adolescents / Women

*Available in 18 different shades. When ordering, please add the colour code according to the 646M3 colour swatches to the end of the order number, e.g., **8S13=7L4**

- For special cleaner for prosthetic gloves, see pages 81, 341

6

7

8

646D49 647G468

Accessories for Myo Prosthetic Gloves

633S2 Procomfort Gel

Acts as a lubricant to aid in donning the prosthetic glove over the inner hand.

Article number	Net contents
633S2	250 ml



640F12 Special Cleaner

Technical Data

Article number	640F12
for	prosthetic gloves



635P15 Colouring Pen

Refillable, for colouring prosthetic gloves

Article number	Colour
635P15=2	red
635P15=5	blue



1

2

3

4

5

6

7

8

1

2

3


4

5

6

7

8



Myo Flexion

1

2

3

4

5

6

7

8

Wrist Units

1
2
3
4
5
6
7
8



647G460

10V51 MyolinoWrist 2000

The 2000 10V51=* MyolinoWrist is a mechanical ball wrist unit with adjustable friction. When combined with the Electric Hand 2000 8E51, the hand can be positioned in any direction, helping to significantly eliminate unnatural compensation with the arm. This additional freedom provides more natural movements, whilst at the same time, the wrist unit promotes physiologically correct posture.

Article number	Wrist connection Ø
10V51=1	40 mm
10V51=2	40 mm

Technical Data

Article number	10V51=1	10V51=2
Weight	47 g	50 g
Overall length	32 mm	40 mm



647G351

10V38 MyoWrist Transcarpal

Enables flexion and extension of a myoelectrically controlled 8E44 Transcarpal Hand. Can be locked in 5 positions from -40° to $+40^\circ$. Passive pronation and supination is possible thanks to the quick-disconnect wrist. In combination with the optional 10S17 Electric Wrist Rotator, electromotive pronation and supination can be offered.

Article number	Side
10V38=L7 1/4	Left (L)
10V38=L7 3/4	Left (L)
10V38=L8 1/4	Left (L)
10V38=R7 1/4	Right (R)
10V38=R7 3/4	Right (R)
10V38=R8 1/4	Right (R)

Technical Data

Reference Number	10V38
Outer Ø	36 mm
Weight	85 – 89 g
Overall length	24.7 mm
for	8E44=6-L7 1/4 bzw. 8E44=7-L7 1/4, 8E44=6-L7 3/4 bzw. 8E44=7-L7 3/4, 8E44=6-L8 1/4 bzw. 8E44=7-L8 1/4, 8E44=6-R7 1/4 bzw. 8E44=7-R7 1/4, 8E44=6-R7 3/4 bzw. 8E44=7-R7 3/4, 8E44=6-R8 1/4 bzw. 8E44=7-R8 1/4
Consists of:	10V38 MyoWrist Transcarpal 9E388 Cable Seal 9S267=* Lock Ring

10V40 MyoWrist 2Act

The MyoWrist 2Act 10V40 was developed especially for all variants of the Otto Bock System Electric Hands with the 8E41 Threaded Stud to maintain a low structural height for the whole system.


Flexion and extension possible with 5 locking positions from -40° to $+40^{\circ}$ in 20° increments. Rotation is performed passively via the wrist lock. Active rotation is also possible in combination with the MyoRotronic 13E205.

Article number	Outside Ø
10V40	36 mm

Technical Data

Article number	10V40
Weight	55 g
Overall length	26 mm



 647G459

1

2

3

4

5

6

7

8

Accessories for 10A38 and 10A40

1



9S266 Chassis with Quick-Disconnect Wrist

Article number	Complete with
9S266	Handgelenkverschluss

2


3



9E397 Coaxial Bushing

Article number	Control unit	Hand connection
9E397=7	DMC plus control	8E41 + 8E44
9E397=8	Digital Twin control	8E41 + 8E44
9E397=9	SensorHand Speed control	8E41
9E397=10	VariPlus Speed® control	8E41

4



11D61 Laminationring Set

Technical Data

Article number	11D61
for	Myolino Wrist 2000

5


6



709S42 Allen Wrench

7

8



1
2
3
4
5
6
7
8

1

2

3

4

5

6

7

8



Myo Rotation

1

2

3

4

5

6

7

8

1



10S16 Lamination Ring

The 10S16=* lamination ring connects the 2000 8E51 Electric Hand to a custom-made forearm socket. The low installation height allows fittings with relatively long below-elbow residual limbs, since the lamination ring is directly laminated into the outer shell. The O-Rings included in the delivery make it possible to ideally adjust the desired friction according to the prosthesis wearer's needs, allowing complete rotation of the of the 2000 8E51 Electric Hand

2

Article number	For hand size	Ø	with
10S16=34	5, 5 1/2	34 mm	627F3 O-Ring
10S16=38	6, 6 1/2	38 mm	627F3 O-Ring

3



10S1 Lamination Ring

Article number	For hand size
10S1=40	7
10S1=45	7, 7 1/4
10S1=50	7 3/4, System-Elektrogreifer 8E33=*
10S1=54	8 1/4

4

Technical Data

Reference Number	10S1
for	8E38 System Electric Hands or 8E33 System Electric Greifers
Complete with	lamination protection cover for bonding

5



9S110 Lamination Ring

For laminating all 8E39 System Electric Hands or 8E34 System Electric Greifers

Article number	For hand size
9S110=50	7, 7 1/4, 7 3/4, System-Elektrogreifer 8E34=*
9S110=54	8 1/4

6

7



9S258 Lamination Plate

Article number	For hand size	Side
9S258=L7 1/4	7 1/4	Left (L)
9S258=L7 3/4	7 3/4	Left (L)
9S258=L8 1/4	8 1/4	Left (L)
9S258=R7 1/4	7 1/4	Right (R)
9S258=R7 3/4	7 3/4	Right (R)
9S258=R8 1/4	8 1/4	Right (R)

8

Technical Data

Reference Number	9S258
for	8E44 Ottobock Transcarpal Hand

9E169 Coaxial Plug

For connecting the two electrodes and the battery, consisting of:

Article number	9E169
Consists of:	Sicherungsring 9E170 501S50=M4x6 Coaxial plug piece



1

2

10S4 Coupling Piece

Article number	10S4
Complete with	11S4 Lock Ring




3

10S17 Electric Wrist Rotator

6/7.2 Volt

Article number	10S17
Weight	96 g
Operating voltage	6/7.2 V
Approx. no-load current	150 mA
Approx. stall current	1000 mA
Rotation range	360 30°
Idle speed	13.5 U/min
Corresponds to a rotation angle of	81 28°/sec
for	electrical rotation of the System Electric Hand or Greifer,
Consists of:	11S4 Lock Ring 11S61 Wrist Drive with Coaxial Plug 11S25 Lock Ring 9E85 Protective Cap 9E363 Drive Unit 9E365 Protection Plug



 647H204

4

5

6

7

8



647G361

13E205 MyoRotronic

Enables proportional or digital control of the 10S17 Electric Wrist Rotator. Compatible with all MYOBOCK System Electric Hands and System Electric Greifers with Quick-Disconnect Wrist.

Not suitable for use with System Electric Hands for other manufacturer's systems.

5 different programmes can be selected using the 757T13 Myo Select and allow individual adaptation to the respective patient.

One or two 13E200 Electrodes or 13E202 Suction Socket Electrodes or a combination of one electrode and a 9X50/9X52 Linear Control Element are used to control the 13E205 MyoRotronic. Delivery includes a 13E183 Buzzer for acoustic feedback on the switchover.

Technical Data

Article number	13E205
Static current	1 mA
Operating temperature	0 bis +70 °C
Power supply	MyoEnergy Integral 757B35 (7,4 V) EnergyPack 757B20 / 757B21 (7,2 V) X-ChangePack 757B15 (6 V) Wechsel-Akkumulator 757B8 (6 V)
Power off	load dependent between 30 ms and 10 s

- The 757T13 MyoSelect is required to select programmes and adjust the 13E205 MyoRotronic.
- When using the 10S17 Electric Wrist Rotator and the 13E205 MyoRotronic, the distance between the residual limb end and wrist must be **at least 65 mm**.
- The rotation can be switched off using the 9X24 Connection Cable.

8R1 Endoskeletal Adapter

Can be used to fit both long and short residual limbs in combination with the 8E44 Transcarpal Hand (not included in the delivery!).

The Endoskeletal Adapter allows the use of the Transcarpal Hand as a particularly short and lightweight prosthetic hand with quick-disconnect mechanism for fitting of both long and short residual limbs.


Particularly lightweight: compared with a conventional 8E38 System Electric Hand, the weight is reduced by approx. 80 – 100 g (approx. -20 %), depending on size and tube length, without compromising function. In addition, the proximal position of the quick-disconnect mechanism is beneficial to the prosthesis wearer.

Particularly short: the shortest version of the Transcarpal Hand with the quick-disconnect mechanism of the Endoskeletal Adapter is as short as 145 mm, whereas the comparable System Electric Hand is 170 mm (= -15%). (We always measure size 7 3/4 with System Inner Hand.)

Article number	Side	for
8R1=L7 1/4	Left (L)	8E44=6-L7 1/4 8E44=7-L7 1/4
8R1=L7 3/4	Left (L)	8E44=6-L7 3/4 8E44=7-L7 3/4
8R1=L8 1/4	Left (L)	8E44=6-L8 1/4 8E44=7-L8 1/4
8R1=R7 1/4	Right (R)	8E44=6-R7 1/4 8E44=7-R7 1/4
8R1=R7 3/4	Right (R)	8E44=6-R7 3/4 8E44=7-R7 3/4
8R1=R8 1/4	Right (R)	8E44=6-R8 1/4 8E44=7-R8 1/4

Technical Data

Reference Number	8R1
Consists of:	9S263=R/L 9S264 Threaded Connector with Tube Adapter 9S265 Threaded Connector 9S267=* Lock Ring 9E167 Cable Connector 9E388 Cable Seal 9S266 Chassis with Quick-Disconnect Wrist 9E397=7 Coaxial Bushing, brown, with cable 9E397=8 636W23 UHU-Plus Endfest 300 Adhesive 1 pc. Coding Ring, brown



647H501

1

2

3

4

5

6

7

8

1

2

3

4

5

6

7

8



Myo Elbow Components

1

2

3

4

5

6

7

8



DynamicArm

The DynamicArm is an electronically controlled elbow joint with continuous gearing. Its vario-gear permits a more natural reproduction of a natural arm movement pattern through the continuous adjustment of gear ratio to the environmental conditions (input signal of the patient, movement phase, load to be lifted). Its harmonious movements combine high precision with speed to provide the user with more independence in everyday life. The low level of operating noise and the slightly dampened and completely silent free-swing phase support the natural, inconspicuous appearance.

The shape and basic colour of the DynamicArm also match quite well with the human body and allow the amputation to slip into the background of daily life. Silicone elements discreetly divided into various colours dampen noises and shocks, such as those that may be caused when supporting oneself on a rigid surface.

Integrated into an appealing and unobtrusive design with highest possible function, the DynamicArm provides users with quality for life.

12K100N DynamicArm Facelift

Electronically Controlled Active Elbow Component

For Myoelectrically Controlled Prostheses


With electric motor, continuously electronically controlled Vario-Gear, integrated Li-Ion Battery, Bluetooth® Interface BionicLink, electronically controlled Automatic Forearm Balance (AFB) and upper arm rotation joint (sickle joint) with adjustable friction. Maximal forearm length 305 mm. The DynamicArm is a microprocessor controlled elbow joint driven by an electric motor that can be controlled by the user with great precision. The speed of flexion and extension of the elbow are proportionately controlled. Pronation and supination of the wrist can also be proportionally controlled via the integrated electronics using the Electric Rotator 10S17 (optional). Weights of up to 6 kg can be actively lifted, depending on the forearm length. At a normal level of user activity, the integrated Li-Ion battery with a capacity of 1800 mAh is sufficient for approx. 1 day.

The Bluetooth® Interface BionicLink is integrated into the DynamicArm. In combination with the ElbowSoft Software and the 60X5 Bluetooth adapter, it allows wireless adjustments using a PC. Therefore, the prosthesis system can be optimised, even while the patient is wearing the prosthesis. Nine different programs allow the system to be adapted to various requirements. Moreover, all essential parameters can be set individually.

Through the integration of the AXON Bus®, a unique internal data exchange technology developed by Ottobock, the DynamicArm is largely unaffected by electric and magnetic interference (e.g. mobile phones, high-tension power lines, shopping centre security systems, etc.) This further increases operating safety and reliability.

Article number	Upper arm connection Ø	For hand size	For lamination ring Ø	Colour
12K100N=45	70 mm	7 - 7 1/4	45 mm	No. 4
12K100N=50	70 mm	7 3/4 - 8 1/4	50 mm	No. 4
12K100N=50-1	70 mm	7 3/4 - 8 1/4	50 mm	No. 11
12K100N=50-2	70 mm	7 3/4 - 8 1/4	50 mm	No. 15

Colour roughly corresponds to glove colour according to 646M3 colour swatch



 646D229
646D120

 647G152
646T3=3.3

1

2

3

4

5

6

7

8

Control Modes for the DynamicArm

Pr.	Control	Switching mode	DynamicArm control type	Electric Wrist Rotator	System Electric Hand and System Electric Greifer	Indication
1	Two electrodes	<ul style="list-style-type: none"> Sequential switching with long or short co-contraction with autom. switch-back to the hand Vibration active 	Proportional	Proportional	All control types for two strong muscle signals	For amputees with two strong muscle signals
2	Two electrodes	<ul style="list-style-type: none"> Short co-contraction with autom. switch-back to the hand Vibration active 	Proportional	Not applicable	All control types for two strong muscle signals	For amputees with two strong muscle signals
3	Two electrodes and one 1 switch	<ul style="list-style-type: none"> Sequential switching with switch impulse With autom. switch-back to the hand Vibration active 	Proportional	Proportional	All control types for two strong muscle signals	For amputees with two strong muscle signals
	Two switches	<ul style="list-style-type: none"> Sequential switching with switch impulse With autom. switch-back to the hand Vibration active 	Digital	Digital	All control types for switch control	For amputees with weak or no muscle signals
4	Two electrodes and one switch	<ul style="list-style-type: none"> Switching back and forth with switch impulse With autom. switch-back to the hand Vibration active 	Proportional	Not applicable	All control types for two strong muscle signals	For amputees with two strong muscle signals
	Two switches	<ul style="list-style-type: none"> Switching back and forth with switch impulse With autom. switch-back to the hand Vibration active 	Digital	Not applicable	All control types for switch control	For amputees with weak or no muscle signals
5	Two electrodes and one 4-Step Transducer	<ul style="list-style-type: none"> Direct switching with the impulse of a 4-Step Transducer With autom. switch-back to the hand Vibration active 	Proportional	Proportional	All control types for two strong muscle signals	For amputees with two strong muscle signals
	One switch and one 4-Step Transducer	<ul style="list-style-type: none"> Direct switching with the impulse of a 4-Step Transducer With autom. switch-back to the hand Vibration active 	Digital	Digital	All control types for switch control	For amputees with weak or no muscle signals
6	Two electrodes and one Linear Transducer	Four-channel control	Position control with Linear Transducer	Digital	All control types for two strong muscle signals	For amputees with two strong muscle signals

Pr.	Control	Switching mode	DynamicArm control type	Electric Wrist Rotator	System Electric Hand and System Electric Greifer	Indication
7	Two electrodes and one Linear Transducer	<ul style="list-style-type: none"> • Co-contraction • With autom. switch-back to the hand • Vibration active 	Position control with Linear Transducer	Proportional	All control types with two electrodes	For amputees with two strong muscle signals
8	Two electrodes and one Linear Transducer	Not applicable	Position control with Linear Transducer	Not applicable	All control types with two electrodes	<ul style="list-style-type: none"> • For amputees with two muscle signals of any strength • Possibility for simultaneous control of DynamicArm and Hand/Greifer component
	One switch and one Linear Transducer	Not applicable	Position control with Linear Transducer	Not applicable	All control types for switch control	For amputees with weak or no muscle signals
9	One electrode and one Linear Transducer	Not applicable	Position control with Linear Transducer	Not applicable	All control types with one electrode	<ul style="list-style-type: none"> • For amputees with one strong muscle signal • Possibility for simultaneous control of DynamicArm and Hand/Greifer component

1

2

3

4

5

6

7

8



ErgoArm

The higher the amputation level, the greater the demands on the fitting technique. The two ErgoArm Electronic plus and ErgoArm Hybrid elbow components facilitate myoelectric fitting of high amputation levels. They offer up-to-date technology with the aim of combining low weight and extraordinary function with appealing design.

Due to its electronically-controlled lock, the **12K50 ErgoArm Electronic plus** can be locked or unlocked with myo-signals as well as with a switch. This happens inconspicuously in fractions of a second and offers reliable loading in any position.

The **12K44 ErgoArm Hybrid plus** is recommended for hybrid prostheses with a myoelectric hand and a cable-controlled elbow joint. Due to the integrated "EasyPlug" cable guide, all electric cables disappear inside the prosthesis so that they are inconspicuous and protected.

To better fit the patient, both elbow joints are now available with two different wrist connections (50 mm and 45 mm) and in 3 colours.

12K50 ErgoArm Electronic plus

With EasyPlug (internal electrical connection) and Automatic Forearm Balance (AFB) for myoelectrically controlled prostheses.

With internal electronic lock and EasyPlug (internal electrical connection), Automatic Forearm Balance (AFB) and upper arm rotation joint (sickle joint), with adjustable friction. Plastic forearm shell, length 305 mm, circumference approx. 260 mm, skin-coloured plastic elbow ball.

The continuously-adjustable electronic lock can be locked or unlocked either with myoelectric signals (2 electrodes: co-contraction) or by means of a switch. Various programs that can be accessed via colour-coded connections permit individual adaptation of the lock control system to the user's individual requirements. The slip-stop function makes it possible to lower the forearm in a controlled manner without having to fully release and reactivate the lock (not included with every model!). The lock can bear a load of up to 230 N at a forearm length of 305 mm.

The electrode and battery connecting cables can be plugged into the elbow cover. Since there are no external cables, its overall durability is increased, reducing the risk of defects caused by broken cables and enhancing cosmetic appearance at the same time. The internal infinite-position lock can be unlocked or locked under load in any position.

Article number	Upper arm connection Ø	For hand size	For lamination ring Ø	Colour
12K50=45	70 mm	7 – 7 1/4	45 mm	No. 4
12K50=45-1	70 mm	7 – 7 1/4	45 mm	No. 11
12K50=45-2	70 mm	7 – 7 1/4	45 mm	No. 15
12K50=50	70 mm	7 3/4 – 8 1/4	50	No. 4
12K50=50-1	70 mm	7 3/4 – 8 1/4	50 mm	No. 11
12K50=50-2	70 mm	7 3/4 – 8 1/4	50 mm	No. 15

Colour roughly corresponds to glove colour according to 646M3 colour swatch



647H437

1

2

3

4

5

6

7

8

Switching Modes 12K50

Switching mode		Lock	SLIP-STOP	Four Channel Processor II	Hand
1	white	Press the switch → release = lock Press the switch → release = unlock	NO	All switching modes	All versions
2	red	Co-contraction = lock Co-contraction = unlock	NO	Only Program 1 (white) Four-channel Control or Program 9 Adjustment Cap	All versions with 2 electrodes Not recommended for hands with Digital or Digital-Twin Control System
3	green	Press and hold down switch = elbow mode Electrode OPEN = unlock Electrode CLOSE = lock Release the switch = Hand mode	YES	All switching modes	
4	blue	Press the switch → release = elbow mode Electrode OPEN = unlock Electrode CLOSE = lock Press the switch → release = hand mode	YES		
5	yellow	Press the switch → release = elbow mode Electrode OPEN = unlock Electrode CLOSE = lock 10 s no electrode signal = hand mode or press the switch → release = hand mode	YES		
6	magenta	Co-contraction = elbow mode Electrode OPEN = unlock Electrode CLOSE = lock Co-contraction = hand mode	YES	Only Program 1 (white) Four-channel Control or Program 9 Adjustment Cap	
7	orange	Co-contraction = elbow mode Electrode OPEN = unlock Electrode CLOSE = lock 10 s no electrode signal = hand mode or co-contraction = hand mode	YES		

Vibration feedback signal for successful switching over between hand and elbow is activated in Programs 3–7.

1x vibration = hand mode (the hand is controlled by electrode signals)
2x vibration = elbow mode (the elbow is controlled by electrode signals)

ErgoArm Hybrid plus

12K44 ErgoArm Hybrid plus

With EasyPlug (internal electrical connection) and Automatic Forearm Balance (AFB) for myoelectrically controlled hybrid prostheses.

With internal infinite-position lock and upper arm rotation joint (sickle joint) with adjustable friction. Plastic forearm shell, length 305 mm, circumference approx. 260 mm, skin-coloured plastic elbow ball.

The slip-stop function makes it possible to lower the forearm in a controlled manner without having to fully release and reactivate the lock. The lock can bear a load of up to 230 N at a forearm length of 305 mm.


The electrode and battery connecting cables can be plugged into the elbow cover. Since there are no external cables, its overall durability is increased, reducing the risk of defects caused by broken cables and enhancing cosmetic appearance at the same time.

The internal infinite-position lock can be unlocked or locked under load in any position.

Article number	Upper arm connection Ø	For hand size	For lamination ring Ø	Colour
12K44=45	70 mm	7 – 7 1/4	45 mm	No. 4
12K44=45-1	70 mm	7 – 7 1/4	45 mm	No. 11
12K44=45-2	70 mm	7 – 7 1/4	45 mm	No. 15
12K44=50	70 mm	7 3/4 – 8 1/4	50 mm	No. 4
12K44=50-1	70 mm	7 3/4 – 8 1/4	50 mm	No. 11
12K44=50-2	70 mm	7 3/4 – 8 1/4	50 mm	No. 15

Colour roughly corresponds to glove colour according to 646M3 colour swatch



 647H437

1

2

3

4

5

6

7

8

Accessories for 12K44 and 12K50

1



743A23 Alignment Aid for ErgoArm

Alignment and foaming aid, allows the alignment of an interim prosthesis for functional training in the rehabilitation phase.

2

3



21A207 Clamp Stopples Set

The clamp stopple set is used to connect the 12K48=* or 12K49=* ErgoArm Forearm Component to the elastic strap of the 21A35=1 Triple-Control Above-Elbow Harness.

Article number	Consists of:
21A207	10 pieces clamp stopples 10 pieces short threaded fittings 1 piece twist drill Ø 5.5 mm

4

5




13Z68 Adapter

Adapter for mounting a 12K50 Ottobock Elbow Component to a Hosmer lamination ring.

6

7

8



1
2
3
4
5
6
7
8

1
2
3
4
5
6
7
8



Myo Battery Manager

1

2

3

4

5

6

7

8

Batteries

1



757B13 Interchangeable Battery

For MyoBock 4.8 Volt children's system

Article number	757B13
Weight	40 g
Capacity	1.5 mAh
Average discharge level	4.8 V
Final discharging level	4.4 V
Dimensions LxWxH	69x26.5x14.5 mm
Technology	Nickel metal hydride
for	integration into all socket shapes
Complete with	on/off switch and replaceable integrated fuse

647G47

646D326

2

3

- For charging the 757B13 Interchangeable Battery, only use a 757L13 Pulse Charger.
- Can only be used for MyoBock 4.8 V children's system

4



757B35 MyoEnergy Integral

For MyoBock System

The MyoEnergy Integral 757B35=* is an integrated power supply system made up of several components.

The **charging receptacle** enables the battery contact, indicates the current charge status and serves to switch the prosthesis on and off and release it in an emergency.

The **communication cable** with a 3-pin receptacle is used for the exchange of data.

The **supply cable** establishes the connection between the battery and the respective prosthesis component.

The **battery** consists of two cells with different capacities.

646D442

647G533

Article number	757B35=1	757B35=3
Capacity	600 mAh	1150 mAh
Approx. output voltage	7.4 V	
Approx. charging time	2.5 h	
Technology	Lithium polymer	Lithium-ion
for	System MyoBock	
Consists of:	Battery Lamination Dummy Charging Receptacle Lamination Dummy Drilling Template for Charging Receptacle	

5

6

- The MyoEnergy Integral can be used for adult and children fittings (in combination with the 9E420=* 7in1 Controller).

7

8

757B20/757B21 EnergyPack


For MyoBock System, 6/7.2 Volt,

Compared to NiCd battery, these batteries have a lower self-discharge level, higher cell voltage and greater capacity. They have no memory effect.

Article number	757B20	757B21
Weight	65 g	51 g
Capacity	900 mAh	800 mAh
Charging time (full charge)	ca. 3.5 h	ca. 3.0 h
Technology	Lithium-ion	
for	integration into all socket shapes, even in cases of long residual limbs	

- Can only be used for MyoBock adult systems



 647H356

757B15 X-ChangePack

For MyoBock System, 6 Volt

In comparison with NiCd batteries, NiMh batteries are the same size but have a considerably higher capacity.


Recommended especially for retrofitting existing 6 Volt prosthetic systems.

For this reason, old prosthetic systems can be operated without energy management (the Ottobock 6 Volt System is not suitable for operation with Li-Ion batteries) without problems.

Article number	757B15
Weight	77 g
Capacity	550 mAh
Nominal voltage	6 V
Dimensions	81x28x16 mm
Technology	Nickel metal hydride
for	integration into all socket shapes, even in cases of long residual limbs

- For charging the 757B15 Interchangeable Battery, only use a 757L14 Pulse Charger.
- Can only be used for MyoBock adult systems



 647H482

Battery Receptacle

1



647H6

757Z149 Battery Receptacle Set

Article number	757Z149
for	the 757B13 Interchangeable Battery
Consists of:	757Z150 Battery Receptacle 757Z151 Cable Retention Piece with 627F1=3.1x1.6 O-Ring 757Z152 Casting Template 757Z153 Spacer Insert

2

3



647H369

757Z184=1 Battery Receptacle Set

Article number	757Z184=1
for	757B20 EnergyPack
Consists of:	757Z185=1 Battery Receptacle 757Z186 Locking Lever 757Z188=1 Spacer Insert 757Z187=1 Casting Template 757Z189=1 Foil Template

4

5



757Z185=1 Battery Receptacle

Battery frame without locking lever

Article number	757Z185=1
for	757B20 EnergyPack

6

7



647H492

757Z184=2 Battery Receptacle Set

with screw clamp connection for quick removal from the socket. Integrated connection sockets for electrodes, hand cables or for the 757P41 and 757P39 Connection Cables for connection to the MyoBoy

Article number	757Z184=2
for	757B20 EnergyPack
Consists of:	757Z185=2 Battery Receptacle 757Z186 Locking Lever 757Z187=1 Casting Template 757Z189=1 Foil Template 757Z202 Plug 506G1=M3x10 Set Screw

8

757Z185=2 Battery Receptacle

Battery Receptacle **with screw clamp connection**, without locking lever

Article number	757Z185=2
for	757B20 EnergyPack
Consists of:	757Z202 Plug 506G1=M3x10 Set Screw



1

757Z190=1 Battery Receptacle Set

Article number	757Z190=1
for	757B21 EnergyPack
Consists of:	757Z191=1 Battery Receptacle 757Z186 Locking Lever 757Z193=1 Spacer Insert 757Z192=1 Casting Template 757Z194=1 Foil Template



647H369

2

3

757Z191=1 Battery Receptacle

Battery receptacle without locking lever

Article number	757Z191=1
for	757B21 EnergyPack



4

5

757Z190=2 Battery Receptacle Set

with screw clamp connection for quick removal from socket. Integrated connection sockets for electrodes, hand cables or for 757P41 and 757P39 Connection Cables for connecting to the MyoBoy

Article number	757Z190=2
for	757B21 EnergyPack
Consists of:	757Z191=2 Battery Receptacle 757Z186 Locking Lever 757Z192=1 Casting Template 757Z194=1 Foil Template 757Z202 Plug 506G1=M3x10 Set Screw



647H492

6

7

757Z191=2 Battery Receptacle

Battery Receptacle **with screw clamp connection**, without locking lever

Article number	757Z191=2
for	757B21 EnergyPack
Consists of:	757Z202 Plug 506G1=M3x10 Set Screw



8

1



757Z103=1 Battery Receptacle Set

Article number	757Z103=1
for	757B15 X-ChangePack
Consists of:	757Z104=1 Battery Receptacle 757Z105=1 Cable Retention Piece with O-Ring 757Z106=1 Casting Template 757Z107=1 Spacer Insert

2

647H6

3



757Z103=2 Battery Receptacle Set

The battery is inserted 2.5 mm into device.

Article number	757Z103=2
for	757B15 X-ChangePack
Consists of:	757Z104=2 Battery Receptacle 757Z105=2 Cable Retention Piece with O-Ring 757Z106=2 Casting Template 757Z107=2 Spacer Insert

4

Accessories for 757Z184 and 757Z190 Battery Receptacles

5



757Z186 Locking Lever

6



757Z195 Locking Lever

Article number	757Z195
for	bilateral amputees
Complete with	enlarged release lever for increased user comfort

7

8


Chargers

757L13 Pulse Charger

For simultaneous charging one or two 757B13 Interchangeable Batteries. The case is made of impact and break-resistant plastic material. The design conforms to recommended standards. Includes 757L16-2 Mains Adapter with interchangeable EU and US plug.

Article number	757L13
Weight	220 g
Operating temperature	0 bis +40 °C
Dimensions LxWxH	130x70x45 mm
Supply voltage range	100-240 automatic adjustment V/AC
Mains frequency	40 -70 Hz
Charging current for 757L13	Pulse current technology, average direct current charge approx. 50 mA, reduction to trickle charge
Charging time (for full charge)	ca. 5 h



 647H169

757L35 MyoCharge Integral

The MyoEnergy Integral integrated into the socket is charged using the 757L35 MyoCharge Integral. This is done by simply connecting the charging plug to the charging receptacle on the outside of the socket. Thanks to an integrated magnet, the charging plug can be easily attached to the charging receptacle. The special contour of the charging receptacle and charging plug assures the quick, reliable positioning of the two components to each other. LEDs indicate the readiness of the charger and the current battery capacity.

Article number	757L35
Operating temperature	0 bis +60 °C
Storage temperature	-20 to +60 °C
Mains frequency	50 - 60 Hz
Supply voltage	100 - 240 V



 647G534


- The MyoCharge Integral can be operated for MyoEnergy Integral (757B25=1, 757B35=1 and 757B35=3) models.

757L20 Li-Ion Charger

For charging one or two 757B20/757B21 EnergyPacks. The case is made of impact-resistant plastic material. Includes 757L16-2 mains adapter with EU and additional US plug. The design conforms to recommended standards, protection class II.

Article number	757L20
Operating temperature	0 bis +50 °C
Mains frequency	47 - 63 Hz
Nominal voltage	100 - 240 (automatic adjustment) V/AC
Approx. charging time for EnergyPacks	3.5 (757B20) 3.0 (757B21) h



 647H357

1



2

 647G260

757L14 Pulse Charger

For simultaneously charging one or two 757B15 X-ChangePacks.
The case is made of impact and break-resistant plastic material.
The design conforms to recommended standards, protection class II.
Includes 757L16-2 Mains Adapter with interchangeable EU and US plug.

Article number	757L14
Weight	220 g
Operating temperature	0 bis +40 °C
Dimensions LxWxH	130x70x45 mm
Supply voltage range	100 - 240 V/AC
Mains frequency	40 - 70 Hz
Charging time (for full charge)	ca. 10 h
Charging current for 757L14	Pulse current technology, average direct current charge approx. 70 mA, reduction to trickle charge

3

4

5

6

7


8

Accessories for chargers

4X74 Car Charging Cable 12 V

This cable allows charging of one 757B35=* or two 757B20/757B21 EnergyPacks, 757B15 X-ChangePacks or 757B13 Interchangeable Batteries using the cigarette lighter of a car.

Article number	4X74
for	MyoBock 757L35, 757L20, 757L14 and 757L13 chargers



1

2

757L16-2 Universal Power Supply

Mains unit for powering the MyoBock 757L35, 757L20, 757L14, 757L13 chargers



3

4

757S1=AUS Adapter for Australia



5

757S1=GB Adapter for the UK



6

7

8

Battery Connection

1




13E51=2 Battery Connecting Cable

For connecting a 757B15 X-ChangePack to 9E169 Coaxial Plug, 13E205 MyoRotronic, 13E190 or 13E190=150 Distributor, 10S17 Electric Wrist Rotator (for switch control) or 12K44 ErgoArm Hybrid plus or 12K50 ErgoArm Electronic plus.

Article number	Length	for
13E51=2	200 mm	757B15

2



3




13E51=3 Battery Connecting Cable

For connecting a 757B15 X-ChangePack to the 13E98 Connecting Cable or various switch systems.

Article number	Length	for
13E51=3	200 mm	757B15

4



5




13E51=4 Battery Connecting Cable

For connecting the 757B15 X-ChangePack to the 9E169 Coaxial Plug, 10S17 Electric Wrist Rotator or 13E205 MyoRotronic

Article number	Length	for
13E51=4	600 mm	757B15

6



7




13E188 Battery Connecting Cable

For connecting the 757Z185=1 or 757Z191=1 Battery Receptacle to 9E169 Coaxial Plug, 10S17 Electric Wrist Rotator or 13E205 MyoRotronic or the 13E190 or 13E190=150 distributor

Article number	Length	for
13E188=200	200 mm	757B20/757B21
13E188=600	600 mm	757B20/757B21

8



13E190 Distributor

For connecting the 757B20 or 757B21 EnergyPack, 13E200 or 13E202 Electrodes and System Electric Hands to **9E53 Hand Cable**.
Without extension cable.

Article number	for
13E190	757B20/757B21



- Apply 633F11 Silicone Grease to the bushings prior to connection.
Protects against corrosion!

13E190=150 Distributor

For connecting the 757B20 or 757B21 EnergyPack, 13E200 or 13E202 Electrodes and System Electric Hands to **9E53 Hand Cable and extension cable**.

Article number	Length	for
13E190=150	150 mm	757B20/757B21



- Apply 633F11 Silicone Grease to the bushings prior to connection.
Protects against corrosion!

1

2

3

4

5

6

7

8

1

2

3

4

5

6

7

8



Myo Control Elements

1

2

3

4


5

6

7

8

Electrodes



647G334

13E202 Suction Socket Electrode

Electrode cable connection with IDC termination, with 13E206 Electrode Accessories, without electrode cable.

This generation of electrodes are based on the well-known 13E200 Electrode. Embedded into a mounting suspension of elastic material, this electrode creates an airtight seal between the inner socket and outer socket.

If used correctly, the 13E202 Suction Socket Electrode also prevents sweat from penetrating between the outer and inner socket, therefore effectively preventing damage to the electrical and mechanical components caused by corrosion.

The Suction Socket Electrode can not only be used for standard fittings, but is also suitable for application in suction sockets. Combining the suction socket electrode with a 12V10 Tube Valve for a suction socket creates a vacuum effect in the socket, ensuring optimal linkage between the residual limb and the socket.

As with the 13E200 Electrode, state-of-the-art shielding and filtering technologies largely protect the 13E202 Suction Socket Electrode against high frequency interference caused, for example, by cell phones, walkie-talkies, computers or anti-theft systems in shopping centres so that the correct function of the myoelectrically controlled prosthesis is not affected.

The electrode contacts are made from **pure titanium** and as they do not contain nickel, they are also suitable for people with allergies.

The frequency filter's full protection effect will only be provided if the mains frequency and filter frequency are identical.

Article number	13E202=50	13E202=60
HZ	50	60
Frequency bandwidth	90 - 450 Hz	
Ambient temperature	-15 bis +60 °C	
Operating voltage U	4.8 - 7.2 V	

- Use 633F11 Silicone Grease to seal the plug connections.
- Remove any excessive grease after connecting the electrode cable.
- For accessories for vacuum forming of inner sockets, see pages 122, 335, 337
- 12V10 Tube Valve, see Page 336

13E200=* Electrode

Electrode cable connection with IDC termination, with **13E153 Electrode Accessories**

These MyoBock electrodes are particularly sensitive in the range of low muscle signals. The change in amplification now takes place logarithmically, which enables enhanced differentiation of the signal level, particularly in the range of high muscle signals.

Thanks to modern frequency shielding and filtering technologies, it is significantly less sensitive to low and high frequency interferences that are emitted, for example, by mobile phones or shopping centre security systems.

The electrode contacts are made from **pure titanium** and as they do not contain nickel, they are also suitable for people with allergies.

The frequency filter's full protection effect will only be provided if the mains frequency and filter frequency are identical.

Article number	13E200=50	13E200=60
Weight	4.5 g	
HZ	50	60
Operating voltage	4.8 - 7.2 V	
Dimensions LxWxH	27x18x9.5 mm	
Frequency bandwidth	90 - 450 Hz	
Ambient temperature	-15 bis +60 °C	

- Use 633F11 Silicone Grease to seal the plug connections. Remove any excessive grease after connecting the electrode cable.
- For accessories for vacuum forming of inner sockets, see pages 122, 335, 337



647H490

1

2

3

4

5

6

7

8

Accessories for electrodes

1



13E206 Electrode Accessory Set

Electrode Accessories

Article number	13E206
for	13E202 Suction Socket Electrode
Consists of:	13E203 Lamination Template for Inner Socket 13E204 Lamination template for Outer Socket 507S15 Lamination Disk, serrated 503F3 Socket Screw with Allen head 13E80 Sensitivity Adjustment Tool

2

3



13E153 Electrode Accessories

Article number	13E153
for	laminated inner sockets with 13E200 Electrode
Consists of:	13E191 Pattern for inner shell 13E192 Lamination template for Outer Socket 507S15 Lamination Disk, serrated 503F3 Socket Screw with Allen head 13E80 Sensitivity Adjustment Tool

4

- Use a 709S10=2 Allen Head Wrench for 503F3 Socket Screw when connecting the inner and outer socket.

5



13E201 Electrode Accessories

Article number	13E201
for	vacuum-formed inner sockets with 13E200 Electrode
Consists of:	Pattern for inner shell 13E172 Electrode mounting bracket (skin colours) Lamination dummy for electrode mounting bracket 29C5=M4x9 Setting Nut

6

- The 13E201 Electrode Accessories are only available as a set
- The 13E172 Electrode Mounting Bracket and 29C5=M4x9 Setting Nut can be ordered individually

7



13E135 Electrode Mounting Bracket Set

For positioning and assembling the MYOBOCK electrodes to the interim plaster socket or the interim socket of Thermolyn (616T52 or 616T53).

Article number	13E135
for	13E200 Electrodes

8

Control Elements

9X50/9X51 Control Element

For assembly in a harness system


The **9X50 Linear Control Element** allows continuous, proportional control of prosthesis components using harnesses.

The **4-step 9X51 Control Element** allows for proportional control of prosthetic components at four different speed levels using harnesses. In combination with the 12K100 DynamicArm, it provides the possibility to selectively switch between the different system components.

Reference Number	9X50/9X51
Weight	6 g
Cable travel	8 mm
Max. operating force	10 N

• For 13E129=G* Connecting Cable, see Page 127



 647H475

9X52/9X53 Control Element

For assembly between the outer and inner prosthetic socket


The **9X52 Linear Control Element** allows continuous, proportional control of prosthesis components using harnesses.

The **4-step 9X53 Control Element** allows for proportional control of prosthetic components at four different speed levels using harnesses. In combination with the 12K100 DynamicArm, it provides the possibility to selectively switch between the different system components.

Reference Number	9X52/9X53
Weight	11 g
Cable travel	8 mm
Max. operating force	10 N

• For 13E129=G* Connecting Cable, see Page 127



 647H485

1

2

3

4

5

6

7

8

Switches

1



647G400

9X14 Harness Pull Switch

To be used in a harness or as a cable pull switch

For controlling the 10S17 Electric Wrist Rotator, the 8E38=7, 8E38=8, 8E38=9 and 8E12 System Electric Hand, or the 8E33=7 and 8E33=9 System Electric Greifer.

The control cables can be plugged in after lifting the cover of the switch.

Functional sequence: off – function 1, – off – function 2.

Article number	9X14
Weight	19 g
Scope of Delivery	501T16=M2x6 Cap Screw (4 pieces)

• For Connecting Cable, see pages 126-127

2

3



647G401

9X18 Cable Pull Switch

With steel cable, bow and wedge lock.

Standard use: the switch housing is screwed onto the socket, the steel cable is connected to the harness or another pull cable for controlling the 10S17 Electric Wrist Rotator, the 8E38=7, 8E38=8, 8E38=9 and 8E12 System Electric Hand or the 8E33=* System Electric Greifer.

Functional sequence: off – function 1, – off – function 2.

Article number	9X18
Weight	17 g
Scope of Delivery	501S46=M3x5 Oval Head Countersunk Screw (2 pieces) 501S46=M3x8 Oval Head Countersunk Screw (2 pieces)

• For Connecting Cable, see pages 126-127

4

5



647G402

9X25 Rocker Switch

With laterally attached four-conductor flat cable with plug connection,

for controlling the 10S17 Electric Wrist Rotator, the

8E38=7, 8E38=8 and 8E12, System Electric Hand, the 8E33=* System Electric Greifer.

Article number	9X25
with	501S46=M2x8 Oval Head Countersunk Screw (2 pieces) 501S46=M2x5 Oval Head Countersunk Screw (2 pieces)

- Use a 13E99 Connection Cable for opening and closing.
- Use 13E50 or 13E97 Connection Cable to control pronation and supination.
- Use the 13E98 Connection Cable to control opening and closing of the 8E12 System Electric Hand.
- For Connecting Cable, see pages 126-127

7

8

9X37 Pressure Switch

The pressure switch makes it possible to control the 8E38=* System Electric Hand, the 8E33=* System Electric Greifer or the 10S17 Electric Wrist Rotator in connection with the 13E205 MyoRotronic.

The switch's special feature allows the separate adjustment of both control points.

Article number	9X37
with	501S75=M2x8 Oval Head Countersunk Screw (2 pieces) 501S75=M2x5 Oval Head Countersunk Screw (2 pieces)

- Use a 13E99 Connection Cable for opening and closing.
- For Connecting Cable, see Page 127



 647H71

1

2

3

4

5


6

7

8

Connection Possibilities

1



13E50 Connecting Cable

To control electric pronation and supination with the 10S17 Electric Wrist Rotator. The connection cable provides the electrical connection between the wrist rotator and the 9X14 Harness Pull Switch, 9X18 Cable Pull Switch, or 9X25 Rocker Switch.


Article number	Length
13E50=250	250 mm
13E50=1200	1200 mm

2



3

4



13E97 Connection Cable with On-Off Switch

To control electric pronation and supination with the 10S17 Electric Wrist Rotator. Pronation and supination function may be turned off using the switch. The connection cable provides the electrical connection between the electric wrist rotator and 9X14 Harness Pull Switch, 9X18 Cable Pull Switch or 9X25 Rocker Switch.

Article number	Length
13E97=250	250 mm
13E97=1200	1200 mm

5



6

7



13E98=1200 Connecting Cable

For controlling opening and closing of the 8E12 System Electric Hand. The connection cable provides the electrical connection between the System Electric Hand and the 9X14 Harness Pull Switch, 9X18 Cable Pull Switch or 9X25 Rocker Switch.

Article number	Length
13E98=1200	1,200 mm

8

13E99=1200 Connecting Cable

For controlling opening and closing of the 8E38=9, 8E38=7, 8E38=8 System Electric Hand or the 8E33=* System Electric Greifer.

The connection cable provides the electrical connection between the coaxial plug or electric wrist rotator and the 9X14 Harness Pull Switch, 9X18 Cable Pull Switch, 9X25 Rocker Switch or 9X37 Pressure Switch.

Article number	Length
13E99=1200	1,200 mm



1

2

9E185 Extension Cable

For elongation of the hand cable of the 8E39, 8E41, 8E44 System Electric Hand or 8E34 System Electric Greifer and connection to the 757Z185=2 or 757Z191=2

Article number	Length
9E185=30	300 mm
9E185=40	400 mm
9E185=50	500 mm



3

4

- Apply 633F11 Silicone Grease to the bushings prior to connection.

757P41 Connection Cable

For connecting the 757Z185=2 or 757Z191=2 Battery Receptacle and the 9E169 Coaxial Plug or the 13E205 MyoRotronic.



5

13E129 Electrode Cable with Straight Plug and 13E121 Plug Connector

For connecting the electrode and the 9E369/9E370 4in1 LS Controller or 9E420 7in1 Controller, the 9E169 Coaxial Plug, 13E190 Distributor or 10S17 Electric Wrist Rotator or 13E205 MyoRotronic

Article number	Length
13E129=G100	100 mm
13E129=G300	300 mm
13E129=G600	600 mm
13E129=G1000	1000 mm



6

7

8

1
2
3
4
5
6
7
8



Myo Master

1

2

3

4

5

6

7

8

Control for systems designed for children



646D442 647G595

7in1 Controller

The 9E420 7in1 Controller is a control system with 7.4 Volt for the Electric Hand 2000 that processes muscle signals and sends them to the prosthetic hand. The 7in1 Controller provides 7 programmes, some of which have been adopted from the adult system.

Children can benefit from the range of possibilities offered by systems for adults, their prosthetic hand can be adapted more individually and they are simultaneously prepared for a subsequent fitting with the adult system.

Numerous Control Programmes

7 programme modes can be used to ideally adjust to the patient's needs.

DMC plus:	Control with 2 electrodes
DMC LowInput:	Control with 2 electrodes
AutoControl LowInput:	Control with 2 electrodes
DigitalControl:	Control with 2 electrodes
VarioControl:	Control with 1 electrode
DoubleChannel:	Control with 1 electrode
EVO Digital:	Control with 1 electrode

Power Supply

Power is supplied to the 7in1 Controller and the prosthetic components by the 757B25=1 or 757B35=* MyoEnergy Integral.

Article number	Side
9E420=L	Left (L)
9E420=R	Right (R)



646D326 647H209

4in1 Controller LS

The 4in1 Controller LS extends fitting possibilities in the area of long forearm residual limbs and wrist disarticulations. Four different controls can be selected by various coloured coded plugs. The 4in1 Controller LS is required for a fitting with an 8E51 Electric Hand 2000.

Article number	For hand size
9E369	5
9E370	5 1/2, 6, 6 1/2

Technical Data

Article number	9E369	9E370
Consists of:	9E371 Contact Plate 9E372 Lamination Pattern 9E373Controller for 9E369 9E374 Controller for 9E370	9E371 Contact Plate 9E372 Lamination Pattern 9E373Controller for 9E369 9E374 Controller for 9E370

Accessories for 4in1 Controller LS

13E182 Coding Plug Set

The 4in1 Controller LS differentiates between 4 functions and between right and left arm prostheses using colour coded plugs:

Article number	Side	Colour	9E370	9E369
13E184=1	Left (L)	white	2 Elektroden-digital	1 Elektrode-digital (EVO)
13E184=2	Right (R)	red	2 Elektroden-digital	1 Elektrode-digital (EVO)
13E184=3	Left (L)	green	2 Elektroden-DMC	2 Elektroden-DMC
13E184=4	Right (R)	blue	2 Elektroden-DMC	2 Elektroden-DMC
13E184=5	Left (L)	yellow	2 Elektroden-DMC LowInput	2 Elektroden-DMC LowInput
13E184=6	Right (R)	purple	2 Elektroden-DMC LowInput	2 Elektroden-DMC LowInput
13E184=7	Left (L)	orange	1 Elektrode-EVO	1 Elektrode-DMC Low (EVO)
13E184=8	Right (R)	black	1 Elektrode-EVO	1 Elektrode-DMC Low (EVO)

Digital: "Open/Close" by two electrodes DMC: proportional control by two electrodes Low Input: as with DMC, only for children with lower muscular voltage from approx. 20 μ V EVO: proportional "opening" by one electrode, automatic "closing" function via motor



- The eight coded plugs are also available separately.

1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



Myo Software

1

2

3

4

5

6

7

8

1

2

3



646C48 PAULA with MyoBoy

PAULA is a comprehensive software package providing essential support for fabricating upper limb prostheses. The software assists certified prosthetists with the planning, socket design and construction of myoelectrically controlled prostheses, hybrid prostheses, cable-controlled prostheses and passive prostheses. The integration into the Ottobock Data Station provides a common platform, offering prosthetists a familiar user interface.

- 646F277 Myo Poster
- 646F286 Movo-Poster

646D315

647G461

4

5

6



646C52 PAULA

PAULA is a comprehensive software package providing essential support for fabricating upper limb prostheses. The software assists certified prosthetists with the planning, socket design and construction of myoelectrically controlled prostheses, hybrid prostheses, cable-controlled prostheses and passive prostheses. The integration into the Ottobock Data Station provides a common platform, offering prosthetists a familiar user interface.

646D315

646G461

7

8

757M11=X-Change MyoBoy

The MyoBoy serves to provide optimum support for the prosthetist and patient when training muscle activity and for the realistic simulation of the MyoBock systems.

The data that are collected allow the most suitable control system for the individual fitting to be selected.

As standard equipment, the MyoBoy is included in the scope of delivery of the Ottobock PAULA software (646C48=*). The MyoBoy can be ordered separately under article number 757M11=X-Change as a replacement or exchange.




Article number	757M11=X-Change
Scope of Delivery	757G265=1 Quick Reference Guide 757Z18 Grounding Cable 757P44 Electrode Adapter 757Z174 Electrode Armband Screwdriver USB Cable Carrying Case
	Optional: 757P23 Test Adapter 757P39 MyoBoy Communication Cable

646C42 ElbowSoft

The "ElbowSoft" software is to be used exclusively for adjusting the DynamicArm and any Ottobock System components connected to it (8E38=* System Electric Hands, 8E33=* System Electric Greifers, 10S17 Electric Wrist Rotator) within the scope of an exoprosthetic fitting. Participation in an Ottobock product training is mandatory prior to using the product. To qualify for software updates, additional product training courses may be necessary.



 647G341



1

2

3

4

5

6

7

8

1



560X3 MyolinoSoft

In order to offer patient-specific controller settings, setup software is also being used for the first time in a hand system for children. Rather than using specified settings, it works with various parameters, such as the patient side, control program, switching threshold and gripping speed that can be individually adjusted. This is especially advantageous for children with high mobility requirements.

The 560X3 MyolinoSoft Software is designed for setting the 7in1 9E420 Controller and the Electric Hand 2000 connected to it. 7 programme versions are available.

2

647G578=V1.0

3

4



60X6 MyolinoLink

The MyolinoLink 60X6 is used for wireless data communication between Ottobock products (e.g. 7in1 Controller) and a PC with Bluetooth® functionality.

The MyolinoLink is equipped with an LED. When the LED is lit, this indicates that the connection between the MyolinoLink and PC has been established successfully. When the LED is flashing, this indicates that data is being exchanged between the MyolinoLink and the PC.

Patient-specific settings can only be administered with the Ottobock MyolinoSoft software.

5



647G594

6



646C50=V1.2 MovoPAULA

MovoPAULA is a comprehensive software package providing essential support for the fabrication of cable-controlled and passive upper limb prostheses. The software assists prosthetists with the planning, socket design and construction of the prosthesis. Thanks to the integration into the Ottobock Data Station, a common platform has been created, offering the prosthetist a familiar user interface.

MovoPAULA is particularly well suited for prosthetists who exclusively carry out fittings with conventional upper limb prostheses and who are not myo-certified.

7

8

647G463

646C51=V1.2 TED & MyoBoy

The TED software guides therapists through patient management, thereby enabling accurate documentation of the entire fitting process. In addition to myo-prosthesis training, the software includes a component selection for the future prosthesis. All steps performed by the therapist can be documented and placed at the disposal of the decision makers who are involved in the fitting. All data can be transferred to the full version of PAULA, thus facilitating cooperation with prosthetists.



 647G462

1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



Movo Terminal Devices

1

2

3


4

5

6

7

8



Otto Bock System Hands

The System Hand, developed by Otto Bock, was designed to meet the functional needs of the patient. The System Hand consists of three parts – the hand chassis with hand mechanism, the shaped inner hand, and the cosmetic glove – providing both functionality and a natural appearance.

The gripping function depends on the construction of the hand mechanism:

Passive system hands are used for **passive arm prostheses**. They are opened by the sound hand and close automatically.

The single-cable and double-cable System Hands belong to the cable-activated prostheses or "active prehensile arms" and are activated by a harness control system.

The **single-cable system** is opened by pulling the cable and closes automatically.

The **double-cable system** is closed by pulling the control cable. Another pull releases the lock and opens the hand. Connection to the forearm is provided by the threaded stud or in the case of long residual limbs by the chassis.

The system inner hand encloses the hand mechanism and provides the shape of the prosthetic glove, which is selected according to the size of the inner hand and is available in various models and colours.

Ottobock System Hand -passive-

With threaded stud and system inner hand,


Hand sizes:

6 3/4	for children
7 1/4	for adolescents and women
7 3/4	for men
8	for men

The passive system hand can be used for all lengths of the residual limb to be fitted with a cosmetic prosthesis. It is opened with the sound hand and closes independently. The system is lightweight and stable.

- The Prosthetic glove must be ordered separately. See pages 172-174, 214-216



 647G444

Article number	Side	Threaded stud	Size	Inner hand
8K18=L6 3/4	Left (L)	M12x1,5	6 3/4	8X14=L6 3/4
8K18=L7 1/4	Left (L)	M12x1,5	7 1/4	8X14=L7 1/4
8K18=L7 3/4	Left (L)	M12x1,5	7 3/4	8X14=L7 3/4
8K18=L8	Left (L)	M12x1,5	8	8X14=L8
8K18=R6 3/4	Right (R)	M12x1,5	6 3/4	8X14=R6 3/4
8K18=R7 1/4	Right (R)	M12x1,5	7 1/4	8X14=R7 1/4
8K18=R7 3/4	Right (R)	M12x1,5	7 3/4	8K21=L8
8K18=R8	Right (R)	M12x1,5	8	8X14=R8

Technical Data

Reference Number	8K18	8K18	8K18	8K18
Size	6 3/4	7 1/4	7 3/4	8
Weight	ca. 185 g	ca. 250 g	ca. 280 g	ca. 290 g

Article number	Side	Threaded stud	Size	Inner hand
8K19=L6 3/4	Left (L)	1/2"-20	6 3/4	8X14=L6 3/4
8K19=L7 1/4	Left (L)	1/2"-20	7 1/4	8X14=L7 1/4
8K19=L7 3/4	Left (L)	1/2"-20	7 3/4	8X14=L6 3/4
8K19=L8	Left (L)	1/2"-20	8	8X14=L8
8K19=R6 3/4	Right (R)	1/2"-20	6 3/4	8X14=R6 3/4
8K19=R7 1/4	Right (R)	1/2"-20	7 1/4	8X14=R7 1/4
8K19=R7 3/4	Right (R)	1/2"-20	7 3/4	8X14=R7 3/4
8K19=R8	Right (R)	1/2"-20	8	8X14=R8

Technical Data

Reference Number	8K19	8K19	8K19	8K19
Size	6 3/4	7 1/4	7 3/4	8
Weight	ca. 185 g	ca. 250 g	ca. 280 g	ca. 290 g

Spare parts for passive system hands

1



8X14 System Inner Hand

Article number	For hand size	Side
8X14=L6 3/4	6 3/4	Left (L)
8X14=L7 1/4	7 1/4	Left (L)
8X14=L7 3/4	7 3/4	Left (L)
8X14=L8	8	Left (L)
8X14=R6 3/4	6 3/4	Right (R)
8X14=R7 1/4	7 1/4	Right (R)
8X14=R7 3/4	7 3/4	Right (R)
8X14=R8	8	Right (R)

Technical Data

Reference Number	8X14
Complete with	finger support in the little and ring finger and 9S187 Retainer Ring

2

3

4



9S15 Retainer Ring, narrow

Article number	For hand size	Ø
9S15=42	6 3/4	42

5



9S187 Retainer Ring, wide

Article number	For hand size
9S187=7 1/4	7 1/4
9S187=7 3/4	7 3/4, 8

Technical Data

Article number	9S187=7 1/4	9S187=7 3/4
for	8X14 System Inner Hands	8X14 System Inner Hands

6

7



9S6 Finger and Thumb Tip

Article number	For hand size
9S6	6 3/4, 7 1/4, 7 3/4 and 8

8

9S67 Chassis

With threaded stud

Article number	Side	Ø	Threaded stud	for
9S67=L40	Left (L)	40 mm	M12x1,5	8K18=L6 3/4
9S67=L44-N	Left (L)	44 mm	M12x1,5	8K18=L7 1/4
9S67=L48-N	Left (L)	48 mm	M12x1,5	8K18=L7 3/4 8K18=L8
9S67=R40	Right (R)	40 mm	M12x1,5	8K18=R6 3/4
9S67=R44-N	Right (R)	44 mm	M12x1,5	8K18=R7 1/4
9S67=R48-N	Right (R)	48 mm	M12x1,5	8K18=R7 3/4 8K18=R8



1

2

9S65 Chassis

With threaded stud

Article number	Side	Ø	Threaded stud	for
9S65=L40	Left (L)	40 mm	1/2"-20	8K19=L6 3/4
9S65=L44-N	Left (L)	44 mm	1/2"-20	8K19=L7 1/4
9S65=L48-N	Left (L)	48 mm	1/2"-20	8K19=L7 3/4 8K19=L8
9S65=R40	Right (R)	40 mm	1/2"-20	8K19=R6 3/4
9S65=R44-N	Right (R)	44 mm	1/2"-20	8K19=R7 1/4
9S65=R48-N	Right (R)	48 mm	1/2"-20	8K19=R7 3/4 8K19=L8



3

4

5

501T52=M3x6 Flange Button Head Socket Screw



6

501T52=M3x8 Flange Button Head Socket Screw



7

8

Accessories for passive system hands

1



9E94 Lamination Ring

Article number	For hand size	Ø
9E94=44	6 3/4	44 mm
9E94=50	7 1/4, 7 3/4, 8	50 mm

2

3



Chassis

For long below-elbow and through-wrist residual limbs

Article number	For hand size	Ø	Suitable for
9S96=40	6 3/4	40 mm	9E94 Lamination Ring
9S184=44-N	7 1/4	44 mm	9E94 Lamination Ring
9S184=48-N	7 3/4, 8	48 mm	9E94 Lamination Ring

4



Technical Data

Article number	9S96=40	9S184=44-N	9S184=48-N
for	replacement of the chassis with threaded stud	replacement of the chassis with threaded stud	replacement of the chassis with threaded stud
Suitable for	9E94 Lamination Ring	9E94 Lamination Ring	9E94 Lamination Ring

5



506G4 Set Screw

(4 set screws are required per chassis)

Article number	For chassis
506G4=M3x5	9S96=40
506G4=M4x5	9S184

6

7

8

8K20 Ottobock System Hand –voluntary opening–

With **perlon cable** through the palm of the hand (**inner pull**),
with threaded stud and system inner hand

These system hands are suitable for all cable-activated prostheses. They are opened by pulling the cable (active) and close independently with concurrent locking.

Hand sizes:

6 3/4 for children
7 1/4 for adolescents and women
7 3/4 for men
8 for men

Article number	Side	Threaded stud	Size	Inner hand
8K20=L6 3/4	Left (L)	M12x1,5	6 3/4	8X14=L6 3/4
8K20=L7 1/4	Left (L)	M12x1,5	7 1/4	8X14=L7 1/4
8K20=L7 3/4	Left (L)	M12x1,5	7 3/4	8X14=L7 3/4
8K20=L8	Left (L)	M12x1,5	8	8X14=L8
8K20=R6 3/4	Right (R)	M12x1,5	6 3/4	8X14=R6 3/4
8K20=R7 1/4	Right (R)	M12x1,5	7 1/4	8X14=R7 1/4
8K20=R7 3/4	Right (R)	M12x1,5	7 3/4	8X14=R7 3/4
8K20=R8	Right (R)	M12x1,5	8	8X14=R8

Technical Data

Reference Number	8K20	8K20	8K20	8K20
Size	6 3/4	7 1/4	7 3/4	8
Weight	ca. 215 g	ca. 300 g	ca. 330 g	ca. 340 g



647G445

◦ The prosthetic glove must be ordered separately. See pages 172-174, 214-216

1

2

3


4


5

6

7

8



 647G445

8K22 Single-Cable System Hand

With **perlon cable** through the back of the hand (**outer pull**), with threaded stud and system inner hand

These system hands are suitable for all cable-activated prostheses. They are opened by pulling the cable (active) and close independently with concurrent locking.

Hand sizes:

6 3/4 for children
7 1/4 for adolescents and women
7 3/4 for men
8 for men

Article number	Side	Threaded stud	Size	Inner hand
8K22=L6 3/4	Left (L)	M12x1,5	6 3/4	8X18=L6 3/4
8K22=L7 1/4	Left (L)	M12x1,5	7 1/4	8X14=L7 1/4
8K22=L7 3/4	Left (L)	M12x1,5	7 3/4	8X14=L7 3/4
8K22=L8	Left (L)	M12x1,5	8	8X14=L8
8K22=R6 3/4	Right (R)	M12x1,5	6 3/4	8X14=R6 3/4
8K22=R7 1/4	Right (R)	M12x1,5	7 1/4	8X14=R7 1/4
8K22=R7 3/4	Right (R)	M12x1,5	7 3/4	8X14=R7 3/4
8K22=R8	Right (R)	M12x1,5	8	8X14=R8

Technical Data

Reference Number	8K22	8K22	8K22	8K22
Size	6 3/4	7 1/4	7 3/4	8
Weight	ca. 215 g	ca. 300 g	ca. 330 g	ca. 340 g

• The prosthetic glove must be ordered separately. See pages 172-174, 214-216

8K21 Single-Cable System Hand

With **steel cable** through the palm of the hand (**inner pull**),
with threaded stud and system inner hand

These system hands are suitable for all cable-activated prostheses. They are opened by pulling the cable (active) and close independently with concurrent locking.

Hand sizes:

7 1/4 for adolescents and women

7 3/4 for men

8 for men

Article number	Side	Threaded stud	Size	Inner hand
8K21=L7 1/4	Left (L)	1/2"-20	7 1/4	8X14=L7 1/4
8K21=L7 3/4	Left (L)	1/2"-20	7 3/4	8X14=L7 3/4
8K21=L8	Left (L)	1/2"-20	8	8X14=L8
8K21=R7 1/4	Right (R)	1/2"-20	7 1/4	8X14=R7 1/4
8K21=R7 3/4	Right (R)	1/2"-20	7 3/4	8X14=R7 3/4
8K21=R8	Right (R)	1/2"-20	8	8X14=R8

Technical Data

Reference Number	8K21	8K21	8K21
Size	7 1/4	7 3/4	8
Weight	ca. 300 g	ca. 330 g	ca. 340 g

• The prosthetic glove must be ordered separately. See pages 172-174, 214-216



 647G445

1

2

3

4

5

6

7

8



647G445

8K23 Single-Cable System Hand

With **steel cable** through the back of the hand (**outer pull**),
with threaded stud and system inner hand

These system hands are suitable for all cable-activated prostheses. They are opened by pulling the cable (active) and close independently with concurrent locking.

Hand sizes:

6 3/4	for children
7 1/4	for adolescents and women
7 3/4	for men
8	for men

Article number	Side	Threaded stud	Size	Inner hand
8K23=L6 3/4	Left (L)	1/2"-20	6 3/4	8X14=L6 3/4
8K23=L7 1/4	Left (L)	1/2"-20	7 1/4	8X14=L7 1/4
8K23=L7 3/4	Left (L)	1/2"-20	7 3/4	8X18=R8 1/4
8K23=L8	Left (L)	1/2"-20	8	8X14=L8
8K23=R6 3/4	Right (R)	1/2"-20	6 3/4	8X14=R6 3/4
8K23=R7 1/4	Right (R)	1/2"-20	7 1/4	8X14=R7 1/4
8K23=R7 3/4	Right (R)	1/2"-20	7 3/4	8X14=R7 3/4
8K23=R8	Right (R)	1/2"-20	8	8X14=R8

Technical Data

Reference Number	8K23	8K23	8K23	8K23
Size	6 3/4	7 1/4	7 3/4	8
Weight	ca. 215 g	ca. 300 g	ca. 330 g	ca. 340 g

• The prosthetic glove must be ordered separately. See pages 172-174, 214-216

8K24 Double-Cable System Hand

With **perlon cable** through the palm of the hand (**inner pull**),
with threaded stud and system inner hand

The hand closes following a deliberate pull of the cable. A subsequent adjustment pull increases the grip force and locks in any grip position (active). Through renewed activation of the cable, the hand becomes released and opens independently (double cable).

Hand sizes:


6 3/4 for children
7 1/4 for adolescents and women
7 3/4 for men
8 for men

Article number	Side	Threaded stud	Size	Inner hand
8K24=L7 1/4	Left (L)	M12x1,5	7 1/4	8X14=L7 1/4
8K24=L7 3/4	Left (L)	M12x1,5	7 3/4	8X14=L7 3/4
8K24=L8	Left (L)	M12x1,5	8	8X14=L8
8K24=R7 1/4	Right (R)	M12x1,5	7 1/4	8X14=R7 1/4
8K24=R7 3/4	Right (R)	M12x1,5	7 3/4	8X14=R7 3/4
8K24=R8	Right (R)	M12x1,5	8	8X14=R8

Technical Data

Reference Number	8K24	8K24	8K24
Size	7 1/4	7 3/4	8
Weight	ca. 340 g	ca. 370 g	ca. 380 g, ca. 380 g.



 647G446

ⓘ The prosthetic glove must be ordered separately. See pages 172-174, 214-216

1

2

3

4


5

6

7

8



 647G446

8K26 Double-Cable System Hand

With **perlon cable** through the back of the hand (**outer pull**), with threaded stud and system inner hand

The hand closes following a deliberate pull of the cable. A subsequent adjustment pull increases the grip force and locks in any grip position (active). Through renewed activation of the cable, the hand becomes released and opens independently (double cable).

Hand sizes:

6 3/4 for children
7 1/4 for adolescents and women
7 3/4 for men
8 for men

Article number	Side	Threaded stud	Size	Inner hand
8K26=L7 1/4	Left (L)	M12x1,5	7 1/4	8X14=L7 1/4
8K26=L7 3/4	Left (L)	M12x1,5	7 3/4	8X14=L7 3/4
8K26=L8	Left (L)	M12x1,5	8	8X14=L8
8K26=R7 1/4	Right (R)	M12x1,5	7 1/4	8X14=R7 1/4
8K26=R7 3/4	Right (R)	M12x1,5	7 3/4	8X14=R7 3/4
8K26=R8	Right (R)	M12x1,5	8	8X14=R8

Technical Data

Reference Number	8K26	8K26	8K26
Size	7 1/4	7 3/4	8
Weight	ca. 340 g	ca. 370 g	ca. 380 g

- The prosthetic glove must be ordered separately. See pages 172-174, 214-216

8K27 Double-Cable System Hand

With **steel cable** through the back of the hand (**outer pull**),
with threaded stud and system inner hand

The hand closes following a deliberate pull of the cable. A subsequent adjustment pull increases the grip force and locks in any grip position (active). Through renewed activation of the cable, the hand becomes released and opens independently (double cable).

Hand sizes:

6 3/4 for children
7 1/4 for adolescents and women
7 3/4 for men
8 for men

Article number	Side	Threaded stud	Size	Inner hand
8K27=L7 1/4	Left (L)	1/2"-20	7 1/4	8X14=L7 1/4
8K27=L7 3/4	Left (L)	1/2"-20	7 3/4	8X14=L7 3/4
8K27=L8	Left (L)	1/2"-20	8	8X14=L8
8K27=R7 1/4	Right (R)	1/2"-20	7 1/4	8X14=R7 1/4
8K27=R7 3/4	Right (R)	1/2"-20	7 3/4	8X14=R7 3/4
8K27=R8	Right (R)	1/2"-20	8	8X14=R8

Technical Data

Reference Number	8K27	8K27	8K27
Size	7 1/4	7 3/4	8
Weight	ca. 340 g	ca. 370 g	ca. 380 g

• The prosthetic glove must be ordered separately. See pages 172-174, 214-216



647G446

1

2

3

4

5

6

7

8

Spare parts for cable system hands

1



8X14 System Inner Hand

Article number	For hand size	Side
8X14=L6 3/4	6 3/4	Left (L)
8X14=L7 1/4	7 1/4	Left (L)
8X14=L7 3/4	7 3/4	Left (L)
8X14=L8	8	Left (L)
8X14=R6 3/4	6 3/4	Right (R)
8X14=R7 1/4	7 1/4	Right (R)
8X14=R7 3/4	7 3/4	Right (R)
8X14=R8	8	Right (R)

Technical Data

Reference Number	8X14
Complete with	finger support in the little and ring finger and 9S187 Retainer Ring

2

3

4



9S15 Retainer Ring, narrow

Article number	For hand size	Ø
9S15=42	6 3/4	42

5



9S187 Retainer Ring, wide

Article number	For hand size
9S187=7 1/4	7 1/4
9S187=7 3/4	7 3/4, 8

Technical Data

Article number	9S187=7 1/4	9S187=7 3/4
for	8X14 System Inner Hands	8X14 System Inner Hands

6

7

8

Chassis

With threaded stud

Article number	Ø	Threaded stud	for
9S10=40	40 mm	M12x1,5	8K18=L/R6 3/4 8K20=L/R6 3/4 8K22=L/R6 3/4 8K24=L/R6 3/4 8K26=L/R6 3/4
9S166=44-N	44 mm	M12x1,5	8K18=L/R7 1/4 8K20=L/R7 1/4 8K22=L/R7 1/4 8K24=L/R7 1/4 8K26=L/R7 1/4
9S166=48-N	48 mm	M12x1,5	8K18=L/R7 3/4 8K18=L/R8 8K20=L/R7 3/4 8K20=L/R8 8K22=L/R7 3/4 8K22=L/R8 8K24=L/R7 3/4 8K24=L/R8 8K26=L/R7 3/4 8K26=L/R8
9S52=40	40 mm	1/2"-20	8K19=L/R6 3/4 8K23=L/R6 3/4 8K27=L/R6 3/4
9S185=44-N	44 mm	1/2"-20	8K19=L/R7 1/4 8K21=L/R7 1/4 8K23=L/R7 1/4 8K27=L/R7 1/4
9S185=48-N	48 mm	1/2"-20	8K19=L/R7 3/4 8K19=L/R8 8K21=L/R7 3/4 8K21=L/R8 8K23=L/R7 3/4 8K23=L/R8 8K27=L/R7 3/4 8K27=L/R8



1

2

3

4

5

9S6 Finger and Thumb Tip

Article number	For hand size
9S6	6 3/4, 7 1/4, 7 3/4 and 8



6

21A8 Suspension Rosette



7

8

1



21A17=3.8x1000 Cable-Control Unit Spiral

Article number	Ø
21A17=3.8x1000	3.8 mm

2

Accessories for cable system hands

3



9E94 Lamination Ring

Article number	For hand size	Ø
9E94=44	6 3/4	44 mm
9E94=50	7 1/4, 7 3/4, 8	50 mm

4



Chassis

For long below-elbow and through-wrist residual limbs

Article number	For hand size	Ø	Suitable for
9S96=40	6 3/4	40 mm	9E94 Lamination Ring
9S184=44-N	7 1/4	44 mm	9E94 Lamination Ring
9S184=48-N	7 3/4, 8	48 mm	9E94 Lamination Ring

5



Technical Data

Article number	9S96=40	9S184=44-N	9S184=48-N
for	replacement of the chassis with threaded stud	replacement of the chassis with threaded stud	replacement of the chassis with threaded stud
Suitable for	9E94 Lamination Ring	9E94 Lamination Ring	9E94 Lamination Ring

6

7



506G4 Set Screw

(4 set screws are required per chassis)

Article number	For chassis
506G4=M3x5	9S96=40
506G4=M4x5	9S184

8

501T52=M3x6 Flange Button Head Socket Screw



1

501T52=M3x8 Flange Button Head Socket Screw



2

10A43 Attachment Plate with stud

For connecting a **cable-activated hook or system hand** to the **10V8 Wrist Unit**

Article number	Ø	Interior thread
10A43=M12x1.5	50 mm	M12x1,5
10A43=1/2"-20	50 mm	1/2"-20

Technical Data

Reference Number	10A43
Material	stainless steel
Complete with	interior thread



3

4

10A44 Threaded Stud

For connecting a hook to the **10V9, 10V10** or **10V25** Wrist Units without ratchet

Article number	Interior thread
10A44=M12x1.5	M12x1,5
10A44=1/2"-20	1/2"-20

Technical Data

Reference Number	10A44
Material	stainless steel
Complete with	interior thread M12x1,5



5

6

7

8

1



10A56 Attachment Plate with short stud

attachment plate Ø 50 mm,

For connecting a **cable-activated hook or system hand** to the **10V30 Wrist Unit**

Article number	Ø	Interior thread
10A56=M12x1.5	20 mm	M12x1,5
10A56=1/2"-20	20 mm	1/2"-20

2

Technical Data

Reference Number	10A56
Material	stainless steel
Complete with	interior thread M12x1,5

3


4

5

6

7

8



Cable-activated Hooks

In cable-activated prostheses, a hook, instead of a prosthetic hand, can also be used as the terminal device. It is opened by pulling a harness and is closed automatically by spring or rubber elements.

In addition to the cable-activated hooks for children and adolescents, the standard hooks for adults as well as all-purpose hooks have been included in this catalogue.

The various elements for attaching the cable-activated hooks to the wrist construction are described at the end of this section.

1

2

3

4

5

6

7

8

Cable-activated Hooks for Children and Adolescents

1



647G443

10A25 Cable-activated Hook for Children

Article number	Side	Threaded stud
10A25=LM12x1.5	Left (L)	M12x1,5
10A25=RM12x1.5	Right (R)	M12x1,5
10A25=L1/2"-20	Left (L)	1/2"-20
10A25=R1/2"-20	Right (R)	1/2"-20

Technical Data

Reference Number	10A25
Material	aluminum alloy

- Connecting elements for cable-activated hooks are listed on pages 153, 155-156, 165-167, 183-184.

2

3

4



647G443

10A37 Cable-Activated Hook for Youths

Article number	Side	Threaded stud
10A37=LM12x1.5	Left (L)	M12x1,5
10A37=RM12x1.5	Right (R)	M12x1,5
10A37=L1/2"-20	Left (L)	1/2"-20
10A37=R1/2"-20	Right (R)	1/2"-20

Technical Data

Reference Number	10A37
Material	aluminum alloy

- Connecting elements for cable-activated hooks are listed on pages 153, 155-156, 165-167, 183-184.

5

6

7

8

Accessories for 10A25 and 10A37

10Y1 Finger Cover

Article number	10Y1
Material	Plastic
Colour	skin colour



10Y8 Finger Cover

Article number	10Y8
Material	Plastic
Colour	skin colour
Scope of Delivery	1 set comprises one cover with or without knobs



10Y13 Rubber Band



1

2

3

4

5

6

7

8

Cable-Activated Hooks for Adults

1
2
3
4
5
6
7
8



647G204

10A70 MovoHook 2Grip

Article number	Side	Threaded stud
10A70=LM12X1.5	Left (L)	M12X1.5
10A70=RM12X1.5	Right (R)	M12X1.5
10A70=L1/2"-20	Left (L)	1/2"-20
10A70=R1/2"-20	Right (R)	1/2"-20

Technical Data

Reference Number	10A70
Weight	120 g
Material	Aluminium
Opening width	100 mm

- Connecting elements for cable-activated hooks are listed on pages 153, 155-156, 165-167, 183-184.



647G204

10A71 MovoHook 2Grip

This coating offers superior gripping security, especially in instances where metal has a tendency to be very slippery on certain materials, such as paper or glass. Furthermore, the coating does not adhere to textiles, so sliding through a sleeve while getting dressed or undressed does not pose a problem. The material is hygienically clean. Resistance to chemicals can be tested on a case-by-case basis upon request. The coating is permanent and inconspicuous. Finger covers are not required, providing yet another advantage.

Article number	Side	Threaded stud
10A71=LM12X1.5	Left (L)	M12X1.5
10A71=RM12X1.5	Right (R)	M12X1.5
10A71=L1/2"-20	Left (L)	1/2"-20
10A71=R1/2"-20	Right (R)	1/2"-20

Technical Data

Reference Number	10A71
Weight	120 g
Material	Aluminium with branch coating for an even better grip
Opening width	100 mm

- Connecting elements for cable-activated hooks are listed on pages 153, 155-156, 165-167, 183-184.

10A80 MovoHook 2Grip

Article number	Side	Threaded stud
10A80=LM12X1.5	Left (L)	M12X1.5
10A80=RM12X1.5	Right (R)	M12X1.5
10A80=L1/2"-20	Left (L)	1/2"-20
10A80=R1/2"-20	Right (R)	1/2"-20

Technical Data

Reference Number	10A80
Weight	250 g
Material	Stainless steel
Opening width	100 mm

- Connecting elements for cable-activated hooks are listed on pages 153, 155-156, 165-167, 183-184.



647G204

10A81 MovoHook 2Grip

The familiar MovoHook 2Grip cable-activated hooks are now also available with coated branches for an even better grip. This coating offers superior gripping security, especially in instances where metal has a tendency to be very slippery on certain materials, such as paper or glass. Furthermore, the coating does not adhere to textiles, so sliding through a sleeve while getting dressed or undressed does not pose a problem. The material is hygienically clean. Resistance to chemicals can be tested on a case-by-case basis upon request. The coating is permanent and inconspicuous. Finger covers are not required, providing yet another advantage.

Article number	Side	Threaded stud
10A81=LM12X1.5	Left (L)	M12X1.5
10A81=RM12X1.5	Right (R)	M12X1.5
10A81=L1/2"-20	Left (L)	1/2"-20
10A81=R1/2"-20	Right (R)	1/2"-20

Technical Data

Reference Number	10A81
Weight	250 g
Material	Stainless steel with branch coating for an even better grip
Opening width	100 mm

- Connecting elements for cable-activated hooks are listed on pages 153, 155-156, 165-167, 183-184.



647G204

1



647G442

10A11 Cable-activated Hook for Adults

Standard hook

Article number	Side	Threaded stud
10A11=LM12x1.5	Left (L)	M12x1,5
10A11=RM12x1.5	Right (R)	M12x1,5
10A11=L1/2"-20	Left (L)	1/2"-20
10A11=R1/2"-20	Right (R)	1/2"-20

Technical Data

Reference Number	10A11
Material	aluminum alloy

- Connecting elements for cable-activated hooks are listed on pages 153, 155-156, 165-167, 183-184.

2

3

4



647G442

10A60 Cable-activated Hook for Adults

Article number	Side	Threaded stud
10A60=LM12x1.5	Left (L)	M12x1,5
10A60=RM12x1.5	Right (R)	M12x1,5
10A60=L1/2"-20	Left (L)	1/2"-20
10A60=R1/2"-20	Right (R)	1/2"-20

Technical Data

Reference Number	10A60
Material	stainless steel

- Connecting elements for cable-activated hooks are listed on pages 153, 155-156, 165-167, 183-184.

5

6



647G442

10A18 Cable-activated Hook for Adults

Standard hook, model 58

Article number	Side	Threaded stud
10A18=LM12x1.5	Left (L)	M12x1,5
10A18=RM12x1.5	Right (R)	M12x1,5
10A18=L1/2"-20	Left (L)	1/2"-20
10A18=R1/2"-20	Right (R)	1/2"-20

Technical Data

Reference Number	10A18
Material	aluminum alloy

- Wrist connections for cable-activated hooks are listed on page

7

8

10A12 All-purpose Hook


Article number	Side	Threaded stud
10A12=L	Left (L)	ohne
10A12=R	Right (R)	ohne
10A12=LM12x1.5	Left (L)	M12x1,5
10A12=RM12x1.5	Right (R)	M12x1,5
10A12=L1/2"-20	Left (L)	1/2"-20
10A12=R1/2"-20	Right (R)	1/2"-20

Technical Data

Reference Number	10A12
Material	stainless steel

- Connecting elements for cable-activated hooks are listed on pages 153, 155-156, 165-167, 183-184.



 647G442

1

2

3

4

5

6

7

8

Spare parts for cable-activated hooks

1



10Y1 Finger Cover

Article number	10Y1
Material	Plastic
Colour	skin colour

2


3



10Y2 Double Spring

4

5



10Y12 Spring

3 pieces are required to complete renew all springs.

Article number	10Y12
for	10A12 All-purpose Hook

6

7

8

Accessories for cable-activated hooks

21A13 Connection Piece for Cable-Activated Hook

Article number	for	Package
21A13=1	21A7	Perlon Cable
21A13=2	21A7	Steel Cable
21A13=3	10Y32	Perlon Cable
21A13=4	10Y32	Perlon Cable

Technical Data

Reference Number	21A13
Complete with	With 21A9 Connector Screw and 21A10 Washer
Suitable for	10A11, 10A12 and 10A60



21A14 Connection Piece

Article number	21A14
for	10A18, 10A25 and 10A37 in conjunction with 21A19/20=*
Complete with	21A8 Clamp Sleeve



1

2

3

4

5

6

7

8

1



21A44 Connection Piece for Hook

Article number	21A44
for	10A12, 10A18, 10A25 and 10A37 in conjunction with 21A35/36=*

2




21A8 Suspension Rosette

3



21A9 Connector Screw

4



21A10 Washer

5



10A43 Attachment Plate with stud

For connecting a **cable-activated hook or system hand** to the **10V8 Wrist Unit**

Article number	Ø	Interior thread
10A43=M12x1.5	50 mm	M12x1,5
10A43=1/2"-20	50 mm	1/2"-20

7

Technical Data

Reference Number	10A43
Material	stainless steel
Complete with	interior thread

8

10A44 Threaded Stud

For connecting a hook to the **10V9, 10V10** or **10V25** Wrist Units without ratchet

Article number	Interior thread
10A44=M12x1.5	M12x1,5
10A44=1/2"-20	1/2"-20

Technical Data

Reference Number	10A44
Material	stainless steel
Complete with	interior thread M12x1,5



1

2

10A56 Attachment Plate with short stud

attachment plate Ø 50 mm,

For connecting a **cable-activated hook or system hand** to the **10V30 Wrist Unit**

Article number	Ø	Interior thread
10A56=M12x1.5	20 mm	M12x1,5
10A56=1/2"-20	20 mm	1/2"-20

Technical Data

Reference Number	10A56
Material	stainless steel
Complete with	interior thread M12x1,5



3

4

5

6

7

8

Implements

1



10A2 All-purpose Clamp

Article number	Complete with	BPL-Pos.
10A2	standard attachment	BPL-Pos. 2.01.14.04

2

3



10A3 All-purpose Hook

Nickel-plated

Article number	Complete with
10A3	standard attachment

4

 647G455

5



10A4 All-purpose Ring

Nickel-plated

Article number	Complete with
10A4	standard attachment

6

 647G455

7




10A8 Movable All-purpose Ring

Nickel-plated

Article number	Complete with
10A8	standard attachment

8

 647G455

10V15 Standard Connector

for accommodating **work tools** with a standard stud (\varnothing 13 mm) and for connecting to wrist units with or without ratchets

Article number	Complete with
10V15	attachment plate and stud



647G455

1

2

11S33 Wrist Lock with Standard Connector

For accommodating work tools with a standard stud (\varnothing 13 mm) and attaching to the Ottobock System Wrist Unit

Article number	11S33=40	11S33=44	11S33=48
Outer \varnothing	40 mm	44 mm	48 mm



647G455

3

4

10A22 Hand Brush

Article number	Complete with	BPL-Pos.
10A22	two suction cups	BPL-Pos. 2.01.14.48



5

6

7

8

1
2
3
4
5
6
7
8



Movo Prosthetic Gloves

1

2

3

4

5

6

7

8

1



8S6=170x65 Prosthetic Glove for Children

For hand sizes 6 3/4, with long sleeve.
Standard prosthetic glove available in 18 different shades

Article number	For hand size	Side	Sleeve length
8S6=170x65L	6 3/4	Left (L)	300 mm
8S6=170x65R	6 3/4	Right (R)	300 mm

When ordering, please include the colour code

2

3

- In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.
- To order them, please enter an **N** before the =:
example 8S11=210x78xL4 standard glove
8S11**N**=210x78xL4 MyoSkin Natural
Please ask our customer service for the available colours.

The entire range of colours for Skin Natural models comprises of six colours. The new 646M47 colour scale helps when selecting the right pattern.

4



8S4=190x76 Prosthetic Glove for Adolescents

For hand sizes 7 1/4, with short sleeve.
Standard prosthetic glove available in 18 different shades

Article number	For hand size	Side	Sleeve length
8S4=190x76L	7 1/4	Left (L)	215 mm
8S4=190x76R	7 1/4	Right (R)	215 mm

When ordering, please include the colour code

5

6

- In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.
- To order them, please enter an **N** before the =:
example 8S11=210x78xL4 standard glove
8S11**N**=210x78xL4 MyoSkin Natural
Please ask our customer service for the available colours.

The entire range of colours for Skin Natural models comprises of six colours. The new 646M47 colour scale helps when selecting the right pattern.

7

8

8S4=210x78 Prosthetic Glove for Men

For hand sizes 7 3/4, with short sleeve.
Standard prosthetic glove available in 18 different shades

Article number	For hand size	Side	Sleeve length
8S4=210x78L	7 3/4	Left (L)	220 mm
8S4=210x78R	7 3/4	Right (R)	220 mm

When ordering, please include the colour code

- In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.
- To order them, please enter an **N** before the =:
example 8S11=210x78xL4 standard glove
8S11**N**=210x78xL4 MyoSkin Natural
Please ask our customer service for the available colours.

The entire range of colours for Skin Natural models comprises of six colours. The new 646M47 colour scale helps when selecting the right pattern.



1

2

3

8S4=220x80 Prosthetic Glove for Men


For hand sizes 8, with short sleeve.
Standard prosthetic glove available in 18 different shades

Article number	For hand size	Side	Sleeve length
8S4=220x80L	8	Left (L)	225 mm
8S4=220x80R	8	Right (R)	225 mm

When ordering, please include the colour code

- In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.
- To order them, please enter an **N** before the =:
example 8S11=210x78xL4 standard glove
8S11**N**=210x78xL4 MyoSkin Natural
Please ask our customer service for the available colours.

The entire range of colours for Skin Natural models comprises of six colours. The new 646M47 colour scale helps when selecting the right pattern.



4

5

6

7

8

1



8S5=195x78 Prosthetic Glove for Women

For hand sizes 7 1/4, with long sleeve.

Standard prosthetic glove available in 18 different shades

Article number	For hand size	Side	Sleeve length
8S5=195x78L	7 1/4	Left (L)	340 mm
8S5=195x78R	7 1/4	Right (R)	340 mm

When ordering, please include the colour code

2

- In addition to the standard prosthetic glove, Otto Bock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.

To order them, please enter an **N** before the =:

example 8S11=210x78xL4 standard glove

8S11**N**=210x78xL4 MyoSkin Natural

Please ask our customer service for the available colours.

3

The entire range of colours for Skin Natural models is covered with six colours. The new 646M47 colour scale helps when selecting the right pattern.


4

5

6

7

8



1
2
3
4
5
6
7
8

1

2

3

4

5

6

7

8



Movo Wrist Units

1

2

3

4

5

6

7

8

1



10V39 MovoWrist Flex

Enables flexion and extension of a cable-controlled or passive terminal device;
Can be locked in 5 positions from -15° to $+45^{\circ}$
and rotation over 360° with 20 different positions
Total length: 33 mm, of which 12 mm is distally visible

647G375

Article number	10V39=45	10V39=50
Outer Ø	45 mm	50 mm

2

Accessories for MovoWrist Flex

3



10A30 Adapter

For connecting a terminal device to a MovoWrist Flex

Article number	10A30=M12x1.5	10A30=1/2"-20	10A30=1
for	Terminal devices with metric thread	Terminal devices with inch thread	Working tools with standard stud

4



5

6



11D1 Laminationring

Please note the diameter of the joint.

11D1=45
11D1=50

Article number	11D1=45	11D1=50
Diameter	45 mm	50 mm
for	MovoWrist Flex 10V39	

7

8

Wrist Units

10V18 Ottobock Wrist Unit

Article number	10V18=28	10V18=34	10V18=40	10V18=45	10V18=50
Outer Ø	28 mm	34 mm	40 mm	45 mm	50 mm
Thread	M12x1,5				
Complete with	interior thread and cylindrical lamination ring				



 647G453

10V36 Ottobock Wrist Unit

Article number	10V36=34	10V36=45	10V36=50
Outer Ø	34 mm	45 mm	50 mm
Thread	1/2"-20		
Complete with	interior thread and cylindrical lamination ring		




 647G453

10V32 Ottobock Ball Wrist Unit

Article number	10V32=40	10V32=45	10V32=50
Outer Ø	40 mm	45 mm	50 mm
Thread	M12x1,5		
Complete with	interior thread and cylindrical lamination ring		



 647G453

10V34 Ottobock Ball Wrist Unit

Article number	10V34=45	10V34=50
Outer Ø	45 mm	50 mm
Thread	1/2"-20	
Complete with	interior thread and cylindrical lamination ring	




 647G453

10V9 Wrist Unit without Ratchet

Joint Piece

Article number	10V9
Material	plastic
Diameter	50 mm



 647G452

- For connecting to a system hand to 10A13=50 attachment with stud or 10A44 Threaded Stud
- For connecting to hook with 10A44 Threaded Stud

1



647G452

10V10 Wrist Unit without Ratchet

Joint Piece

Article number	10V10
Material	light metal
Diameter	50 mm

- For connecting to a system hand with 10A44 Threaded Stud
- For connecting to hook with 10A44 Threaded Stud

2



647G452

10V25 Wrist Unit without Ratchet

With flexion device

Article number	10V25
Diameter	50 mm
Complete with	Flexion Device
Consists of:	10V10 Wrist Unit without ratchet + 10V26 Flexion Device

- For connecting to the system hand with 10A44 Threaded Stud
- For connecting to hook with 10A44 Threaded Stud

3

4



647G451

10V8 Wrist Unit – ratchet type rotation

Article number	10V8
Diameter	50 mm
Complete with	lock lever

- Connection to hand with 10A43 Attachment Plate with threaded stud
- Connection to hook with 10A43 Attachment Plate with threaded stud

5

6



647G451

10V30 Wrist Unit – ratchet type rotation, short

Article number	10V30
Diameter	50 mm
for	long below-elbow residual limbs
Complete with	lock lever

- Connection to hand with 10A56 Attachment Plate with short threaded stud
- Connection to hook with 10A56 Attachment Plate with short threaded stud

7

8

Spare parts for wrist units


11D20 Cylindrical Lamination Ring

Article number	11D20=28	11D20=34	11D20=40	11D20=45	11D20=50
Outer Ø	28 mm	34 mm	40 mm	45 mm	50 mm
for	10V18=28	10V18=34	10V18=40 10V32=40	10V18=45 10V32=45 10V34=45 10V36=45	10V8 10V9 10V10 10V18=50 10V32=50 10V34=50 10V36=50



501S40 Oval Head Screw (sheet metal screw)

Article number	For lamination ring Ø
501S40=3.5x9.5	28, 34 mm
501S40=3.5x13	40, 45, 50 mm



11D27 Rubber Friction Ring

Article number	For lamination ring Ø	Ø
11D27=25	28 mm	25 mm
11D27=32	34, 40, 45, 50 mm	32 mm

Technical Data

Article number	11D27=25	11D27=32
for	10V18 and 10V36	10V18 and 10V36



501S27 Slotted Oval Head Screw

Article number	501S27=M3x8
Material	
for	10V8 10V9 10V10 10V25



2:1

Accessories for wrist units

1



10V26 Flexion Device

For subsequent installation into 10V9/10V10 Wrist Unit without ratchet or 10V8 Wrist Unit with ball ratchet (not for 10V30)

2

Article number	10V26
for	10V8 10V9 10V10

3

4

5

6

7

8

Connecting elements for system hands

10A43 Attachment Plate with stud

For connecting a cable-activated hook or system hand to the 10V8 Wrist Unit

Article number	Ø	Interior thread
10A43=M12x1.5	50 mm	M12x1,5
10A43=1/2"-20	50 mm	1/2"-20

Technical Data

Reference Number	10A43
Material	stainless steel
Complete with	interior thread



1

2

10A56 Attachment Plate with short stud

attachment plate Ø 50 mm,

For connecting a cable-activated hook or system hand to the 10V30 Wrist Unit

Article number	Ø	Interior thread
10A56=M12x1.5	20 mm	M12x1,5
10A56=1/2"-20	20 mm	1/2"-20

Technical Data

Reference Number	10A56
Material	stainless steel
Complete with	interior thread M12x1,5



3

4

10A44 Threaded Stud

For connecting a hook to the 10V9, 10V10 or 10V25 Wrist Units without ratchet

Article number	Interior thread
10A44=M12x1.5	M12x1,5
10A44=1/2"-20	1/2"-20

Technical Data

Reference Number	10A44
Material	stainless steel
Complete with	interior thread M12x1,5



5

6

7

8

Connecting elements for cable-activated hooks and implements

1



2

10A56 Attachment Plate with short stud

attachment plate Ø 50 mm,
For connecting a **cable-activated hook or system hand** to the **10V30 Wrist Unit**

Article number	Ø	Interior thread
10A56=M12x1.5	20 mm	M12x1,5
10A56=1/2"-20	20 mm	1/2"-20

3

Technical Data

Reference Number	10A56
Material	stainless steel
Complete with	interior thread M12x1,5

4



10A43 Attachment Plate with stud

For connecting a **cable-activated hook or system hand** to the **10V8 Wrist Unit**

Article number	Ø	Interior thread
10A43=M12x1.5	50 mm	M12x1,5
10A43=1/2"-20	50 mm	1/2"-20

5

Technical Data

Reference Number	10A43
Material	stainless steel
Complete with	interior thread

6



10V15 Standard Connector

for accommodating **work tools** with a standard stud (Ø 13 mm) and for connecting to wrist units with or without ratchets

Article number	10V15
Complete with	attachment plate and stud

7

647G455

8

The page features a large grid of small dots for taking notes. On the right side, there is a vertical sidebar with a light beige background for the top three sections and a teal background for the bottom five sections. The numbers 1 through 8 are printed in white on this sidebar, corresponding to the sections of the dot grid.

1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



Movo Elbow Joints

1

2

3

4

5

6

7

8

Elbow Components for Children

1
2
3
4
5
6
7
8



647G469

12K19=40 Elbow Set-Up for Children

With manual elbow lock (10 locking positions in 8° increments) and upper arm rotation joint (sickle joint) with lamination ring.
Plastic forearm shell, skin colour, length approx. 250 mm, circumference 210 mm.
The elbow ball is made of skin-coloured plastic.

Article number	Upper arm connection Ø	For hand size	Colour	Length
12K19=40	54 mm	6 3/4	skin-coloured	

Technical Data

Article number	12K19=40
Material	Plastic
Ø-Lamination ring	40 mm



Practical recommendation:

- This joint can also be fitted in conjunction with the 2000 8E51=* Childrens Hand.



647G570

12K12 MovolinoArm Friction

The 12K12 MovolinoArm Friction is available in one size. It is the perfect complement to the current Ottobock product portfolio, since it allows fittings on the right and left side for children 3 to 5 years old. It is compatible with passive, cable-controlled and myoelectric arm prosthesis systems. It weighs only 182g.

The elbow has one friction setting for humeral rotation and one for flexion or extension of the forearm. Parents can easily set this friction setting.

Yet another advantage is that the elbow is compatible with the components of the 7.4 Volt Children System.

The MovolinoArm Friction features an impressive design that has a very realistic natural appearance.

The hand adapter of the 12K12 MovolinoArm Friction is naturally compatible with the Electric Hand 2000. The 10A40 Wood Adapter is available for passive prosthetic fittings. The 10V18=34 or 10V36=34 Wrist Unit is compatible with body-powered fittings.

Article number	Upper arm connection Ø	Wrist connection Ø
12K12	43.5 mm	34 mm



Practical recommendation:

- This joint can also be fitted in conjunction with the 2000 8E51=* Childrens Hand.

Spare parts for 12K12 and 12K19

13G21 Lamination Ring Set

Article number	13G21
Consists of:	one lamination ring with lamination dummies and O-Ring

10S18=40 Lamination Ring for Children's Forearm

The 10S18=40 Lamination Ring for Child's Forearm enables an 2000 8E51 Electric Hand to be combined with the 12K19=40 Elbow Set-up for Children. To achieve this, the lamination ring is directly glued into the elbow set-up using 636K18=1 Orthocryl sealing resin compact glue.

Article number	Ø	Elbow Set-Up
10S18=40	33 mm	12K19=40



1

709Z9 Special Key

The 709Z9 Special Key is used to tighten the counter nut. The counter nut may become damaged when tightening in the vice or, for example, with pliers.



2

3

13G8=54 Lamination Ring



4

5

21A18 Nylon wire

Article number	Ø	Length	Order by
21A18=2X1	2 mm	1 m	lfm
21A18=2X5	2 mm	5 m	lfm
21A18=2X10	2 mm	10 m	lfm
21A18=2X25	2 mm	25 m	lfm



6

Accessories for 12K12

11D12=34 Wrist Unit Lamination Ring

For connection with the Children's Hand 2000

Article number	Complete with
11D12=34	copper friction ring, O-Rings and protective cover



7

8

1



10A40 Wood Hand Adapter

for connecting the passive inner hand with the forearm socket

Article number	Complete with
10A40	threaded stud made of plastic, M12x1.5

2

10V18=34 Ottobock Wrist Unit

Article number	10V18=34
Outer Ø	34 mm
Thread	M12x1,5

3



10V36=34 Ottobock Wrist Unit

Article number	10V36=34
Outer Ø	34
Thread	1/2"-20

4

5

6

7

8

12K42 ErgoArm plus

For cable-controlled prostheses.

With internal infinite-position lock, Automatic Forearm Balance (AFB) and upper arm rotation joint

(sickle joint), with adjustable friction.

Plastic forearm shell, plastic elbow ball, skin colour, length 305 mm, circumference 250 mm.

The slip-stop function makes it possible to lower the forearm in a controlled manner without having to fully release and reactivate the lock. The lock can bear a load of up to 230 N at a forearm length of approximately 305 mm. The internal infinite-position lock can be unlocked or locked under load in any position.

The Automatic Forearm Balance supports flexion of the joint and enables natural swing characteristics.

Article number	Upper arm connection Ø	For hand size	For lamination ring Ø	Colour
12K42=45	70 mm	6 3/4 – 7 1/4	45 mm	No. 4
12K42=45-1	70 mm	6 3/4 – 7 1/4	45 mm	No. 11
12K42=45-2	70 mm	6 3/4 – 7 1/4	45 mm	No. 15
12K42=50	70 mm	7 3/4 – 8 1/4	50 mm	No. 4
12K42=50-1	70 mm	7 3/4 – 8 1/4	50 mm	No. 11
12K42=50-2	70 mm	7 3/4 – 8 1/4	50 mm	No. 15

Colour roughly corresponds to glove colour according to 646M3 colour swatch



674H437

12K41 ErgoArm

For cable-controlled prostheses.

With internal infinite-position lock and upper arm rotation joint (sickle joint) with adjustable friction.

Plastic forearm shell, plastic elbow ball, skin colour, length 305 mm, circumference 250 mm.

The slip-stop function makes it possible to lower the forearm in a controlled manner without having to fully release and reactivate the lock. The lock can bear a load of up to 230 N at a forearm length of 305 mm. The internal infinite-position lock can be unlocked or locked under load in any position.

Article number	Upper arm connection Ø	For hand size	For lamination ring Ø	Colour
12K41=45	70 mm	6 3/4 – 7 1/4	45 mm	No. 4
12K41=45-1	70 mm	6 3/4 – 7 1/4	45 mm	No. 11
12K41=45-2	70 mm	6 3/4 – 7 1/4	45 mm	No. 15
12K41=50	70 mm	7 3/4 – 8 1/4	50 mm	No. 4
12K41=50-1	70 mm	7 3/4 – 8 1/4	50 mm	No. 11
12K41=50-2	70 mm	7 3/4 – 8 1/4	50 mm	No. 15

Colour roughly corresponds to glove colour according to 646M3 colour swatch



674H438

- Please note that the variants -1 and -2 are fabricated only after receipt of an order and that, for this reason, delivery of such articles will take longer.

Spare parts for 12K42 + 12K41

1



12K48 Forearm

Article number	Ø	Suitable for
12K48=45	45 mm	12K42=45 12K44=45 12K50=45
12K48=50	50 mm	12K42=50 12K44=50 12K50=50

2

3



12K49 Forearm

Article number	Ø	Suitable for
12K49=45	45 mm	12K41=45
12K49=50	50 mm	12K41=50

4

5



12A13 Elbow Joint

Article number	Colour
12A13-1	No. 11
12A13-15	No. 15

• Marked articles are available in colour shade no. 11 (please add -1 to the end of the order number) or in colour shade no. 15 (add -2 to the end of the order number).

6



13G68 Clamp Plate

7



13Z47 Lamination Ring

8

13Z48 Ball Cap

Article number	Colour
13Z48-1	No. 11
13Z48-2	No. 15



- Marked articles are available in colour shade no. 11 (please add -1 to the end of the order number) or in colour shade no. 15 (add -2 to the end of the order number).

13Z50 Thread Segment

Article number	Colour
13Z50-1	No. 11
13Z50-2	No. 15



- Marked articles are available in colour shade no. 11 (please add -1 to the end of the order number) or in colour shade no. 15 (add -2 to the end of the order number).

13Z51 Hole Covering

Article number	Colour
13Z5-1	No. 11
13Z5-2	No. 15



- Marked articles are available in colour shade no. 11 (please add -1 to the end of the order number) or in colour shade no. 15 (add -2 to the end of the order number).

13Z52 Spring Telescope



13Z53 Switch Cable



1

2

3

4

5

6

7

8

1



13Z55 Lamination Protection Cover



2



13Z56 Tongue Cap

Article number	Colour
13Z56-1	No. 11
13Z56-2	No. 15

3

4



13Z57 Pressure Piece

Article number	Colour
13Z57-1	No. 11
13Z57-2	No. 15

5



13Z58 Eccentric

Article number	Colour
13Z58-1	No. 11
13Z58-2	No. 15

6

7



13Z59 Lamination Protection Cover

Article number	Colour
13Z59-1	No. 11
13Z59-2	No. 15

8



501S84=M4x20 Countersunk Head Screw

501S101=M4x12 Countersunk Head Screw



1

627F13=60x2.5 O-Ring



2

3

4

5


6

7

8

1

Accessories for 12K42 + 12K41



743A23 Alignment Aid for ErgoArm

Alignment and foaming aid, allows the alignment of an interim prosthesis for functional training in the rehabilitation phase.

2

3



21A207 Clamp Stopples Set

The clamp stopple set is used to connect the 12K48=* or 12K49=* ErgoArm Forearm Component to the elastic strap of the 21A35=1 Triple-Control Above-Elbow Harness.

Article number	Consists of:
21A207	10 pieces clamp stopples 10 pieces short threaded fittings 1 piece twist drill Ø 5.5 mm

4

5



13Z68 Adapter

Adapter for mounting a 12K50 Ottobock Elbow Component to a Hosmer lamination ring.

6

7

8

Elbow Component with Cable Lock

12K27 Elbow Component

With single-sided, rotatable cable lock (18 locking positions in about 7.2° increments) and upper arm rotation joint (sickle joint) with lamination ring.

Plastic forearm shell, skin colour, length 280 mm, circumference approx. 260 mm or 300 mm, skin-coloured plastic elbow ball.

Thanks to the exterior joint construction and special technology for interior and exterior rotation of the forearm, this joint allows the distal residual limb to reach up to the joint ball. It is therefore suitable for all amputation levels and for the wrist disarticulation.

Article number	Upper arm connection Ø	For hand size	For lamination ring Ø	Interior joint Ø
12K27=58x45	70 mm	6 3/4 – 7 1/4	45 mm	58 mm
12K27=58x50	70 mm	7 3/4 – 8	50 mm	58 mm
12K27=78x45	90 mm	6 3/4 – 7 1/4	45 mm	78 mm
12K27=78x50	90 mm	7 3/4 – 8	50 mm	78 mm



647G470

Spare parts for 12K27

16Y26 Pull Cable

Article number	for	Complete with
16Y26	16X12 and 16X13	screw connection



21A18 Nylon wire

Article number	Ø	Length	Order by
21A18=2X1	2 mm	1 m	lfm
21A18=2X5	2 mm	5 m	lfm
21A18=2X10	2 mm	10 m	lfm
21A18=2X25	2 mm	25 m	lfm



1

Elbow Component with Cable Lock and Automatic Forearm Balance (AFB)



647G470

12K33 Elbow Component

With single-sided, rotatable cable lock (18 locking positions in about 7.2° increments), automatic forearm balance (AFB) and upper arm rotation joint (sickle joint).

Plastic forearm shell, skin colour, length 280 mm, circumference approx. 260 mm, skin-coloured plastic elbow ball.

Thanks to the exterior joint construction and innovative technology for interior and exterior rotation of the forearm, this joint allows the distal residual limb to reach up to the joint ball. It is therefore suitable for all amputation levels and for elbow disarticulation.

Article number	Upper arm connection Ø	For hand size	For lamination ring Ø	Interior joint Ø
12K33=58x45	70 mm	6 3/4 - 7 1/4	45 mm	58 mm
12K33=58x50	70 mm	7 3/4 - 8	50 mm	58 mm

2

3

Spare parts for 12K33



16Y26 Pull Cable

Article number	for	Complete with
16Y26	16X12 and 16X13	screw connection

4

5



12K35 Automatic Forearm Balance (AFB)

Automatic Forearm Balance

6



13G50 Strap Clip

7

8

Elbow Component with passive lock

12K6 Elbow Component

With manual elbow lock (13 locking positions in approximately 8° increments).
Plastic forearm shell, skin colour, length approx. 260 mm, circumference approx. 245 mm, skin-coloured plastic elbow ball, with wood adapter.

Article number	Upper arm connection Ø	For hand size	For lamination ring Ø
12K6=45	64 mm	6 3/4 – 7 1/4	45 mm
12K6=50	64 mm	7 3/4 – 8	50 mm



 647G469

12K5 Elbow Component

With manual elbow lock (13 locking positions in approximately 8° increments) and upper arm rotation joint (sickle joint).

Plastic forearm shell, skin colour, length approx. 260 mm.
Circumference approx. 245 mm, elbow ball is made of skin coloured plastic.

Article number	Upper arm connection Ø	For hand size	For lamination ring Ø
12K5=45	67 mm	6 3/4 – 7 1/4	45 mm
12K5=50	67 mm	7 3/4 – 8	50 mm



 647G469

12K20 Elbow Component

With manual elbow lock (13 locking positions in approximately 8° increments) and upper arm rotation joint (sickle joint).
Plastic forearm shell, **dark brown**, length approx. 260 mm, circumference 245 mm. plastic elbow ball, **dark brown**.

Article number	Upper arm connection Ø	For hand size	For lamination ring Ø
12K20=45	67 mm	6 3/4 – 7 1/4	45 mm
12K20=50	67 mm	7 3/4 – 8	50 mm



 647G469

Spare parts for 12K5, 12K6 and 12K20

13G8=67 Lamination Ring

Article number	Outside Ø	for
13G8=67	67 mm	12K5 12K20 12R1 12R2 12R5



1



21A18 Nylon wire

Article number	Ø	Length	Order by
21A18=2X1	2 mm	1 m	lfm
21A18=2X5	2 mm	5 m	lfm
21A18=2X10	2 mm	10 m	lfm
21A18=2X25	2 mm	25 m	lfm

2

3

4

5

6

7

8

The page features a large grid of small dots for taking notes. On the right side, there is a vertical sidebar with a light beige background for the top three sections and a teal background for the bottom five sections. The numbers 1 through 8 are printed in white on this sidebar, corresponding to the sections of the dot grid.

1
2
3
4
5
6
7
8



Movo Shoulder Joints

1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



647G349

12S6 MovoShoulder Swing

This shoulder joint opens new opportunities for prosthetic fittings in the shoulder area. MovoShoulder Swing is ideal for prostheses with basic functions or in combination with high-tech components. The free swing of up to 40° reduces pressure from the prosthesis socket and allows even bilateral amputees to achieve natural, harmonious movements.. Locking at 30° anteversion and unlocking is controlled by specific upper body movements or with the sound hand. No additional control elements such as switches or a body harness are required. Abduction of up to 20° facilitates more comfortable movement patterns during many activities of daily life. This especially applies to activities done close to the body or while sitting. The joint weighs 242 g and has a total length of 230 mm.

Article number	Side
12S6=L	Left (L)
12S6=R	Right (R)


Technical Data

Reference Number	12S6
Weight	242 g
Overall length	230 mm



Practical recommendation:

- This joint can also be fitted in conjunction with the 12K100=* DynamicArm.



647G473

12S4 Ottobock Shoulder Joint

The shoulder joint has two axes with separate friction and two arm bars.

Article number	12S4
Weight	134 g
Ø-Lamination ring	43 mm
Length	220 mm



647G473

12S5 Ottobock Shoulder Joint

The shoulder joint features two axes, separate friction, a shoulder bracket and two upper arm bars.

Article number	Side
12S5=L	Left (L)
12S5=R	Right (R)

Technical Data

Reference Number	12S5
Weight	149 g
Overall length	220 mm
Ø-Lamination ring	43 mm

Spare parts for 12S4 and 12S6

13Z16=43 Lamination Ring



13D1 O-Ring set

Technical Data

Article number	13D1
for	12S6

13D2 Shoulder Bracket Set

Article number	Side	for
13D2=L	Left	12S6=L
13D2=R	Right	12S6=R

1

2

3

4

5

6

7

8



1

2

3

4

5

6

7

8

Movo Terminal Devices

1

2

3

4

5

6

7

8



 647G444

Ottobock System Hand -passive-

With threaded stud and system inner hand,

Hand sizes:

6 3/4

for children

7 1/4

for adolescents and women

7 3/4

for men

8

for men

The passive system hand can be used for all lengths of the residual limb to be fitted with a cosmetic prosthesis. It is opened with the sound hand and closes independently. The system is lightweight and stable.

- The Prosthetic glove must be ordered separately. See pages 172-174, 214-216

Article number	Side	Threaded stud	Size	Inner hand
8K18=L6 3/4	Left (L)	M12x1,5	6 3/4	8X14=L6 3/4
8K18=L7 1/4	Left (L)	M12x1,5	7 1/4	8X14=L7 1/4
8K18=L7 3/4	Left (L)	M12x1,5	7 3/4	8X14=L7 3/4
8K18=L8	Left (L)	M12x1,5	8	8X14=L8
8K18=R6 3/4	Right (R)	M12x1,5	6 3/4	8X14=R6 3/4
8K18=R7 1/4	Right (R)	M12x1,5	7 1/4	8X14=R7 1/4
8K18=R7 3/4	Right (R)	M12x1,5	7 3/4	8K21=L8
8K18=R8	Right (R)	M12x1,5	8	8X14=R8

Technical Data

Reference Number	8K18	8K18	8K18	8K18
Size	6 3/4	7 1/4	7 3/4	8
Weight	ca. 185 g	ca. 250 g	ca. 280 g	ca. 290 g

Article number	Side	Threaded stud	Size	Inner hand
8K19=L6 3/4	Left (L)	1/2"-20	6 3/4	8X14=L6 3/4
8K19=L7 1/4	Left (L)	1/2"-20	7 1/4	8X14=L7 1/4
8K19=L7 3/4	Left (L)	1/2"-20	7 3/4	8X14=L6 3/4
8K19=L8	Left (L)	1/2"-20	8	8X14=L8
8K19=R6 3/4	Right (R)	1/2"-20	6 3/4	8X14=R6 3/4
8K19=R7 1/4	Right (R)	1/2"-20	7 1/4	8X14=R7 1/4
8K19=R7 3/4	Right (R)	1/2"-20	7 3/4	8X14=R7 3/4
8K19=R8	Right (R)	1/2"-20	8	8X14=R8

Technical Data

Reference Number	8K19	8K19	8K19	8K19
Size	6 3/4	7 1/4	7 3/4	8
Weight	ca. 185 g	ca. 250 g	ca. 280 g	ca. 290 g

Spare parts for passive system hands

8X14 System Inner Hand

Article number	For hand size	Side
8X14=L6 3/4	6 3/4	Left (L)
8X14=L7 1/4	7 1/4	Left (L)
8X14=L7 3/4	7 3/4	Left (L)
8X14=L8	8	Left (L)
8X14=R6 3/4	6 3/4	Right (R)
8X14=R7 1/4	7 1/4	Right (R)
8X14=R7 3/4	7 3/4	Right (R)
8X14=R8	8	Right (R)

Technical Data

Reference Number	8X14
Complete with	finger support in the little and ring finger and 9S187 Retainer Ring



9S15 Retainer Ring, narrow

Article number	For hand size	Ø
9S15=42	6 3/4	42



9S187 Retainer Ring, wide

Article number	For hand size
9S187=7 1/4	7 1/4
9S187=7 3/4	7 3/4, 8

Technical Data

Article number	9S187=7 1/4	9S187=7 3/4
for	8X14 System Inner Hands	8X14 System Inner Hands



9S6 Finger and Thumb Tip

Article number	For hand size
9S6	6 3/4, 7 1/4, 7 3/4 and 8



1

2

3

4

5

6

7

8

1



9S67 Chassis

With threaded stud

Article number	Side	Ø	Threaded stud	for
9S67=L40	Left (L)	40 mm	M12x1,5	8K18=L6 3/4
9S67=L44-N	Left (L)	44 mm	M12x1,5	8K18=L7 1/4
9S67=L48-N	Left (L)	48 mm	M12x1,5	8K18=L7 3/4 8K18=L8
9S67=R40	Right (R)	40 mm	M12x1,5	8K18=R6 3/4
9S67=R44-N	Right (R)	44 mm	M12x1,5	8K18=R7 1/4
9S67=R48-N	Right (R)	48 mm	M12x1,5	8K18=R7 3/4 8K18=R8

2

3



9S65 Chassis

With threaded stud

Article number	Side	Ø	Threaded stud	for
9S65=L40	Left (L)	40 mm	1/2"-20	8K19=L6 3/4
9S65=L44-N	Left (L)	44 mm	1/2"-20	8K19=L7 1/4
9S65=L48-N	Left (L)	48 mm	1/2"-20	8K19=L7 3/4 8K19=L8
9S65=R40	Right (R)	40 mm	1/2"-20	8K19=R6 3/4
9S65=R44-N	Right (R)	44 mm	1/2"-20	8K19=R7 1/4
9S65=R48-N	Right (R)	48 mm	1/2"-20	8K19=R7 3/4 8K19=L8

4



5

6



501T52=M3x6 Flange Button Head Socket Screw

7



501T52=M3x8 Flange Button Head Socket Screw

8

The page features a large grid of dots for taking notes. On the right side, there is a vertical sidebar with eight numbered tabs (1-8). The tabs are arranged vertically, with tabs 1-4 in a light brown color and tabs 5-8 in a teal color. The grid of dots covers the majority of the page area, starting below the header and ending above the footer.

1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



Movo Prosthetic Gloves

1

2

3

4

5

6

7

8

1



8S6=170x65 Prosthetic Glove for Children

For hand sizes 6 3/4, with long sleeve.
Standard prosthetic glove available in 18 different shades

Article number	For hand size	Side	Sleeve length
8S6=170x65L	6 3/4	Left (L)	300 mm
8S6=170x65R	6 3/4	Right (R)	300 mm

When ordering, please include the colour code

2

3

- In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.
- To order them, please enter an **N** before the =:
example 8S11=210x78xL4 standard glove
8S11**N**=210x78xL4 MyoSkin Natural
Please ask our customer service for the available colours.

The entire range of colours for Skin Natural models comprises of six colours. The new 646M47 colour scale helps when selecting the right pattern.

4



8S4=190x76 Prosthetic Glove for Adolescents

For hand sizes 7 1/4, with short sleeve.
Standard prosthetic glove available in 18 different shades

Article number	For hand size	Side	Sleeve length
8S4=190x76L	7 1/4	Left (L)	215 mm
8S4=190x76R	7 1/4	Right (R)	215 mm

When ordering, please include the colour code

5

6

- In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.
- To order them, please enter an **N** before the =:
example 8S11=210x78xL4 standard glove
8S11**N**=210x78xL4 MyoSkin Natural
Please ask our customer service for the available colours.

The entire range of colours for Skin Natural models comprises of six colours. The new 646M47 colour scale helps when selecting the right pattern.

7

8

8S4=210x78 Prosthetic Glove for Men

For hand sizes 7 3/4, with short sleeve.
Standard prosthetic glove available in 18 different shades

Article number	For hand size	Side	Sleeve length
8S4=210x78L	7 3/4	Left (L)	220 mm
8S4=210x78R	7 3/4	Right (R)	220 mm

When ordering, please include the colour code

- In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.
- To order them, please enter an **N** before the =:
example 8S11=210x78xL4 standard glove
8S11**N**=210x78xL4 MyoSkin Natural
Please ask our customer service for the available colours.

The entire range of colours for Skin Natural models comprises of six colours. The new 646M47 colour scale helps when selecting the right pattern.



1

2

3

8S4=220x80 Prosthetic Glove for Men


For hand sizes 8, with short sleeve.
Standard prosthetic glove available in 18 different shades

Article number	For hand size	Side	Sleeve length
8S4=220x80L	8	Left (L)	225 mm
8S4=220x80R	8	Right (R)	225 mm

When ordering, please include the colour code

- In addition to the standard prosthetic glove, Ottobock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.
- To order them, please enter an **N** before the =:
example 8S11=210x78xL4 standard glove
8S11**N**=210x78xL4 MyoSkin Natural
Please ask our customer service for the available colours.

The entire range of colours for Skin Natural models comprises of six colours. The new 646M47 colour scale helps when selecting the right pattern.



4

5

6

7

8

1



8S5=195x78 Prosthetic Glove for Women

For hand sizes 7 1/4, with long sleeve.

Standard prosthetic glove available in 18 different shades

Article number	For hand size	Side	Sleeve length
8S5=195x78L	7 1/4	Left (L)	340 mm
8S5=195x78R	7 1/4	Right (R)	340 mm

When ordering, please include the colour code

2

- In addition to the standard prosthetic glove, Otto Bock also offers additional models of the **Skin Natural** series. The multilayer structure of the gloves gives them a depth effect so that they have a very vivid appearance. The outer translucent layer allows the coloured fibres of the material to shine through. This simulates the natural vein structure of human skin.

To order them, please enter an **N** before the =:

example 8S11=210x78xL4 standard glove

8S11**N**=210x78xL4 MyoSkin Natural

Please ask our customer service for the available colours.

3

The entire range of colours for Skin Natural models is covered with six colours. The new 646M47 colour scale helps when selecting the right pattern.

4

5

6

7

8

The page features a large grid of dots for taking notes. On the right side, there is a vertical sidebar with eight numbered tabs (1-8). The tabs are arranged vertically, with tabs 1-4 in a light brown color and tabs 5-8 in a teal color. The grid of dots covers the majority of the page area.

1

2

3

4

5

6

7

8

1

2

3

4

5

6

7

8



Movo Wrist Units

1

2

3

4

5

6

7

8

Adapter

1



10R1 Adapter

For connecting an inner hand or prosthetic hand to the modular arm components

Article number	Complete with
10R1	With threaded stud M12x1.5

647G471

2

3



10R2 Adapter

Light metal

For connecting a hook or system hand to modular arm components

Article number	Interior thread
10R2=M12x1.5	M12x1,5
10R2=1/2"-20	1/2"-20

Technical Data

Reference Number	10R2
Connection-Ø	20 mm
Complete with	inner thread

647G471

4

5



10R3 Adapter with Flexion

Light metal

For connecting a hook or system hand to modular arm components

Article number	Interior thread
10R3=M12x1.5	M12x1,5
10R3=1/2"-20	1/2"-20

Technical Data

Reference Number	10R3
Connection-Ø	20 mm
Complete with	inner thread

647G471

6

7

8

Modules

10R4 Knurled Plate


For connecting an inner hand or prosthetic hand to the 10R2 and 10R3 adapter

Article number	Threaded stud
10R4=M12x1.5	M12x1,5 / M12x1,5
10R4=1/2"-20	M12x1,5 / 1/2"-20

Technical Data

Article number	10R4=M12x1.5	10R4=1/2"-20
for	10R2=M12x1.5 10R3=M12x1.5	10R2=1/2"-20 10R3=1/2"-20
Complete with	Threaded stud on both ends	




 647G471

13R11 Modular Adapter for ArmLiners

For connecting Ottobock 14Y1 ArmLiner and 14A1 Lock Set

Article number	Complete with
13R11	Ottobock Modular Arm Component



 647G164

1

2

3

4

5

6

7

8

1

2

3

4

5

6

7

8



Movo Elbow Components

1

2

3

4

5

6

7

8

Modular arm components for above-elbow residual limb

1



12R1 Ottobock Modular Arm Component

separate friction

Article number	Side
12R1=L	Left (L)
12R1=R	Right (R)

2



Technical Data

Reference Number	12R1
for	above-elbow residual limb,
Complete with	passive elbow lock, upper arm and forearm rotation

3

647G471

4

Spare parts for 12R1

5



13R1 Modular Elbow Joint

Article number	Side
13R1=L	Left (L)
13R1=R	Right (R)

Technical Data

Reference Number	13R1
Complete with	passive lock

6

647G471

7



13G8=67 Lamination Ring

Article number	Outside Ø	for
13G8=67	67 mm	12K5 12K20 12R1 12R2 12R5

8

13R3 Light Metal Tube

Gold anodised

Article number	Length	Complete with
13R3	250 mm	reinforced ring and 11D58 Glide Bushing



1

Accessories for 12R1

15K1=30 Pre-Shaped Foam Block

skin coloured Perlon Prosthetic Stockinette

Article number	15K1=30
Colour	skin-coloured
Forearm circumference	300 mm
Complete with	proximally attached 13R8 Connection Flange and 1 m 623T8=9



 647G471

2

3

4

13Y1 Flexion Cable Retainer

Technical Data

Article number	13Y1
for	forearm tube



5

6

13R9 Connection Ring with groove

Plastic

For attaching the perlon stockinette to the wrist unit

Article number	For hand size
13R9=45	6 3/4
13R9=50	7 1/4
13R9=55	7 3/4, 8
13R9=65	prosthetic hands (adaptable)



7

Technical Data

Article number	13R9=45	13R9=50	13R9=55	13R9=65
Outer Ø	45 mm	50 mm	55 mm	65 mm

8

1



 647G471

12R2 Ottobock Modular Arm Component

separate friction

Article number	Side
12R2=L	Left (L)
12R2=R	Right (R)

Technical Data

Reference Number	12R2
for	above-elbow residual limb
Complete with	cable lock suitable for left and right side (18 locking positions in about 7.2° increments)

2

3



12R5 Ottobock Modular Arm Component

separate friction

Article number	Side
12R5=L	Left (L)
12R5=R	Right (R)

Technical Data

Reference Number	12R5
for	above-elbow residual limb
Complete with	cable lock suitable for left and right side (18 locking positions in about 7.2° increments)

4

5

6

7


8

Spare parts for 12R2 and 12R5

13R1 Modular Elbow Joint

Article number	Side	Complete with
13R1=L	Left (L)	passive lock
13R1=R	Right (R)	passive lock



 647G471

13G8=67 Lamination Ring

Article number	Outside Ø	for
13G8=67	67 mm	12K5 12K20 12R1 12R2 12R5



13R3 Light Metal Tube

Gold anodised

Article number	Length	Complete with
13R3	250 mm	enforced ring and 11D58 Glide Bushing



16Y27 Joint Piece with Cable Lock

(18 ratchet settings of approximately 7.2°),
suitable for both sides

Article number	for	Complete with
16Y27	16X12 and 16X13	joint plate and flat head screws



Accessories for 12R2 and 12R5

1



15K2=30 Pre-Shaped Foam Block

preformed

Article number	15K2=30
Colour	skin-coloured
Forearm circumference	300 mm
Complete with	separate 13R8 Connection Flange and 1.30 m 623T8=9 skin coloured Perlon Stockinette

2

3



13Y1 Flexion Cable Retainer

Technical Data

Article number	13Y1
for	forearm tube

4



13R9 Connection Ring with groove

Plastic

For attaching the perlon stockinette to the wrist unit

Article number	For hand size
13R9=45	6 3/4
13R9=50	7 1/4
13R9=55	7 3/4, 8
13R9=65	prosthetic hands (adaptable)

5


Technical Data

Article number	13R9=45	13R9=50	13R9=55	13R9=65
Outer Ø	45 mm	50 mm	55 mm	65 mm

6

Modular arm components for shoulder disarticulation

7



12R3 Ottobock Modular Arm Component

For connecting a shoulder joint

Article number	Side
12R3=L	Left (L)
12R3=R	Right (R)

647G471

8

Reference Number	12R3
for	shoulder disarticulation
Complete with	passive elbow lock and 13R7 adapter

Accessories for 12R3

15K3=30 Pre-Shaped Foam Block

preformed, two-piece

Article number	15K3=30
Colour	skin-coloured
Forearm circumference	300 mm
Complete with	approx. 2 m 623T8=9 skin coloured Perlon Stockinette



1


2

13R9 Connection Ring with groove

Plastic

For attaching the perlon stockinette to the wrist unit

Article number	For hand size
13R9=45	6 3/4
13R9=50	7 1/4
13R9=55	7 3/4, 8
13R9=65	prosthetic hands (adaptable)



3

4

Technical Data

Article number	13R9=45	13R9=50	13R9=55	13R9=65
Outer Ø	45 mm	50 mm	55 mm	65 mm

13Y1 Flexion Cable Retainer

Technical Data

Article number	13Y1
for	forearm tube




5


6

12S4 Ottobock Shoulder Joint

The shoulder joint has two axes with separate friction and two arm bars.

Article number	12S4
Weight	134 g
Ø-Lamination ring	43 mm
Length	220 mm



 647G473

7

8

1



12S5 Ottobock Shoulder Joint

The shoulder joint features two axes, separate friction, a shoulder bracket and two upper arm bars.

Article number	Side
12S5=L	Left (L)
12S5=R	Right (R)

2

 647G473

Technical Data

Reference Number	12S5
Weight	149 g
Overall length	220 mm
Ø-Lamination ring	43 mm

3

4

5

6

7

8

12R4 Ottobock Modular Arm Component

For connecting a 12S7 shoulder joint

Article number	Side
12R4=L	Left (L)
12R4=R	Right (R)

Technical Data

Reference Number	12R4
for	shoulder disarticulation
Complete with	passive elbow lock



 647G471

Accessories for 12R4

15K3=30 Pre-Shaped Foam Block

preformed, two-piece

Article number	15K3=30
Colour	skin-coloured
Forearm circumference	300 mm
Complete with	approx. 2 m 623T8=9 skin coloured Perlon Stockinette



13R9 Connection Ring with groove

Plastic

For attaching the perlon stockinette to the wrist unit

Article number	For hand size
13R9=45	6 3/4
13R9=50	7 1/4
13R9=55	7 3/4, 8
13R9=65	prosthetic hands (adaptable)



Technical Data

Article number	13R9=45	13R9=50	13R9=55	13R9=65
Outer Ø	45 mm	50 mm	55 mm	65 mm

1

2

3

4


5

6

7

8

1



12S7 Ottobock Ball Shoulder Joint

Article number	12S7
Material	Light metal
Weight	78 g
Complete with	threaded stud M12x1.5

647G473

2

Spare parts for 12S7

3



13X5 Joint Ball

4



13X4 Joint Clamp

5



13X3 Friction Ring

Article number	13X3
Material	plastic

6

7

647G454

Accessories for 12S7

10R2=M12x1.5 Adapter

for connecting modular arm components to 12S7 Ball Shoulder Joint

Article number	Interior thread	Complete with
10R2=M12x1.5	M12x1,5	internal thread M12x1.5

8

Modular Elbow Joint

13R1 Modular Elbow Joint

Article number	Side	Complete with
13R1=L	Left (L)	passive lock
13R1=R	Right (R)	passive lock



 647G471

Spare parts for 13R1

13G37=N Joint Protector



1

2

3

4

5

6

7

8

Modules

1



13R3 Light Metal Tube

Gold anodised

Article number	Length	Complete with
13R3	250 mm	reinforced ring and 11D58 Glide Bushing

2

3



13R4 Light Metal Tube

Gold anodised

Article number	Length	Complete with
13R4	120 mm	reinforced ring and 11D58 Glide Bushing

4

5



10R5 Adapter

For connecting the 13R3/13R4 Tubes and the 13R6/13R7 Adapters to the elbow joint

Article number	for
10R5	hand and upper arm rotation

6



13R7 Adapter

Article number	for
13R7	12S4 and 12S5 Shoulder Joints

7

8



13R5 Adapter

For connecting the 13G8 Lamination Ring to the aluminium alloy tube

13R6 Adapter

For connecting the 13G8 Lamination Ring to the elbow joint using the 10R5 Adapter



1

13G8=67 Lamination Ring

Article number	Outside Ø	for
13G8=67	67 mm	12K5 12K20 12R1 12R2 12R5




2

3

501S35=M4x6 Truss Head Screw

Article number	for
501S35=M4x6	13R5, 13R6 and 13R7



4

501Z2=M4x18 Cap Screw

Article number	for
501Z2=M4x18	13R5



5

6

506G1=M6x10 Set Screw

Article number	for
506G1=M6x10	10R5



7

8

1



13R8=64 Mounting Flange

for attaching the foam cover to the lamination ring

Article number	13R8=64
-----------------------	----------------

Material	Plastic
-----------------	---------

2



13Y1 Flexion Cable Retainer

Article number	for
-----------------------	------------

13Y1	forearm tube
-------------	--------------

3

4

5

6

7

8

The page features a large grid of dots for taking notes. On the right side, there is a vertical sidebar with eight numbered tabs (1-8). The tabs are arranged vertically, with tabs 1-4 in a light beige color and tabs 5-8 in a teal color. The numbers are white and centered on each tab.

1
2
3
4
5
6
7
8



Movo Arm Bars

1

2

3

4

5

6

7

8

Elbow Joint Bars

1




16X12 Elbow Joint Bars

One joint bar features rotatable cable lock (18 locking positions in 7.2° increments). Second joint bar without lock. Flat bar profile, suitable for both sides

Article number	Head Ø	Upper/lower bar length	Bar width / thickness
16X12	30 mm	320/320 mm	14/2 mm


2



Technical Data

Article number	16X12
Material	Stainless steel
for	orthoses and prostheses

3

 647G472

4




16X13 Elbow Joint Bar with Cable Lock

Joint features rotatable cable lock (18 locking positions in 7.2° increments). Flat bar profile, suitable for both sides

Article number	Head Ø	Upper/lower bar length	Bar width / thickness
16X13	30 mm	320/320 mm	14/2 mm

5

 647G472

Technical Data

Article number	16X13
Material	Stainless steel
for	orthoses and prostheses

6

7

8

16X14 Elbow Joint Bar without Lock


Free motion joint, flat bar profile

Article number	Head Ø	Upper/lower bar length	Bar width / thickness
16X14	30 mm	320/320 mm	14/2 mm

Technical Data

Article number	16X14
Material	Stainless steel



 647G472

Spare parts for 16X12, 16X13, 16X14

16Y27 Joint Piece with Cable Lock

(18 ratchet settings of approximately 7.2°),
suitable for both sides

Article number	16Y27
for	16X12 and 16X13
Complete with	joint plate and flat head screws



16Y31 Joint Piece without Lock

Article number	16Y31
for	16X12 and 16X14
Complete with	joint plate and Flat Head Screws




16Y26 Pull Cable

Article number	16Y26
for	16X12 and 16X13
Complete with	screw connection



1



21Y79 Bow with wedge lock

Article number	21Y79
for	16X12 and 16X13
Complete with	gib and cotter

2

3

4

5

6

7

8

Arm Bars


16U4 Double Joint Arm Bars

Free-motion joints, concave half-round bar profile

Article number	Upper/lower bar length	Bar width	Bar length
16U4	270 / 270 mm	14 mm	2 mm

Technical Data

Article number	16U4
Material	Stainless steel
for	orthoses and prostheses



16U5 Double Joint Arm Bars

Free motion joints, with step-up, concave half-round bar profile

Technical Data

Article number	16U5
Material	stainless steel
for	short below-elbow residual limb
Upper/lower bar length	270 / 270 mm
Bar width	14 mm
Bar thickness	2 mm
Length of the step-up bar	105 mm



1

2

3

4

5

6

7

8

1



16X4 Arm Bars

Joints with automatic ratchet lock, locks only in extension, releases at a max. flexion of approx. 135°, concave half-round bar profile


Article number	Head Ø	Upper/lower bar length	Bar width	Bar thickness
16X4	26 mm	210/230 mm	14 mm	2 mm

2

Technical Data

Article number	16X4
Material	Stainless steel

3



4



16U8 Arm Bar

Article number	Width	Material thickness
16U8	14 mm	2.45 mm

5

Technical Data

Article number	16U8
Material	Stainless steel
Overall length	300 mm
Package content	2 pieces pc.

6



Practical recommendation:

- 16U7 and 16U8 together form a bilateral joint bar.

7

8

16U7 Arm Bar with Joint

Free motion joint in flexion direction

Article number	Width	Material thickness	Length from joint centre
16U7	14 mm	2.45 mm	307 mm

Technical Data

Article number	16U7
Material	Stainless steel
Package content	2 pieces pcs.



Practical recommendation:

- 16U6 and 16U7 together form an above-elbow bar with two-part joint. To obtain the required length, the ends of the 16U6 are cut.
- 16U7 and 16U8 together form a bilateral joint bar.

16U6 Arm Frame

Article number	Width	Material thickness
16U6	14 mm	2.45 mm

Technical Data

Article number	16U6
Material	Stainless steel
Overall length	600 mm



Practical recommendation:

- The 16U6 Arm Frame can be combined with the 16U7 Arm Bar with Elbow Joint to form an above-elbow bar with a two-part joint. To obtain the required length, the ends of the 16U6 can be cut.

17F48=5 Joint



1



17Y17=6x8x2.4 Bushing

Article number	17Y17=6x8x2.4
Material	Brass

2



17Y93=6x5.2xM4 Bearing Nut

3



4



506A12 Rivet Pin

Article number	506A12=4x22	506A12=4x30
Material	Stainless steel	Stainless steel

5



Phillips Countersink Screw Head

Article number	Thread	Length
501T19=M5x6	M5	6 mm

6



16Y9 Joint Plate

Article number	16Y9
Material	Stainless steel

7

Spare parts for 16X4

8



16Y7 Locking Latch

Article number	16Y7
Material	Stainless steel

16Y8 Lock Washer

Article number	16Y8
Material	Stainless steel



1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



Movo Shoulder Joints

1

2

3

4

5

6

7

8

MovoShoulder Swing

1
2
3
4
5
6
7
8



 647G349

12S6 MovoShoulder Swing

This shoulder joint opens new opportunities for prosthetic fittings in the shoulder area. MovoShoulder Swing is ideal for prostheses with basic functions or in combination with high-tech components. The free swing of up to 40° reduces pressure from the prosthesis socket and allows even bilateral amputees to achieve natural, harmonious movements.. Locking at 30° anteversion and unlocking is controlled by specific upper body movements or with the sound hand. No additional control elements such as switches or a body harness are required. Abduction of up to 20° facilitates more comfortable movement patterns during many activities of daily life. This especially applies to activities done close to the body or while sitting. The joint weighs 242 g and has a total length of 230 mm.

Article number	Side
12S6=L	Left (L)
12S6=R	Right (R)

Technical Data

Reference Number	12S6
Weight	242 g
Overall length	230 mm



Practical recommendation:

- This joint can also be fitted in conjunction with the 12K100=* DynamicArm.

Shoulder Joints


12S4 Ottobock Shoulder Joint

The shoulder joint has two axes with separate friction and two arm bars.

Article number	Length
12S4	220 mm

Technical Data

Article number	12S4
Weight	134 g
Ø-Lamination ring	43 mm



 647G473

12S5 Ottobock Shoulder Joint


The shoulder joint features two axes, separate friction, a shoulder bracket and two upper arm bars.

Article number	Side
12S5=L	Left (L)
12S5=R	Right (R)

Technical Data

Reference Number	12S5
Weight	149 g
Overall length	220 mm
Ø-Lamination ring	43 mm



 647G473

Spare parts for 12S4 and 12S5

13Z16=43 Lamination Ring



1



647G473

12S7 Ottobock Ball Shoulder Joint

Article number	12S7
Material	Light metal
Weight	78 g
Complete with	threaded stud M12x1.5

2

Accessories for 12S7

647G454

10R2=M12x1.5 Adapter


for connecting modular arm components to 12S7 Ball Shoulder Joint

Article number	Interior thread	Complete with
10R2=M12x1.5	M12x1,5	internal thread M12x1.5

3

4

Spare parts for 12S7



13X5 Joint Ball

5



13X3 Friction Ring

Article number	13X3
Material	plastic

6



13X4 Joint Clamp

7

8

The page features a large grid of dots for taking notes. On the right side, there is a vertical sidebar with eight numbered tabs (1-8). The tabs are arranged vertically, with tabs 1-4 in a light brown color and tabs 5-8 in a teal color. The grid of dots covers the majority of the page area, starting below a horizontal line and extending to the bottom of the page.

1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



Physo

Physo Terminal Devices

1

2

3

4

5

6

7

8

1

The following pages, containing photographs of every passive prosthetic hand available from Otto Bock, are intended to assist you and your patients in selecting the most appropriate prosthetic hand.




2

Otto Bock passive prosthetic hands consist of an inner hand and a prosthetic glove. The standard inner hand is molded from a light-weight foam material with steel wire reinforcement in each finger. Various adapters are available to provide connection to the prosthesis.

3





Shape, color and detailed surface structure of the prosthetic glove give a lifelike appearance to the inner hand. The 43 different models for children, women and men are illustrated for the right hand and are marked for identification as in the table below:

4

		Prosthetic glove	Inner hand
	Children	8S6=	8S9=
	Women	8S5=	8S8=
	Men	8S4=	8S7=

5

The order numbers appear on each page above the passive prosthetic hand. The prosthetic glove (line 1) and inner hand (line 2) each have unique numbers for the left and right side.

	8S5=165x72L	40	8S5=167x72R	
	8S8=165x72L		8S7=167x72R	


6

Composing the order number

To precisely determine the correct size of the prosthetic glove and inner hand, take actual measurements of the „metacarpal circumference“ (B) and „middle finger length“ (C) of your patient guided by the following sketch.

Comparison of the actual „metacarpal circumference“ (B) and „middle finger length“ (C) measurements from your patient with those indicated in the tables on pages 3.3 and 3.4 will help to ensure proper size selection.

7



8

Metacarpal circumf. (B)	Middle finger length (C)	Wrist circumf. (A)	Below-elbow circumf. (D)	Length to elbow (E)
160	68	140	200	384

These measurements lead for example to the following order number for the prosthetic glove and the corresponding inner hand. Actual measurements of all inner hands differ from those shown, as they are reduced by the wall thickness of the cosmetic glove.

8S5=	165 x 72	L
8S8=		

The final step in the selection process is to select the color shade, i.e. the color number. The standard gloves are available in 18 shades which are reproduced on page 3.5 as accurately as possible, together with their corresponding code number. Color shade guides on actual glove material will be sent free of charge on request under order number 646M3. After selecting the color, the order number for prosthetic glove and inner hand should appear as in the following example:

8S5=	165 x 72	L	6
8S8=	165 x 72	L	

If the inner hand is to be used for a partial hand replacement a „T“ must be added to the order number as follows:

8S8=	165 x 72	L	T
------	----------	---	---

All measurements are stated in mm and are subject to a tolerance of ± 5 mm.

Prosthetic Gloves and Inner Hands for Children

Article Number for Ordering		Measurements for the left hand					Measurements for the right hand					
Prosthetic glove	Inner hand	Article No.	metacarpal circumf. (B)	middle finger length (C)	wrist circumf. (A)	below-elbow circumf. (D)	length to elbow (E)	metacarpal circumf. (B)	middle finger length (C)	wrist circumf. (A)	below-elbow circumf. (D)	length to elbow (E)
8S6=		115 x 38L	110	38	107	147	211	113	37	107	144	215
		130 x 51L	126	46	112	150	262	133	47	113	156	263
		142 x 50L	140	46	124	148	294	141	49	125	158	277
		151 x 58L	140	54	123	164	294	142	54	126	168	298
		158 x 54L	158	54	133	188	341	159	53	134	185	352
	165 x 68L	163	59	135	181	380	160	63	140	181	355	
	168 x 70L	164	64	145	180	376	161	63	142	176	369	
	170 x 65L	163	57	146	182	397	165	61	142	180	370	

Prosthetic Gloves and Inner Hands for Women

Article Number for Ordering		Measurements for the left hand					Measurements for the right hand					
Prosthetic glove	Inner hand	Article No.	metacarpal circumf. (B)	middle finger length (C)	wrist circumf. (A)	below-elbow circumf. (D)	length to elbow (E)	metacarpal circumf. (B)	middle finger length (C)	wrist circumf. (A)	below-elbow circumf. (D)	length to elbow (E)
8S5=		165 x 72L	160	68	140	200	384	165	66	149	204	394
		174 x 74L	172	67	151	190	408	170	67	148	185	410
		180 x 80L	167	72	159	225	488	173	71	156	215	470
		182 x 84L	179	73	154	204	290	187	74	157	210	285
		184 x 75L	180	66	149	226	461	181	69	153	221	407
		184 x 78L	182	70	155	221	400	184	73	153	215	426
		185 x 75L	180	70	152	202	402	180	69	157	205	405
		188 x 79L	185	77	163	225	433	185	76	160	235	433
		190 x 77L	188	66	148	210	444	186	66	152	215	438
		190 x 93L	185	82	182	222	467	184	88	178	218	464
		192 x 78L	186	69	154	223	403	189	72	163	224	395
		194 x 82L	188	73	157	215	468	185	73	156	214	488
		195 x 78L	195	69	165	214	460	194	73	167	223	437
	195 x 79L	190	70	166	228	408	196	74	168	235	405	
	208 x 89L	197	77	175	234	456	201	79	180	235	458	

Prosthetic Gloves and Inner Hands for Men

Article Number for Ordering		Measurements for the left hand					Check no.	Measurements for the right hand					
Prosthetic glove	Inner hand	Article No.	metacarpal circumf. (B)	middle finger length (C)	wrist circumf. (A)	below-elbow circumf. (D)		length to elbow (E)	metacarpal circumf. (B)	middle finger length (C)	wrist circumf. (A)	below-elbow circumf. (D)	length to elbow (E)
		202 x 74L	195	73	160	226	480	206 x 76R	200	73	167	212	468
		203 x 83L	200	70	166	233	420	199 x 82R	192	72	160	227	422
		203 x 85L	195	76	162	232	481	206 x 85R	196	76	166	236	465
		205 x 81L	202	72	171	180	196	205 x 81R	199	76	175	180	190
		206 x 80L	206	80	184	250	483	205 x 80R	205	80	185	255	474
		206 x 87L	200	80	160	220	420	215 x 88R	207	81	166	240	468
		207 x 86L	197	80	166	227	483	209 x 86R	205	82	168	231	471
		208 x 85L	205	75	176	236	475	212 x 83R	208	76	186	240	475
		211 x 88L	202	81	170	243	460	212 x 86R	215	81	174	246	463
		212 x 93L	201	81	175	245	510	215 x 93R	206	85	175	239	502
		213 x 85L	213	73	187	243	462	218 x 85R	211	76	187	245	460
		214 x 82L	210	74	179	250	450	215 x 83R	210	71	185	246	444
		218 x 85L	218	75	180	233	498	218 x 83R	214	75	180	240	470
		220 x 91L	210	77	189	236	453	214 x 90R	208	78	181	230	453
		221 x 81L	218	68	182	245	468	225 x 82R	225	76	183	239	460
		228 x 84L	224	73	180	245	445	222 x 84R	218	79	173	244	444
		228 x 88L	223	78	186	250	345	228 x 89R	220	86	194	247	419
		232 x 94L	228	81	184	270	533	230 x 93R	225	80	184	250	540
		238 x 92L	223	76	187	255	465	244 x 94R	223	76	189	255	455

8S4= 8S7=


Examples for ordering:

8S4= 202 x 74 L 6 bzw. 202 x 74 L


If the inner hand is to be used for partial hand replacement a „T“ must be added to the article number, for example:

8S7= 202 x 74 L T





Prosthetic gloves in the colors 2, 4, 6, 8, 10 and 14 are available as stock. For all other colors, please note that the gloves are fabricated only after receipt of an order and that, for this reason, delivery of such gloves will take longer.



1

2

3

4


5


6



7

8





1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6




7

8





- 1
- 2
- 3
- 4
- 5
- 6
- 7
- 8



1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6




7

8





- 1
- 2
- 3
- 4
- 5
- 6
- 7
- 8



1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6




7

8





- 1
- 2
- 3
- 4
- 5
- 6
- 7
- 8



1

2

3

4


5

6




7

8





- 1
- 2
- 3
- 4
- 5
- 6
- 7
- 8



1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8





1


2

3

4

5

6



7

8

1

2

3

4


5

6



7

8




8K5 Physolino Babyhand

For babies and infants, made of medical silicone; hand circumference approx. 110 mm, finger length approx. 40 mm; the colour approximately corresponds to colour no. 4, 11, or 16 of Ottobock's colour scale for prosthetic gloves, with the 11D31 lamination ring, Ø 24 mm

Article number	Side	Approximately corresponds to colour no.
8K5=1L1	Left (L)	1
8K5=1L4	Left (L)	4
8K5=1L11	Left (L)	11
8K5=1L16	Left (L)	16
8K5=1R1	Right (R)	1
8K5=1R4	Right (R)	4
8K5=1R11	Right (R)	11
8K5=1R16	Right (R)	16

- 743Y42 Lamination Dummy for Physolino Babyhand
For matching prosthetic gloves, see Page 306



 647G360

1

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



Physo Wrist Units

1

2

3

4

5

6

7

8

1



11D31 Lamination Ring for Physolino Baby Hand

Article number	Ø	for
11D31	24 mm	socket replacements

2



10A40 Wood Hand Adapter

for connecting the passive inner hand with the forearm socket

Article number	Complete with
10A40	threaded stud made of plastic, M12x1.5

3

4



10A41 Threaded Stud

For connecting inner hand to wrist units without ratchet

Article number	Complete with
10A41	threaded stud M12x1.5

5

6

7

8

A large grid of dotted lines for taking notes, spanning most of the page width and height.

1
2
3
4
5
6
7
8

1
2
3
4
5
6
7
8



Liner

1

2

3

4

5

6

7

8

1

ArmLiner

2



 647H323

14Y1 Silicon ArmLiner

Residual limb encasement is of the utmost importance for the quality and comfort of an upper limb prosthesis. Using Ottobock ArmLiners, which were developed especially for the needs of arm prosthetics, clearly improves wear comfort, ensures good residual limb adhesion and reduces frictional forces. With an adequate residual limb length, it is possible to do without an elbow envelopment. This enables unrestricted pronation and supination.

For the fitting of below-elbow and above-elbow residual limbs.

3

Silicone offers particular advantages:

- Temperature resistant from -60 °C to +200 °C
- Resistant to water, sweat and weather
- Neutral smell and taste
- High gas permeability
- Good adhesion

4

Article number	Length	Scope 1	Scope 2
14Y1=110	200 mm	110 mm	150 mm
14Y1=140	200 mm	140 mm	160 mm
14Y1=160	200 mm	160 mm	180 mm
14Y1=180	200 mm	180 mm	200 mm
14Y1=200	200 mm	200 mm	220 mm
14Y1=220	200 mm	220 mm	240 mm

5

- To select the right ArmLiner, the residual limb circumference is measured 3 cm before the distal and proximal end of the residual limb

6



14A1 Lock Set

For attaching a 14Y1 Ottobock Silicone ArmLiner to the prosthetic socket

7

Spare parts for 14A1

8



14A110 Lock with Release Pin

14A107 Pin

Article number	Length
14A107	23 mm



1

Accessoires for 14Y1

453H10 Derma Clean

- For the skin, prosthesis, orthosis and liner
- Cleans gently and reliably
- Ph-neutral, free of alkali and phosphates
- Anti-bacterial formula
- 453H10=1-N (1 bottle)

Article number	Order by	Content
453H10	6 bottles	300 ml



2

3

453H12 Derma Prevent

- Prevents chafing
- Inhibits contact with external allergens
- Covers the skin with a protective coating and leaves it soft and supple
- Inhibits perspiration and odour formation through the individual release of an active substance
- 453H12=1 (1 bottle)

Article number	Order by	Content
453H12	6 bottles	100 ml



4

5

453H14 Derma Repair

- Moisturises and promotes the regeneration of dry, irritated skin
- Reduces the effects of excessive strain and soothes irritated skin
- Antibacterial formula, restores the skin's defences against harmful environmental influences
- Regulates moisture and makes the skin noticeably more supple and elastic
- Enhances skin function, promotes blood circulation in the skin and supports cell growth
- Contains panthenol and vitamin E
- 453H14=1 (1 bottle)

Article number	Order by	Content
453H14	6 bottles	200 ml



6

7

8

Accessoires for 14A1

1



14A107 Pin

Article number	Length
14A107=1	28 mm
14A107=2	200 mm

2

3



14A111 Dummy Set

For aligning a prosthesis with Ottobock Silicone ArmLiner

Article number	Consists of:
14A111	Pin dummy with thread Pin dummy without thread Shaped dummy for lamination ring

4



640F18 Spray for Donning Silicone Liners

Article number	Content
640F18	45 ml

5

6

640F18=900 Refill Package

Article number	Content
640F18=900	900 ml

7

8

A large grid of dotted lines for taking notes, spanning most of the page width and height.

1
2
3
4
5
6
7
8

1

2

3

4

5

6

7

8



Body harnesses and accessories

1

2

3

4

5

6

7

8

Triple-Control Above-Elbow Harness

1



21A35=1 Triple-Control Above-Elbow Harness

For fixing the prosthetic socket as well as controlling body-powered prostheses, for right and left side fittings, with perlon cable

Article number	21A35=1
Complete with	Perlondraht

2



21A13=3 Connection Piece for Cable-Activated Hook

Article number	for	Package
21A13=3	10Y32	Perlon Cable

3

4

Spare parts for 21A35=1



21A37=1 Bowden cable

Article number	21A37=1
Spiral length	500 mm

5



10Y31 Ball-Shaft Adapter

Forms a connection element with 10Y32 coupler

Article number	Suitable for
10Y31=1	21A18=2 Perlon cable
10Y31=2	651D4=2 Steel cable
10Y31=7	3/64" Steel cable
10Y31=8	1/16" Steel cable
10Y31=9	3/32" Steel cable

6

7

8

10Y32 Coupler

Forms a connection element with 10Y31=* ball-shaft adapter

Article number	Suitable for
10Y32=1	21A18=2 Perlon cable
10Y32=2	651D4=2 Steel cable



21Y194 Ring

Harness ring with integrated cable guide



21Y195=25 Stainless Steel Buckle

For positioning and fixing the harness strap



21Y197=1 Connecting Bracket

Article number	Colour
21Y197=1	white



21A38 Axilla Pad Set

Article number	Colour	Content
21A38	white	5 pc.



1

2

3

4

5

6

7

8

1



21Y199 Strap Buckle

Article number	for	Consists of:
21Y199	seamless connection of the 623G23 elastic harness strap to the prosthetic socket	Covering Lower part Upper part Spacer sleeve

2



3



21Y203 Spacer Sleeve

Technical Data

Article number	21Y203
for	21Y199 Strap Buckle

4



503F3 Socket Screw with Allen head

Socket Head Screw (2 pieces)

Article number	Thread	Head Ø	Thread length
503F3	M4	8 mm	7 mm

2:1

5



623G23 Elastic Harness Strap

Article number	Colour	Length	Complete with
623G23	white	RM	tunnel-shaped cable guide

6



623H23 Harness Strap

Article number	Colour	Length
623H23	white	RM

7

8

21A18 Nylon wire

Article number	Ø	Length	Order by
21A18=2X1	2 mm	1 m	lfm
21A18=2X5	2 mm	5 m	lfm
21A18=2X10	2 mm	10 m	lfm
21A18=2X25	2 mm	25 m	lfm



29C5 Setting Nut

Article number	Thread	Length	Head Ø	Shoulder Ø
29C5=M4x9	M4	3.6 mm	9 mm	5.5 mm



Technical Data

Article number	29C5=M4x9
Material	Stainless steel

🔍 2:1

1

2

3

4

5

6

7

8

Below-Elbow Harness



21A36=1 Below-Elbow Harness


For fixing the prosthetic socket as well as controlling body-powered prostheses, for right and left side fittings, with perlon cable

Spare parts for 21A36



21A37=1 Bowden cable

Article number	21A37=1
Spiral length	500 mm



10Y31 Ball-Shaft Adapter

Forms a connection element with 10Y32 coupler

Article number	Suitable for
10Y31=1	21A18=2 Perlon cable
10Y31=2	651D4=2 Steel cable
10Y31=7	3/64" Steel cable
10Y31=8	1/16" Steel cable
10Y31=9	3/32" Steel cable



10Y32 Coupler

Forms a connection element with 10Y31=* ball-shaft adapter

Article number	Suitable for
10Y32=1	21A18=2 Perlon cable
10Y32=2	651D4=2 Steel cable



21Y194 Ring

Harness ring with integrated cable guide



1

21Y195=25 Stainless Steel Buckle

For positioning and fixing the harness strap



2

21Y197=2 Connecting Bracket

Article number	Colour
21Y197=2	white



3

21A38 Axilla Pad Set

Article number	Colour	Content
21A38	white	5 pc.



4

5

21Y199 Strap Buckle

Article number	for	Consists of:
21Y199	seamless connection of the 623G23 elastic harness strap to the prosthetic socket	Covering Lower part Upper part Spacer sleeve



6



7

8

1



21Y203 Spacer Sleeve

Technical Data

Article number	21Y203
for	21Y199 Strap Buckle

2



503F3 Socket Screw with Allen head

Socket Head Screw (2 pieces)

- + For universal fastening applications
- + Stainless steel
- + Flat head so that there is little bulging
- + Allen head ensure there are no sharp edges as is the case with slotted screws

2:1

Article number	Thread	Head Ø	Thread length
503F3	M4	8 mm	7 mm

3

4



623G23 Elastic Harness Strap

Article number	Colour	Length	Complete with
623G23	white	RM	tunnel-shaped cable guide

5

6



623H23 Harness Strap

Article number	Colour	Length
623H23	white	RM

7



21A18 Nylon wire

Article number	Ø	Length	Order by
21A18=2X1	2 mm	1 m	lfm
21A18=2X5	2 mm	5 m	lfm
21A18=2X10	2 mm	10 m	lfm
21A18=2X25	2 mm	25 m	lfm

8

29C5=M4x9 Setting Nut

Article number	Thread	Length	Head Ø	Shoulder Ø
29C5=M4x9	M4	3.6 mm	9 mm	5.5 mm

1

2

3

4

5

6

7

8

Body Harnesses 21A20 / 21A19

1



21A20 Below-Elbow Harness

Complete

Article number	Harness design
21A20=1	Nylon
21A20=2	Steel rope

2

3



21A19 Above-Elbow Harness

Triple-Control Harness, complete

Article number	Harness design
21A19=1	Nylon
21A19=2	Steel rope

4



5

Spare parts for body harnesses

6



10Y3 Cable clamp

7



10Y25 Clamp Sleeve, short

To be clamped to the 651D4=2 steel cable

8



10Y26 Clamp Sleeve, long


To be clamped to the 651D4=2 steel cable

21A3 Above-Elbow Cable Control Unit, complete

Article number	Harness design
21A3=1	Nylon
21A3=2	Steel rope

Technical Data

Reference Number	21A3
Spiral length	500 mm
Complete with	21A7 coupling piece and 21A13 connection piece



1

2

21A4 Below-Elbow Cable Control Unit, complete

Article number	Harness design
21A4=1	Nylon
21A4=2	Steel rope

Technical Data

Reference Number	21A4
Spiral length	300 mm
Complete with	21A7 coupling piece and 21A13 connection piece

3

4

21A5 Eyelet Cable Anchor, small



5

21A6 Spiral Nut



6

21A7 Coupling Piece

Technical Data

Article number	21A7
for	Nylon rope



7

8

1



21A11 Threaded Fitting, long

For screwing on the 21A18=2 perlon cable

2



21A12 Threaded Fitting, short

For screwing on the 21A18=2 perlon cable

Article number	Package content
21A12	2 pcs.

3



21A17=3.8x1000 Cable-Control Unit Spiral

Article number	Ø
21A17=3.8x1000	3.8 mm

4



21A18 Nylon wire

Article number	Ø	Length	Order by
21A18=2X1	2 mm	1 m	lfm
21A18=2X5	2 mm	5 m	lfm
21A18=2X10	2 mm	10 m	lfm
21A18=2X25	2 mm	25 m	lfm

5

6




21A23 Wire Coupling

7



21A24 Eyelet Cable Anchor, large

8



21A25 Spring

Technical Data	
Article number	21A25
for	21A7 coupling piece

21Y37 Cable Strap Connector

Inside width 25 mm

Article number	Inside width
21Y37	25 mm



1

21Y120=47 Leather Lamella

Technical Data

Article number	21Y120=47
for	straps and harnesses



2

623G4 Double Rubber Belt

For 29R127 Clip Fastener

Article number	Colour	Length	Width
623G4=1	gray with white stripes	1 m	30 mm
623G4=5	gray with white stripes	5 m	30 mm
623G4=10	gray with white stripes	10 m	30 mm



3

4



5



623H2=25x5 Carrying Strap

Article number	Colour	Length	Width
623H2=25x5	white	5 m	25 mm

6

651D4=2 Steel Cable

Article number	651D4=2
Complete with	plastic cover



7

8

Accessories for body harnesses

1



21A21=* Above-Elbow Cable Control Unit,

Article number	Harness design
21A21=1	Nylon
21A21=2	Steel rope

Technical Data

Reference Number	21A21=*
Spiral length	500 mm
Complete with	21A23 wire coupling

2

3



21A22=* Below-Elbow Cable Control Unit, complete

Article number	Harness design
21A22=1	Nylon
21A22=2	Steel rope

Technical Data

Reference Number	21A22=*
Spiral length	300 mm
Complete with	21A23 wire coupling

4

5



16H1 Flexible Elbow Joint Connection

"V"-shape

Article number	Complete with
16H1	three-point suspension

6



16H2 Flexible Joint Splints

pair

Technical Data

Article number	16H2
for	below-elbow prostheses

7



501S28 Flat head screw

nickel-plated

Technical Data

Article number	501S28=M3.5x5
for	16H1 and 16H2

8

21A5 Eyelet Cable Anchor, small



1

21A16 “D” Ring

Article number	Inside width
21A16	13 mm



2

Technical Data

Article number	21A16
for	lastic strap

3

503F3 Socket Screw with Allen head

Socket Head Screw (2 pieces)

Article number	Thread	Head Ø	Thread length
503F3	M4	8 mm	7 mm



🔍 2:1

4

507S15 Lamination Disk, serrated

Lamination Plate with hole (2 pieces)

Article number	Outside Ø	Hole Ø
507S15	13.8 mm	3 mm



🔍 2:1

5

516S3 Pad Screw Head

Brass, nickel plated,

Article number	Thread	Thread length
516S3	M4	5.5 mm



🔍 2:1

6

514Z3=25 Loop

Article number	Inside width
514Z3=25	25 mm



8

1



514K3=27 Clamp Buckle

2

504H3=11 Hollow Rivet

nickel-plated, open end

Article number	Head Ø	Material
504H3=11	11 mm	Brass

3



10Y19=1 Screw Coupling

Connection piece between perlon cable and steel cable as well as steel cable and steel cable.

Article number	Consists of:
10Y19=1	coupling sleeve fixed coupling screw, to be attached to the 651D4=2 steel cable.

4



10Y19=2 Screw Coupling

Connection piece between steel cable and perlon cable or spectra cable and between perlon cable and perlon cable or spectra cable.

Article number	Consists of:
10Y19=2	coupling sleeve coupling screw, to be screwed to the 21A18=* perlon cable or through which a spectra cable is to be threaded and then knotted

5

6



736Y6 Clamping Tool

For attaching the coupling screw and clamp sleeve as well as 10Y31=2 Ball-shaft Adapter and 10Y32=2 Coupler to the 651D4=2 Steel Cable

7

8

21A29 Axillary Protector

Article number	Width
21A29=18	18 mm
21A29=25	25 mm

Technical Data

Article number	21A29=18	21A29=25
for	body harnesses	body harnesses



15Y1 Sleeve Protection Pad

Rubber with leather cover

Article number	Length	Width
15Y1=16	16 cm	85 mm
15Y1=18	18 cm	85 mm



1

2

3

4

5

6

7

8



1

2

3

4

5

6

7

8

Socket Comfort, Materials & Workplace Arm Prosthetics

1

2

3

4

5

6

7

8

1



704B5 Mounting Vice

Tool to assist in the disassembly and reassembly

Technical Data

Article number	704B5
for	finger thumb group

2



709S10=2 Allen Wrench

Technical Data

Article number	709S10=2
for	503F3 Socket Screw

3

4



711M17 Pulling Unit

For easy disassembly of the 11S61 Gear Unit

5

6




711M18 Mounting Wrench

For accurately tightening the 11S2 Locking Unit and for releasing the disengaged ball bearings

Article number	Complete with
711M18	503F3 preset torque

7



711M20 Alignment Rod

In combination with the 711M50 Mounting Adapter for installing and removing the 9E369 and 9E370 4in1 Controller LS or 9E420 7in1 Controller or for setting and releasing the 9E349 and 9E350 Electronic Control.

Article number	Complete with
711M20	Allen head

8

711M53 Vacuum Forming Tool

For arm prostheses

Article number	711M53
Outer Ø	26 cm



1

711M56 Vice Clamping Jaws

Two clamping jaws are required to equip the 704B5 Mounting Vice.



2

726W7=28.5 Crown Cutter

To cut electrode mounting holes

Technical Data

Article number	726W7=28.5
for	flexible mounting of electrode, in thin-walled plastic sockets



4

726W9 Conical Plastic Drill

made of high speed steel

Article number	Ø
726W9=14	14 mm
726W9=20	20 mm
726W9=30	30 mm



5

743A19 Alignment Tool

Article number	Thread	Ø	Length
743A19	M8	12 mm	169 mm

Technical Data

Article number	743A19
for	10S16 Lamination Rings
Consists of:	743Y167 Alignment Rod 743Y42=34 Foam Insert for sizes 5 and 5 1/2 743Y42=38 Foam Insert for sizes 6 and 6 1/2



7

8

1



743A18 Alignment Tool

Article number	Thread	Ø	Length
743A18	118	12 mm	280 mm

Technical Data

Article number	743A18
Consists of:	743Y41 Schäumeinsatz, Ø 40 mm 743Y42=45 743Y42=50 743Y42=54

- When ordering replacement inserts please state the diameter (corresponds to the outer diameter of the 10S1 Lamination Ring)

3



11S12 Socket Attachment Piece

For fastening the 743Y41 Alignment Rod,
 For fabricating an arm prosthesis

4

5



12V10 Tube Valve for Suction Socket

In combination with the 13E202 MYOBOCK Electrode, the 12V10 Tube Valve for Suction Socket creates an airtight seal of the socket.

6



99B13 PVC Connection Tube

Formed to 90° angle, acts as a connection channel between the inner and outer socket

Article number	Colour	Ø
99B13=16	skin-coloured	16 mm
99B13=16-7	black	16 mm
99B13=21	skin-coloured	21 mm
99B13=21-7	black	21 mm

7

8

99B83 Tube Dummies

For fabricating vacuum formed inner sockets

Article number	Ø
99B83=16	16 mm
99B83=21	21 mm



1

743Y42=24 Foam Insert

can be used with alignment rod of the
 743A18 Alignment Tool

Technical Data

Article number	743Y42=24
for	Physolino Babyhand



2

3

641T10=M Plaster Shirt

The Plaster Shirt can be used for taking a plaster cast as part of the prosthetic arm fitting process. The plaster shirt is suitable for upper arm residual limbs and shoulder disarticulation.



4

641T9=M Plaster Sock

The Plaster Sock can be used for taking a plaster cast as part of the prosthetic arm fitting process. The plaster sock is used for forearm residual limbs.



5

6

7

8



616T69 ThermoLyn soft (EVA), skin-coloured


Now also available with antibacterial features (616T269)!

For fabricating flexible arm prosthesis sockets, high shrinkage if cooled down too quickly, translucent (colour depends on skin colouring), processing temperature 160 °C in convection and infrared ovens

Article number	Colour	Length	Width
616T69=6	skin-coloured	400 mm	400 mm
616T69=8	skin-coloured	400 mm	400 mm
616T69=10	skin-coloured	400 mm	400 mm
616T69=12	skin-coloured	400 mm	400 mm



 646F265=D

 646A230=D
646D300=D
646D119=D



Practical recommendation:

- We recommend using the 503F3 socket screw with Allen head in combination with the 29C3 or 29C5 Setting Nut
- The 633D5 Double Sided Adhesive Tape can be used to fasten the flexible definitive inner socket within the container socket



633F11 Silicone Grease 400, medium

For the planetary gear set of the drive unit as well as for all cable plug connections and other places that need to be protected to prevent penetration of sweat. Can be applied for cast isolation.

Article number	Net contents
633F11	0.1 kg



633F14=* Special Lubricant

(Molycote Paste DX, White)

Article number	Colour	Delivery	Net contents
633F14=0.050	white	Tube	6 g
633F14=1	white	Dose	9,5 g

Technical Data

Reference Number	633F14=*
Colour	white
for	all accessible gears and axes in the System Electric Hands and Greifers



633F37=0.02 Variotrac 68

For lubricating the 9E79 Transmission,

Article number	Delivery	Net contents
633F37=0.02	Flasche	20 g



1

633S2 Procomfort Gel

Acts as a lubricant to aid in donning the prosthetic glove over the inner hand.

Article number	Net contents
633S2	250 ml



2

634A58 Isopropyl Alcohol

For cleaning sensitive plastics such as PVC, PS, ABS, acrylic, PC

Article number	Colour	Net contents
634A58	transparent	1 l



i 646F297=D

4

636W23 UHU-Plus Endfest 300 Adhesive

Two-component adhesive for strong adhesions

Article number	Colour
636W23	honey



i 646F297=D

5

6



Practical recommendation:

- The higher the curing temperature (up to approximately +180°C), the higher the strength of the adhesion
- Using the 642B2 Measuring Cup has proven useful for mixing
- Can be applied with the 699Y3 Wooden Spatula

7

8

1



636W60 Loctite® 243

For thread lock of the electric hand, also hardens on brass, up to thread M36

Article number	Net contents
636W60	50 ml

2



Practical recommendation:

- Can be activated and loosened by applying heat

3



636K11 Cyamet Rapid Adhesive (Superglue)

For adhering silicone rubber to acrylic resin laminates

Article number	Colour	Net contents
636K11	transparent	20 g

4



 646F297=D

5



636K13 Loctite® 241

Article number	Colour	Net contents
636K13	blue	50 ml

6



Practical recommendation:

- Can be activated and thus removed by applying heat

8

Loctite® is a registered trademark of Loctite.

636K14 Loctite® 601

For locking the threaded stud,
Bottle with 50 ml net contents

Article number	Colour	Delivery	Net contents
636K14	green	Flasche	50 ml



Practical recommendation:

- Can be activated and loosened by applying heat

640F12 Special Cleaner

Technical Data

Article number	640F12
for	prosthetic gloves



640F13 Spray Pump Bottle

Empty container holds 90 ml

Technical Data

Article number	640F13
for	liquids such as 640F12 Special Cleaner



9E161 Screw Cap

For retaining the transmission in the gear housing



Loctite® is a registered trademark of Loctite.

1



636W22 UHU, hard

for locking the motor,
 Tube with 35 ccm net contents

Article number	Colour	Delivery	Net contents
636W22	colourless	Tube	35 ml



2

i 646F297=D

3




636W34 Rubber Adhesive

Article number	Colour	Net contents
636W34	beige	60 g



i 646F297=D

4




637L7 Resin Rapid Plumb

60% SN, Ø 1 mm,
 Coil with 1/2 kg net weight

Article number	Ø	Weight
637L7	1 mm	0.5 kg

5

6



640Z6 Work Mat

7



642G1 Glass Bottle

Article number	for	Scope of Delivery
642G1	lubrication materials, solvents, etc	2 pcs.

8

706F1 Flat Nose Pliers

Article number	Length
706F1	125 mm



1

706R3 Universal Pliers, large

Technical Data

Article number	706R3
for	drive unit



2

3

706R4 Universal Pliers, small

Article number	Length
706R4	115 mm



4

5

706Z2 Side Cutting Pliers

Article number	Length
706Z2	120 mm



6

799P1 Tweezers

Article number	Length
799P1	120 mm



7

8

1



799P2 Tweezers

Article number	Length
799P2	155 mm

2



799Y1 Tool Board

3

4



709G1=7x7 Open-end Wrench

For adjustment of the slip clutch

5

6



709H1 Hook Wrenches

To assemble the 8E39 System Electric Hands with the corresponding lamination ring.
(May also be used to reset the ball bearings in the quick disconnect locking unit.)

7



709Z6=1.8 Special Wrench

For adjusting the drive unit

8

709Z7 Special Wrench

For replacing the hand cable



1

710H3 Precision Screwdriver Set

Blade widths: 1/1.5/1.8/2.3/2.9/3.6 mm



2

3

711M64 Donning Tool for Michelangelo

The Donning Tool is essential for the correct application of the AxonSkin glove in order to protect the mechanism of the Michelangelo Hand and for easier handling by the technician.



4

5

710H4 Screwdriver

Article number	Blade width
710H4=3	3.0 mm
710H4=4.5	4.5
710H4=6	6.0



6

711M1 Assembly Tool

One side with M12x1.5 outer thread and other side with M12x1.5 outer thread for Otto Bock System Electric Hands and Michelangelo Hand



7

8

1



711M2 Mounting Plate

for attaching to a workbench

Technical Data

Article number	711M2
for	711M1 Assembly Tool

2



711M3 Assembly Tool

For accommodating System Electric Hands

Article number	Complete with	Consists of:
711M3	quick disconnect wrist	10S4 Coupling Piece 11S4 Lock Ring 501Z2=M6x25 Allen screw

3

4



711M5 Mounting Aid

For exchanging the DC motor and adjusting the drive unit

5

6



711M7 Mounting Aid

For attaching the 8E39 System Electric Hands to 711M1 Assembly Tool.
 Makes putting on and removing the inner hand and cosmetic glove easier

7



711M12 Soldering Jig

For soldering 9E167 Plug to the coaxial bushing or the hand cable as well as for soldering the plug bushing to the motor cable

8

711M16 Assembly Tool

For assembling the 10S17 Electric Wrist Rotator or 10S4 Coupling Piece in the 10S1 Lamination Ring



1



2

3

Article number	Length
718H5	160 mm

Technical Data

Article number	718H5
Weight	0.05 kg
for	synthetic materials
Complete with	protective cap und bar

4

718Y1 Spare Blade for Deburring Knife

Article number	718Y1
Weight	0.001 kg
for	718H5 Deburring Knife



5

743F1 Pinch Gauge

For testing the grip force of System Electric Hands



6

7

799L1=* Soldering Iron

Complete with 799L3 Soldering Iron Tip

Article number	Volt	Watt
799L1=220	220	16



8

1



799L2 Stand for Soldering Iron

2

3

4

5

6

7

8

1
2
3
4
5
6
7
8



Test, Inspection and Adjustment Equipment

1

2

3

4

5

6

7

8

1



757P2=4 Test Cable

Technical Data

Article number	757P2=4
for	drive units and DC motors

2



3

4



757P2=6 Test Cable

Technical Data

Article number	757P2=6
for	757B15 X-ChangePack

5



6

7

8

757P2=8 Test Cable

Technical Data

Article number	757P2=8
for	13E51=2 and 13E132 Battery Connection Cables



1

2

3

757P2=9 Test Cable

Technical Data

Article number	757P2=9
for	10S17 Electric Wrist Rotator and 8X15 or 8X19 DC Motors



4

5

6

7

8

1



757P2=10 Test Cable

Technical Data

Article number	757P2=10
for	757B13 Interchangeable Battery

2



3

4



757P23 Test Adapter

Technical Data

Article number	757P23
for	8E38, 8E39, 8E41 and 8E44 System Electric Hands or 8E33 or 8E34 System Electric Greifers

5



6

7

8

757P28 Test Plug

When used with the 757T16 System Tester, this test plug may be used to test the function of DMC, DMC plus and DMC VariPlus control.

Technical Data

Article number	757P28
for	use with the DMC system



647H149

757P35 Test Plug

Technical Data

Article number	757P35=1	757P35=2
for	EnergyPack 757B20	EnergyPack 757B21



1

2

3

4

5

6

7

8

1



757P39 MyoBoy Communication Cable

Connects the MyoBoy and 757Z185=2 or 757Z191=2 Battery Receptacle. Permits signal measurement while the prosthesis is in use.

2



3

4



757T16 System Tester

For MyoBock-System,
without test cable, with 757L16-1 AC adapter with EU and US plug.
The system tester is used to test the function of individual components and assemblies in the MyoBock System.
System Electric Hands, System Electric Greifers, connection cables, battery connection cables, batteries and assemblies may be tested either individually or in combinations.

⦿ Please order a test cable separately according to requirements

5



757T13 MyoSelect

The 757T13 MyoSelect is used to identify and adjust Myobock® components, such as System Electric Hands, System Electric Greifers, MyoRotronic and ErgoArm Electronic plus.

The MyoSelect is connected to the MyoBock component and then displays information on the type of component and on the currently selected control mode. The multi-function button can be used to select and adjust alternative control modes.

With the 757T13 MyoSelect, it is even possible to adjust the speed of the MyoHand VariPlus Speed, SensorHand Speed and the System Electric Greifer.

Please note that in order to make adjustments with the MyoSelect, the components must first be equipped with a black coding plug !

New components, such as the MyoHand VariPlus Speed and System Electric Greifer DMC VariPlus are only designed for use with the 757T13 MyoSelect!

6

7

8

1
2
3
4
5
6
7
8



Myo Service Parts

1

2

3

4

5

6

7

8

Service Parts for System Electric Hands

1



9E53 Hand Cable, flat

Article number	for	Complete with
9E53	System Electric Hands	9E70 Cable Bushing

2

3



9E68 Tube Chassis

Article number	For hand size
9E68=7 1/4	7 1/4
9E68=7 3/4	7 3/4
9E68=8 1/4	8 1/4

Reference Number	9E68
for	System Electric Hands with threaded stud M12X1.5, in the sizes 7 1/4, 7 3/4 and 8 1/4
Complete with	Threaded Stud M12X1,5

4

5




9E83 Lock Nut

Technical Data

Article number	9E83
for	9E53 and 9E125 Hand Cable

6



9E125=200 Hand Cable, twin-wire

Article number	for	Complete with
9E125=200	8E12 Electric Hand	9E48 Cable Bushing

7



9E167 Cable Connector

Technical Data

Article number	9E167
for	9E53 and 9E168 Coaxial Bushing

8

9E168 Coaxial Bushing

Article number	9E168=8	9E168=7
for	MyoBock Digital Twin System	MyoBock DMC plus-System
Complete with	9E167 Cable Connector and 501S17=M3×5 Countersunk Head Screw	



1

9E188 Distributor

Article number	9E188
for	System Electric Hands without electronic control, for example, 8E12 or 8E37
Complete with	501T8=M2x4 Cap Screw (2 pieces) 9E189=1 Plug Housings (2 pieces) 519S5=0.8 Plug Bushings (4 pieces)



2

3

9E189=1 Plug Housing, dark grey


Article number	9E189=1
for	DC motors



4

9E194 Drive Unit

Article number	9E194
Weight	80 g
Approx. max current	2.6 A
Approx. no-load current	150 mA
Supply voltage	6/7.2 V
for	Electric Hands in sizes 7 1/4, 7 3/4 and 8 1/4
Complete with	ball bearing and lock element
Consists of:	9E202 Stainless Steel Spacer 9E250=1 Housing with Drive Element 9E205 Lock Element 9E206=2 Planet Set 9E207=1 Drive Ring 9E251 Spacer 9E209=4 Wrapped Spring 616T95=3 ThermoLyn PE 200 8X16 DC Motor 9E56 Adjustment Ring 9E208 Plastic Disk 9E219=3 Spring Support Ring



5

6

7

- ⦿ Not suitable for DMC and DMC plus systems.
- ⦿ 9E205 Lock Element: lubricate the front surface (friction surface for the slip ring) with 633F14 Molykote.
- ⦿ 9E206=2 Planet Set: lubricate with 633F11 Silicone Grease.
- ⦿ 9E56 Motor: secure with 636W22 Adhesive.
- ⦿ For repairs, please state hand type (e.g. 8E38=1) or whether or not an electronic control has been built into the System Electric Hand.

8

1



9E255=1 Slip Clutch

Article number	for	with
9E255=1	Electric Hands in sizes 7 1/4, 7 3/4 and 8 1/4	9E162 Axis with Lock 501S24=M3x4 Flat Head Screw

2



9E274 Transmission

Article number	For hand size	Side
9E274=L	7 1/4, 7 3/4 und 8 1/4	Left (L)
9E274=R	7 1/4, 7 3/4 und 8 1/4	Right (R)

Technical Data

Reference Number	9E274
for	8E38, 8E39 and 8E41 Electric Hands in sizes 7 1/4, 7 3/4 and 8 1/4,
Consists of:	9E79=L/R Transmission 9E46 Wire Ring 513S1 Coil Spring 9E43 Double Spring (2 pieces) 9E33 Planetary Gear Wheel (3 pieces) 9E31 Inner Gear Ring 519S47 Spacer Washer (3 pieces) 9S197 Housing with Roller Bearing 509Y1=2.0 Ball (10 pieces) 9E278 Bevel Gear, complete 9E277 Bevel Gear, short 507U3=4.5x1.8x0.3 Washer 501T8=M1.7x3 Cap screw 9S188 Support Bearing 501T8=M2x10 Cap screw 9E161 Screw Cap

- Lubricate the transmission with 10 drops of 633F37=0.02 Silicone Oil.
- For repairs, please state hand type (e.g. 8E38=1-L8 1/4)

5

6



9E380 Electronic Control

For System Electric Hands DMC plus

The DMC plus control features a DMC and DMC plus control mode.

Article number	Side
9E380=L	Left (L)
9E380=R	Right (R)

647H336

Technical Data

Reference Number	9E380
for	System Electric Hands in sizes 7 1/4, 7 3/4 and 8 1/4 with 501T8=M2.6x6 Screws (2 pieces)
Complete with	electronic switch-off 13E185 Function Plug integrated hand switch

7

8

9E381 Electronic Control

For System Electric Hands Digital Twin


The Digital Twin control features a digital and dual channel control mode.

Article number	Side
9E381=L	Left (L)
9E381=R	Right (R)

Technical Data

Reference Number	9E381
for	System Electric Hands in sizes 7 1/4, 7 3/4 and 8 1/4 with 501T8=M2.6x6 Screws (2 pieces)
Complete with	electronic switch-off 13E185 Function Plug ntegrated hand switch



 647H337

9E385 Electronic Control


For Transcarpal Hand DMC plus

Article number	Side
9E385=L	Left (L)
9E385=R	Right (R)

Technical Data

Reference Number	9E385
Complete with	9E387 Clamp Piece (for electronic control) 13E185 Function Plug 501S122=M2.5x22 Flat Head Screw



 647H403

9E386 Electronic Control

For Transcarpal Hand Digital Twin

Article number	Side
9E386=L	Left (L)
9E386=R	Right (R)

Technical Data

Reference Number	9E386
Complete with	9E387 Clamp Piece (for electronic control) 13E185 Function Plug 501S122=M2.5x22 Flat Head Screw



 647H404

9E388 Cable Seal

Technical Data

Article number	9E388
for	9E53 Hand Cable (for the Transcarpal Hand)



1



9S72 Hand Chassis

Article number	For hand size
9S72=7 1/4	7 1/4
9S72=7 3/4	7 3/4
9S72=8 1/4	8 1/4

Technical Data

Reference Number	9S72
for	System Electric Hands with quick disconnect wrist, in the sizes 7 1/4, 7 3/4 and 8 1/4
Complete with	quick disconnect wrist
with	11S8 Disconnect Piece 11S9 Ball Bearing Cage 11S10 Disconnect Piece (10 pieces) 509Y1 Ball
Consists of:	9S74 Hand Chassis with quick disconnect wrist, 11S6 Ratchet Ring with Inner Gear Teeth 11S27 Spacing Washer 11S7 Pressure Ring 11S2 Locking Unit Quick disconnect piece with ball bearings for locking and coupling to lamination ring

2

3

4



9S74 Hand Chassis with quick disconnect wrist,

Article number	For hand size
9S74=7 1/4	7 1/4
9S74=7 3/4	7 3/4
9S74=8 1/4	8 1/4

Technical Data

Reference Number	9S74
for	System Electric Hands with quick disconnect wrist, in the sizes 7 1/4, 7 3/4 and 8 1/4
Complete with	Tube Chassis with threaded adapter and 11S5 Ratchet Ring

5

6



9S83 Thumb

Article number	For hand size	Side
9S83=L	7 1/4, 7 3/4	Left (L)
9S83=R	7 1/4, 7 3/4	Right (R)

Technical Data

Reference Number	9S83
with	9S8 Finger Thumb Tip 507U42=4.0 Lock Washer

7

8

9S85 Finger Group

Article number	For hand size	Side
9S85=L	7 1/4, 7 3/4	Left (L)
9S85=R	7 1/4, 7 3/4	Right (R)

Technical Data

Reference Number	9S85
Consists of:	9S84 Finger Tip (2 pieces) 502S16=BM 3.5 Hexagon nut 9E74 Steel axis 9E65=7 3/4 Finger Connection Bar




1

2

9S103=7 1/4 Tube chassis with lamination ring

Article number	9S103=7 1/4
for	System Electric Hands for wrist disarticulation In size 7 1/4
Consists of:	9S74=7 1/4 Tube Chassis with threaded adapter and ratchet ring 9S105 Cable Guide 9S110=50 Lamination Ring 11S30 Friction Ring with Inner Gear Teeth 11S27 Spacing Washer 11S7 Pressure Ring



3

4

9S103=7 3/4 Tube chassis with lamination ring

Article number	9S103=7 3/4
for	System Electric Hands for wrist disarticulation In size 7 3/4
Consists of:	9S199=7 3/4 Tube Chassis with threaded adapter and 11S5 Ratchet Ring 9S105 Cable Guide 9S110=50 Lamination Ring 11S30 Friction Ring with Inner Gear Teeth 11S27 Spacing Washer 11S7 Pressure Ring



5

9S103=8 1/4 Tube chassis with lamination ring

Article number	9S103=8 1/4
for	System Electric Hands for wrist disarticulation In size 8 1/4
Consists of:	9S199=8 1/4 9S105 Cable Guide 9S110=54 Lamination Ring 11S30 Friction Ring with Inner Gear Teeth 11S27 Spacing Washer 11S7 Pressure Ring



7

8

1



9S203 Finger-Thumb Group

Article number	For hand size	Side
9S203=L	8 1/4	Left (L)
9S203=R	8 1/4	Right (R)

Technical Data

Reference Number	9S203
for	System Electric Hands in size 8 1/4
Consists of:	9S204 Finger Group 9S211 Thumb 9S92 Finger Axis 501S70=M3.5x6 Flat Head Screw (4 pieces) 9S91 Thumb Axis 9E78 Bow Spring 9S194 Protective Cap

2

3



9S204 Finger Group

Article number	For hand size	Side
9S204=L	8 1/4	Left (L)
9S204=R	8 1/4	Right (R)

Technical Data

Reference Number	9S204
Consists of:	9S8 Finger Thumb Tip 502S16=BM 3.5 Hexagon nut 9E74 Steel axis 9E65=7 3/4 Finger Connection Bar

4

5



9S211 Thumb

Article number	For hand size	Side
9S211=L	8 1/4	Left (L)
9S211=R	8 1/4	Right (R)

Technical Data

Reference Number	9S211
with	9S236 Thumb Tip 507U42=4.0 Lock Washer

6

7


8

9S255 Finger-Thumb Group

Article number	For hand size	Side
9S255=R	7 1/4, 7 3/4	Right (R)
9S255=L	7 1/4, 7 3/4	Left (L)

Technical Data

Reference Number	9S255
for	Transcarpal Hands in sizes 7 1/4, 7 3/4
Consists of:	9S85 Finger Group 9S83 Thumb 9S257 Finger Chassis 9S92 Finger Axis 501S70=M3.5x6 Flat Head Screw (4 pieces) 9S91 Thumb Axis 9S194 Protective Cap



1


2

9S256 Finger-Thumb Group

Article number	For hand size	Side
9S256=L	8 1/4	Left (L)
9S256=R	8 1/4	Right (R)

Technical Data

Reference Number	9S256
for	Transcarpal Hands in size 8 1/4
Consists of:	9S204 Finger Group 9S211 Thumb 9S257 Finger Chassis 9S92 Finger Axis 501S70=M3.5x6 Flat Head Screw (4 pieces) 9S91 Thumb Axis 9S194 Protective Cap



3

4

501S7=M3x5 Countersunk Head Screw

For connecting the 9S197 Housing with roller bearing to the finger chassis



5

6

501S17=M3x5 Countersunk Head Screw

Technical Data

Article number	501S17=M3x5
for	9E168 Coaxial Bushing



7

501T8=M2x4 Cap screw

For securing the distributor to the hand switch



8

1



502S97=M5 Hexagon nut

Article number	Complete with
502S97=M5	flange

2



519S5=0.8 Plug Bushing

Technical Data

Article number	519S5=0.8
for	plug housing

3



520E30 Solder Lug

Technical Data

Article number	520E30
for	System Electric Hands in sizes 7 1/4, 7 3/4 and 8 1/4, for connecting the system ground connection

4

5

6


7

8

Service Parts for System Electric Greifers

9E221 Coaxial Bushing

Article number	for	for	with
9E221=7	MyoBock DMC plus System	9E264=6	501S17=M3x5 Countersunk Head Screw
9E221=8	MyoBock Digital Twin System	8E33=7	501S17=M3x5 Countersunk Head Screw
9E221=9	MyoBock DMC VariPlus System	8E33=9	501S17=M3x5 Countersunk Head Screw



9E235 Connecting Cable

Article number	for	Complete with
9E235	System Electric Greifers for 8E34=9 and 8E34=7 wrist disarticulation	501S17=M3x5 Countersunk Head Screw



9E236=2 Distributor

Technical Data

Article number	9E236=2
for	8E32=6 System Electric Greifers



1

2

3

4

5

6

7

8

1



9E264=7 Electronic Control

Technical Data

Article number	9E264=7
for	Digital Twin system 8E33=7 and 8E34=7 System Electric Greifers

2

 647H74

3



9E264=9 Electronic Control

Technical Data

Article number	9E264=9
for	DMC VariPlus System 8E33=9 and 8E34=9 System Electric Greifers

4



9E346 Drive Unit for System Electric Greifer

Article number	Complete with
9E346	electronic battery-saving circuit

5

6



9S117 Reduction Gear

7



9S127=1 Finger Cover

Article number	with
9S127=1	501S54=M3x10 Flat Head Screw (2 pieces) 9S147=2PAA Rubber pad (pair)

8

9S136=1 Finger Housing Set

Article number	with
9S136=1	501S27=M2x6 Oval Head Countersunk Screw (4 pieces) 9S147=1PAA Rubber pad (pair)



1

11S2 Locking unit

Quick disconnect piece with ball bearings for locking and coupling to lamination ring

Article number	with
11S2	11S8 Disconnect Piece 11S9 Ball Bearing Cage 11S10 Disconnect Piece (10 pieces) 509Y1=3 Ball (20 pieces)



2

3

11S6 Ratchet Ring with Inner Gear Teeth



4

11S7 Pressure Ring



5

11S27 Spacing Washer



6

11S30 Friction Ring with Inner Gear Teeth



7

13E29 Collector Brush

Technical Data

Article number	13E29
for	Greifer electronic control and distributor



8

Service Parts for Electric Wrist Rotator

1



11S61 Wrist Drive with Coaxial Plug

For 10S17 Electric Wrist Rotator, with 9E365 Protection Plug

2

3

4

5

6

7

8

A			
Above-Elbow Cable Control Unit	325, 328	Cable-activated Hook for Children	158
Above-Elbow Harness	324	Cable-Activated Hook for Youths	158
Adapter	104, 178, 196, 220, 232, 234 f, 252	Cable-Control Unit Spiral	154, 326
Adapter for Australia	115	Cable clamp	324
Adapter for the UK	115	Cable Connector	360
Adapter with Flexion	220	Cable Pull Switch	124
Alignment Aid for ErgoArm	104, 196	Cable Seal	363
Alignment Rod	334	Cable Strap Connector	327
Alignment Tool	335 f	Cap screw	367
All-purpose Clamp	168	Cap Screw	235
All-purpose Hook	163, 168	Car Charging Cable 12 V	115
All-purpose Ring	168	Carrying Strap	327
Allen Wrench	86, 334	Chassis	143 f, 153 f, 210
Arm Bar	244	Chassis with Quick-Disconnect Wrist	86
Arm Bars	244	Clamp Buckle	330
Arm Bar with Joint	245	Clamping Tool	330
Arm Frame	245	Clamp Plate	192
Assembly Tool	30, 345 ff	Clamp Sleeve	324
Attachment Plate with short stud	156, 167, 183 f	Clamp Stopple Set	104, 196
Attachment Plate with stud	155, 166, 183 f	Coaxial Bushing	86, 361, 369
Automatic Forearm Balance	198	Coaxial Plug	91
Axilla Pad Set	317, 321	Coding Plug Set	131
Axillary Protector	331	Collector Brush	371
AxonCharge Integral	28	Colouring Pen	81
AxonEnergy Integral	27	Conical Plastic Drill	335
AxonMaster	28	Connecting Bracket	317, 321
AxonSkin Natural	25	Connecting Cable	126 f, 369
AxonSkin Visual / AxonSkin Black	26	Connection Cable	127
AxonSoft	29	Connection Cable with On-Off Switch	126
AxonWrist	27	Connection Piece	165
		Connection Piece for Cable-Activated Hook	165, 316
		Connection Piece for Hook	166
		Connection Ring with groove	225, 228 f, 231
		Connector Screw	166
		Control Element	123
		Controller 7in1	130
		Controller LS 4in1	130
		Cosmetic Glove for Children	261 ff
		Cosmetic Glove for Men	284 ff
		Cosmetic Glove for Women	269 ff, 282 f
		Countersunk Head Screw	194 f, 367
		Coupler	317, 320
		Coupling Piece	91, 325
		Crown Cutter	335
		Cyamet Rapid Adhesive	340
		Cylindrical Lamination Ring	181
B		D	
Ball-Shaft Adapter	316, 320	Derma Clean	311
Ball Cap	193		
Battery Connecting Cable	116		
Battery Receptacle	110 f		
Battery Receptacle Set	110 ff		
Bearing Nut	246		
Below-Elbow Cable Control Unit	325, 328		
Below-Elbow Harness	320, 324		
Bowden cable	316, 320		
Bow with wedge lock	242		
Bushing	246		
C			
Cable-activated Hook for Adults	162		

1	Derma Prevent	311	Finger Group	365 f		
	Derma Protection ArmComfort	65	Finger Housing Set	371		
	Derma Repair	311	Finger Tip	72 f		
	Distributor	117, 361, 369	Finger Tip Blank Set	72		
	Donning Tool for Michelangelo	30, 345	Finger Tip Set	72		
	Double-Cable System Hand	149 ff	Flange Button Head Socket Screw	143, 155, 210		
	Double Joint Arm Bars	243	Flat head screw	328		
	Double Rubber Belt	327	Flat Nose Pliers	343		
	Double Spring	164	Flexible Elbow Joint Connection	328		
	Drive Unit	361	Flexible Joint Splints	328		
2	Drive Unit for System Electric Greifer	370	Flexion Cable Retainer	225, 228 f, 236		
	Dummy Set	312	Flexion Device	182		
	DynamicArm Facelift	97	Foam Insert	337		
	3	E		Forearm	192	
		Eccentric	194	Friction Ring	232, 252	
		Elastic Harness Strap	318, 322	Friction Ring with Inner Gear Teeth	371	
		Elbow Component	197 ff	4	G	
		Elbow Joint	192		Glass Bottle	342
		Elbow Joint Bars	240		5	H
		Elbow Joint Bar with Cable Lock	240			Hand Brush
Elbow Joint Bar without Lock		241	Hand Cable			360
Elbow Set-Up for Children		188	Hand Chassis			364
ElbowSoft		135	Hand Chassis with quick disconnect wrist			364
Electric Greifers for other manufacturer's systems	71	Harness Pull Switch	124			
Electric Hand 2000	35	Harness Strap	318, 322			
Electric Wrist Rotator	91	Hexagon nut	368			
Electrode	121	Hole Covering	193			
Electrode Accessories	122	Hollow Rivet	330			
Electrode Accessory Set	122	Hook Wrenches	344			
Electrode Cable	127	6	I			
Electrode Mounting Bracket Set	122		Interchangeable Battery	108		
Electronic Control	362 f, 370		Isopropyl Alcohol	339		
Endoskeletal Adapter	93		7	J		
EnergyPack	109			Joint	245	
ErgoArm	191			Joint Ball	232, 252	
ErgoArm Electronic plus	101			Joint Bolt	246	
ErgoArm Hybrid plus	103			Joint Clamp	232, 252	
ErgoArm plus	191			Joint Piece with Cable Lock	227, 241	
Extension Cable	127			Joint Piece without Lock	241	
Eyelet Cable Anchor	325 f, 329	Joint Plate		246		
8	F			Joint Protector	233	
	Finger-Thumb Group	366 f		8	K	
	Finger and Thumb Tip	142, 153, 209	Knurled Plate		221	
	Finger Cover	36, 159, 164, 370				
	Finger Filler	36				

L			
Lamination Disk	329	MyoEnergy Integral	108
Lamination Plate	90	MyoHand VariPlus Speed	43 ff
Lamination Protection Cover	194	MyolinoLink	136
Laminationring	178	MyolinoSkin Natural	76
Lamination Ring 90, 144, 154, 189, 192, 199, 205, 224, 227, 235, 251		MyolinoSoft	136
Lamination Ring for Children's Forearm	189	MyolinoWrist 2000	84
Lamination Ring for Physolino Baby Hand	306	MyoRotronic	92
Lamination Ring Set	188	MyoSelect	356
Laminationring Set	86	MyoSkin Natural	77 f
Leather Lamella	327	MyoWrist 2Act	85
Li-Ion Charger	113	MyoWrist Transcarpal	84
Light Metal Tube	225, 227, 234		
Locking Latch	246	N	
Locking Lever	112	Nylon wire	189, 197, 200, 319, 322, 326
Locking unit	371	O	
Lock Nut	360	O-Ring	195
Lock Set	310	O-Ring set	205
Lock Washer	247	Open-end Wrench	344
Lock with Release Pin	310	Ottobock Ball Shoulder Joint	232, 252
Loctite® 241	340	Ottobock Ball Wrist Unit	179
Loctite® 243	340	Ottobock Modular Arm Component	224, 226, 228, 231
Loctite® 601	341	Ottobock Shoulder Joint	204, 229 f, 251
Loop	329	Ottobock System Hand -passive-	141, 208
M		Ottobock System Hand -voluntary opening-	145
Michelangelo Hand	25	Ottobock Wrist Unit	179, 190
Michelangelo Hand with Axon-Bus Prosthetic System	24	Oval Head Screw	181
Modular Adapter for ArmLiners	221	P	
Modular Elbow Joint	224, 227, 233	Pad Screw Head	329
Mounting Aid	346	PAULA	134
Mounting Flange	236	PAULA with MyoBoy	134
Mounting Plate	346	Phillips Countersink Screw Head	246
Mounting Tool	334	Physolino Babyhand	303
Mounting Vice	334	Pin	311 f
Mounting Wrench	334	Pinch Gauge	64, 347
Movable All-purpose Ring	168	Plaster Shirt	337
MovoHook 2Grip	160 f	Plaster Sock	337
MovolinoArm Friction	188	Plug Bushing	368
MovoPAULA	136	Plug Housing	361
MovoShoulder Swing	204, 250	Pre-Shaped Foam Block	225, 228 f, 231
MovoWrist Flex	178	Precision Screwdriver Set	345
MyoBoy	135	Pressure Piece	194
MyoBoy Communication Cable	356	Pressure Ring	371
MyoCharge Integral	113	Pressure Switch	125
		Procomfort Gel	81, 339
		Prosthetic Glove for Adolescents	172, 214

1	Prosthetic Glove for Adolescents and Men	79	Special Cleaner	81, 341
	Prosthetic Glove for Children	79, 172, 214	Special Key	189
	Prosthetic Glove for Children and Adolescents	80	Special Lubricant	338
	Prosthetic Glove for Men	173, 215	Special Thread	64
	Prosthetic Glove for Women	80, 174, 216	Special Wrench	344 f
	Pull Cable	197 f, 241	Spiral Nut	325
	Pulling Unit	334	Spray for Donning Silicone Liners	312
2	Pulse Charger	113 f	Spray Pump Bottle	341
	PVC Connection Tube	336	Spring	164, 326
R				
3	Ratchet Ring with Inner Gear Teeth	371	Spring Telescope	193
	Reduction Gear	370	Stainless Steel Buckle	317, 321
	Refill Packge	312	Standard Connector	169, 184
	Resin Rapid Plumb	342	Stand for Soldering Iron	348
	Retainer Ring	142, 152, 209	Steel Cable	327
	Retaining ring	63	Strap Buckle	318, 321
	Ring	317, 321	Strap Clip	198
4	Rivet Pin	246	Suction Socket Electrode	120
	Rocker Switch	124	Suspension Rosette	153, 166
	Rubber Adhesive	342	Switch Cable	193
	Rubber Band	159	System Electric Greifer Digital Twin	70
	Rubber Friction Ring	181	System Electric Greifer DMC VariPlus	67
	Rubber Gripping Pad Set	72	System Electric Hand Digital Twin	55 ff
			System Electric Hand DMC plus	51 ff
5			System Electric Hand for other systems	61 f
			System Inner Hand	63 f, 142, 152, 209
			System Tester	356
			T	
	Screw Cap	341	TED & MyoBoy	137
	Screw Coupling	330	Test Adapter	354
	Screwdriver	345	Test Cable	352 ff
SensorHand Speed	39 ff	Test Plug	355	
6	Set Screw	144, 154, 235	ThermoLyn soft	338
	Setting Nut	319, 323	Threaded Fitting	326
	Shoulder Bracket Set	205	Threaded Stud	155, 167, 183, 306
	Side Cutting Pliers	343	Thread Segment	193
	Silicon ArmLiner	310	Thumb	364, 366
	Silicone Grease	338	Tongue Cap	194
	Single-Cable System Hand	146 ff	Tool Board	344
7	Sleeve Protection Pad	331	Transcarpal-Hand Digital Twin	60
	Slip Clutch	362	Transcarpal-Hand DMC plus	59
	Slotted Oval Head Screw	181	Transmission	362
	Socket Attachment Piece	336	Triple-Control Above-Elbow Harness	316
	Socket Screw with Allen head	318, 322, 329	Truss Head Screw	235
	Soldering Iron	347	Tube Chassis	360
	Soldering Jig	346	Tube chassis with lamination ring	365
8	Solder Lug	368	Tube Dummies	337
	Spacer Sleeve	318, 322		
	Spacing Washer	371		

Tube Valve for Suction Socket	336
Tweezers	64, 343 f

U

UHU	342
UHU-Plus Endfest 300 Adhesive	339
Universal Pliers	343
Universal Power Supply	115

V

Vacuum Forming Tool	335
Variotrac 68	339
Vice Clamping Jaws	335

W

Washer	166
Wire Coupling	326
Wood Hand Adapter	190, 306
Work Mat	342
Wrist Drive with Coaxial Plug	372
Wrist Lock with Standard Connector	169
Wrist Unit Lamination Ring	189
Wrist Unit without Ratchet	179 f
Wrist Unit – ratchet type rotation	180

X

X-ChangePack	109
--------------	-----

“

“D” Ring	329
----------	-----

1

2

3

4

5

6

7

8

	Article/Reference No	Page	Article/Reference No	Page	Article/Reference No	Page
1	1-10		8K500	24	8S4=228 x 84L	299
	4X74	115	8R1	93	8S4=228 x 88L	300
2	8E12	61	8S4=190x76	172, 214	8S4=228 x 89R	300
	8E32=6	71	8S4=199 x 82R	285	8S4=230 x 93R	301
	8E33=7	70	8S4=202 x 74L	284	8S4=232 x 94L	301
	8E33=9	67	8S4=203 x 83L	285	8S4=238 x 92L	302
	8E34=7	70	8S4=203 x 85L	286	8S4=244 x 94R	302
3	8E34=9	67	8S4=205 x 80R	288	8S5=165 x 72L	269
	8E37	62	8S4=205 x 81L	287	8S5=167 x 72R	269
	8E38=6	51	8S4=205 x 81R	287	8S5=174 x 74L	270
	8E38=7	55	8S4=206 x 76R	284	8S5=175 x 76R	270
	8E38=8	39	8S4=206 x 80L	288	8S5=176 x 80R	271
4	8E38=9	43	8S4=206 x 85R	286	8S5=180 x 80L	271
	8E39=6	52	8S4=206 x 87L	289	8S5=181 x 75R	275
	8E39=7	56	8S4=207 x 86L	290	8S5=182 x 84L	272
	8E39=8	40	8S4=208 x 85L	291	8S5=184 x 75L	273
	8E39=9	44	8S4=209 x 86R	290	8S5=184 x 78L	274
5	8E41=6	53	8S4=210x78	173, 215	8S5=184 x 78R	274
	8E41=7	57	8S4=211 x 88L	292	8S5=185 x 75L	275
	8E41=8	41	8S4=212 x 83R	291	8S5=186 x 92R	278
	8E41=9	45	8S4=212 x 86R	292	8S5=187 x 74R	273
	8E44=6	59	8S4=212 x 93L	293	8S5=187 x 79R	276
6	8E44=7	60	8S4=213 x 85L	294	8S5=188 x 79L	276
	8E51	35	8S4=214 x 82L	295	8S5=189 x 84R	280
	8E500	25	8S4=214 x 90R	297	8S5=190 x 77L	277
	8K5	303	8S4=215 x 83R	295	8S5=190 x 77R	277
	8K18	141, 208	8S4=215 x 88R	289	8S5=190 x 84R	272
7	8K19	141, 208	8S4=215 x 93R	293	8S5=190 x 93L	278
	8K20	145	8S4=218 x 83R	296	8S5=191 x 78R	279
	8K21	147	8S4=218 x 85L	296	8S5=192 x 78L	279
	8K22	146	8S4=218 x 85R	294	8S5=194 x 82L	280
	8K23	148	8S4=220 x 91L	297	8S5=195 x 79L	282
8	8K24	149	8S4=220x80	173, 215	8S5=195x78	174, 216
	8K26	150	8S4=221 x 81L	298	8S5=195x78L	281
	8K27	151	8S4=222 x 84R	299	8S5=195x78R	281
			8S4=225 x 82R	298	8S5=200 x 79R	282

Article/Reference No	Page	Article/Reference No	Page	Article/Reference No	Page
8S5=208 x 89L	283	8S7=212 x 93L	293	8S8=189 x 84R	280
8S5=210 x 89R	283	8S7=213 x 85L	294	8S8=190 x 77L	277
8S6=115 x 37R	261	8S7=214 x 82L	295	8S8=190 x 77R	277
8S6=115 x 38L	261	8S7=214 x 90R	297	8S8=190 x 84R	272
8S6=130 x 51L	262	8S7=215 x 83R	295	8S8=190 x 93L	278
8S6=134 x 52R	262	8S7=215 x 88R	289	8S8=191 x 78R	279
8S6=139 x 51R	263	8S7=215 x 93R	293	8S8=192 x 78L	279
8S6=142 x 50L	263	8S7=218 x 83R	294, 296	8S8=194 x 82L	280
8S6=151 x 58L	264	8S7=218 x 85L	296	8S8=195 x 78L	281
8S6=151 x 59R	264	8S7=220 x 91L	297	8S8=195 x 78R	281
8S6=158 x 54L	265	8S7=221 x 81L	298	8S8=195 x 79L	282
8S6=158 x 68R	266	8S7=222 x 84R	299	8S8=200 x 79R	282
8S6=159 x 53R	265	8S7=225 x 82R	298	8S8=208 x 89L	283
8S6=165 x 68L	266	8S7=228 x 84L	299	8S8=210 x 89R	283
8S6=166 x 70R	267	8S7=228 x 88L	300	8S9=115 x 37R	261
8S6=168 x 70L	267	8S7=228 x 89R	300	8S9=115 x 38L	261
8S6=170x65	172, 214	8S7=230 x 93R	301	8S9=130 x 51L	262
8S6=170x65L	268	8S7=232 x 94L	301	8S9=134 x 52R	262
8S6=170x65R	268	8S7=238 x 92L	302	8S9=139 x 51R	263
8S7=199 x 82R	285	8S7=244 x 94R	302	8S9=142 x 50L	263
8S7=202 x 74L	284	8S8=165 x 72L	269	8S9=151 x 58L	264
8S7=203 x 83L	285	8S8=167 x 72R	269	8S9=151 x 59R	264
8S7=203 x 85L	286	8S8=174 x 74L	270	8S9=158 x 54L	265
8S7=205 x 80R	288	8S8=175 x 76R	270	8S9=158 x 68R	266
8S7=205 x 81L	287	8S8=176 x 80R	271	8S9=159 x 53R	265
8S7=205 x 81R	287	8S8=180 x 80L	271	8S9=165 x 68L	266
8S7=206 x 76R	284	8S8=181 x 75R	275	8S9=166 x 70R	267
8S7=206 x 80L	288	8S8=182 x 84L	272	8S9=168 x 70L	267
8S7=206 x 85R	286	8S8=184 x 75L	273	8S9=170 x 65L	268
8S7=206 x 87L	289	8S8=184 x 78L	274	8S9=170 x 65R	268
8S7=207 x 86L	290	8S8=184 x 78R	274	8S11	79
8S7=208 x 85L	291	8S8=185 x 75L	275	8S11N	77
8S7=209 x 86R	290	8S8=186 x 92R	278	8S12	80
8S7=211 x 88L	292	8S8=187 x 74R	273	8S12N	78
8S7=212 x 83R	291	8S8=187 x 79R	276	8S13	80
8S7=212 x 86R	292	8S8=188 x 79L	276	8S13N	78

	Article/Reference No	Page	Article/Reference No	Page	Article/Reference No	Page
1	8S20	79	9E385	363	9S211	366
	8S20N	76	9E386	363	9S234=PAA	72
2	8S500	26	9E388	363	9S255	367
	8S501	25	9E397	86	9S256	367
	8S502	25, 26	9E420=L	130	9S258	90
	8X14	142, 152, 209	9E420=R	130	9S266	86
	8X18	63	9S6	142, 153, 209	9S278=PAA	73
3	8X24	64	9S10=40	153	9X14	124
	8Y1	64	9S15	142, 152, 209	9X18	124
	9E53	360	9S52=40	153	9X25	124
	9E68	360	9S65	143, 210	9X37	125
	9E83	360	9S67	143, 210	9X50/9X51	123
4	9E94	144, 154	9S72	364	9X52/9X53	123
	9E125=200	360	9S74	364	10A2	168
	9E161	341	9S83	364	10A3	168
	9E167	360	9S85	365	10A4	168
	9E168	361	9S96=40	144, 154	10A8	168
5	9E169	91	9S103=7 1/4	365	10A11	162
	9E185	127	9S103=7 3/4	365	10A12	163
	9E188	361	9S103=8 1/4	365	10A18	162
	9E189=1	361	9S110	90	10A22	169
	9E194	361	9S117	370	10A25	158
6	9E221	369	9S127=1	370	10A30	178
	9E235	369	9S136=1	371	10A37	158
	9E236=2	369	9S138	72	10A40	190, 306
	9E255=1	362	9S145	72	10A41	306
	9E264=7	370	9S149	72	10A43	155, 166, 183, 184
7	9E264=9	370	9S166=44-N	153	10A44	155, 167, 183
	9E274	362	9S166=48-N	153	10A56	156, 167, 183, 184
	9E342	36	9S184=44-N	144, 154	10A60	162
	9E346	370	9S184=48-N	144, 154	10A70	160
	9E347	36	9S185=44-N	153	10A71	160
8	9E369	130	9S185=48-N	153	10A80	161
	9E370	130	9S187	63, 142, 152, 209	10A81	161
	9E380	362	9S203	366	10R1	220
	9E381	363	9S204	366	10R2	220

Article/Reference No	Page	Article/Reference No	Page	Article/Reference No	Page
10R2=M12x1.5	232, 252	10Y26	324	12R2	226
10R3	220	10Y31	316, 320	12R3	228
10R4	221	10Y32	317, 320	12R4	231
10R5	234	11-20		12R5	226
10S1	90	11D1	178	12S4	204, 229, 251
10S4	91	11D12=34	189	12S5	204, 230, 251
10S16	90	11D20	181	12S6	204, 250
10S17	91	11D27	181	12S7	232, 252
10S18=40	189	11D31	306	12V10	336
10V8	180	11D61	86	13D1	205
10V9	179	11S2	371	13D2	205
10V10	180	11S6	371	13E29	371
10V15	169, 184	11S7	371	13E50	126
10V18	179	11S12	336	13E51=2	116
10V18=34	190	11S27	371	13E51=3	116
10V25	180	11S30	371	13E51=4	116
10V26	182	11S33	169	13E97	126
10V30	180	11S61	372	13E98=1200	126
10V32	179	12A13	192	13E99=1200	127
10V34	179	12K5	199	13E129	127
10V36	179	12K6	199	13E135	122
10V36=34	190	12K12	188	13E153	122
10V38	84	12K19=40	188	13E182	131
10V39	178	12K20	199	13E188	116
10V40	85	12K27	197	13E190	117
10V51	84	12K33	198	13E190=150	117
10V500	27	12K35	198	13E201	122
10Y1	159, 164	12K41	191	13E202	120
10Y2	164	12K42	191	13E205	92
10Y3	324	12K44	103	13E206	122
10Y8	159	12K48	192	13E500	28
10Y12	164	12K49	192	13G8=54	189
10Y13	159	12K50	101	13G8=67	199, 224, 227, 235
10Y19=1	330	12K100N	97	13G21	188
10Y19=2	330	12R1	224	13G37=N	233
10Y25	324			13G50	198

	Article/Reference No	Page	Article/Reference No	Page	Article/Reference No	Page	
1	13G68	192	16H1	328	21A16	329	
	13R1	224, 227, 233	16H2	328	21A17=3.8x1000	154, 326	
	13R3	225, 227, 234	16U4	243	21A18	189, 197, 200, 319, 322, 326	
	13R4	234	16U5	243	21A19	324	
	13R5	234	16U6	245	21A20	324	
	13R6	235	16U7	245	21A23	326	
	13R7	234	16U8	244	21A24	326	
2	13R8=64	236	16X4	244	21A25	326	
	13R9	225, 228, 229, 231	16X12	240	21A29	331	
	13R11	221	16X13	240	21A35=1	316	
	13X3	232, 252	16X14	241	21A36=1	320	
	13X4	232, 252	16Y5	246	21A37=1	316, 320	
	13X5	232, 252	16Y7	246	21A38	317, 321	
	13Y1	225, 228, 229, 236	16Y8	247	21A44	166	
3	13Z16=43	205, 251	16Y9	246	21A207	104, 196	
	13Z47	192	16Y26	197, 198, 241	21Y37	327	
	13Z48	193	16Y27	227, 241	21Y79	242	
	13Z50	193	16Y31	241	21Y120=47	327	
	13Z51	193	17F48=5	245	21Y194	317, 321	
	13Z52	193	17Y17=6x8x2.4	246	21Y195=25	317, 321	
	13Z53	193	17Y93=6x5.2xM4	246	21Y197=1	317	
4	13Z55	194	21-35			21Y197=2	321
	13Z56	194	21A3	325	21Y199	318, 321	
	13Z57	194	21A4	325	21Y203	318, 322	
	13Z58	194	21A5	325, 329	29C5	319	
	13Z59	194	21A6	325	29C5=M4x9	323	
	13Z68	104, 196	21A7	325	36-89		
	14A1	310	21A8	153, 166	60X5	135	
5	14A107	311, 312	21A9	166	60X6	136	
	14A110	310	21A10	166	90-99		
	14A111	312	21A11	326	99B13	336	
	14Y1	310	21A12	326	99B83	337	
	15K1=30	225	21A13	165			
	15K2=30	228	21A13=3	316			
	15K3=30	229, 231	21A14	165			
6	15Y1	331					
7							
8							

Article/Reference No	Page	Article/Reference No	Page	Article/Reference No	Page
100-499		610-619		700-799	
453A1	65	616T69	338	646C48	134
453H10	311	620-629		646C50=V1.2	136
453H12	311	623G4	327	646C51=V1.2	137
453H14	311	623G23	318, 322	646C52	134
500-599		623H2=25x5	327	651D4=2	327
501S7=M3x5	367	623H23	318, 322	700-799	
501S17=M3x5	367	624Z12	64	704B5	334
501S27	181	627F13=60x2.5	195	706F1	343
501S28	328	630-639		706R3	343
501S35=M4x6	235	633F11	338	706R4	343
501S40	181	633F37=0.02	339	706Z2	343
501S84=M4x20	194	633S2	81, 339	709G1=7x7	344
501S101=M4x12	195	634A58	339	709H1	344
501T8=M2x4	367	635P15	81	709S10=2	334
501T19=M5x6	246	636K11	340	709S42	86
501T52=M3x6	143, 155, 210	636K13	340	709Z6=1.8	344
501T52=M3x8	143, 155, 210	636K14	341	709Z7	345
501Z2=M4x18	235	636W22	342	709Z9	189
502S97=M5	368	636W23	339	710H3	345
503F3	318, 322, 329	636W34	342	710H4	345
504H3=11	330	636W60	340	711M1	30, 345
506A12	246	637L7	342	711M2	346
506G1=M6x10	235	640-699		711M3	346
506G4	144, 154	640F12	81, 341	711M5	346
507S15	329	640F13	341	711M7	346
514K3=27	330	640F18	312	711M12	346
514Z3=25	329	640F18=900	312	711M16	347
516S3	329	640Z6	342	711M17	334
519S5=0.8	368	641T9=M	337	711M18	334
520E30	368	641T10=M	337	711M20	334
560X3	136	642G1	342	711M50	334
560X500=V1.0	29	646C42	135	711M53	335
				711M56	335
				711M64	30, 345
				718H5	347

	Article/Reference No	Page	Article/Reference No	Page	Article/Reference No	Page
1	718Y1	347	757Z103=2	112		
	726W7=28.5	335	757Z149	110		
	726W9	335	757Z184=1	110		
2	736Y6	330	757Z184=2	110		
	743A18	336	757Z185=1	110		
	743A19	335	757Z185=2	111		
	743A23	104, 196	757Z186	112		
	743F1	64, 347	757Z190=1	111		
3	743Y42=24	337	757Z190=2	111		
	757B13	108	757Z191=1	111		
	757B15	109	757Z191=2	111		
	757B20/757B21	109	757Z195	112		
	757B35	108	799L2	348		
4	757B500	27	799P1	343		
	757L13	113	799P2	344		
	757L14	114	799Y1	344		
	757L16-2	115				
	757L20	113				
	757L35	113				
	757L500	28				
5	757M11=X-Change	135				
	757P2=4	352				
	757P2=6	352				
	757P2=8	353				
6	757P2=9	353				
	757P2=10	354				
	757P23	354				
	757P28	355				
	757P35	355				
7	757P39	356				
	757P41	127				
	757S1=AUS	115				
8	757S1=GB	115				
	757T13	356				
	757T16	356				
	757Z103=1	112				

1

2

3

4

5

6

7

8

Otto Bock HealthCare GmbH
Max-Näder-Str. 15 · 37115 Duderstadt/Germany
T +49 5527 848-1706 · F +49 5527 72330
export@ottobock.de · www.ottobock.com